Rapid Antibiotic Test

Microbiological test with spores of *Geobacillus stearothermophilus*
for antibiotic residue detection in milk, meat and other samples.

DESCRIPTION

Rapid Antibiotic Test is a cultural tube for broad spectrum antibiotic residue detection in foods such as milk and derivatives, meat, fish and other samples.

CONTENT OF THE PACKAGES

Each package contains 50 tubes containing 1.5 ml nutritive medium.

PRINCIPLE OF THE METHOD

Antibiotics are given to animals to prevent diseases and promote growth but are found as residue accumulated in the tissues of muscles, kidneys, liver and in biological fluids.

Rapid Antibiotic Test tubes contain a purple colored nutritive medium with spores of *G. stearothermophilus* ATCC® 10419.

After incubation at $60 \pm 2^\circ\text{C}$ for 4-6 hours the microorganism growth is indicated by the color change of the medium from purple to yellow. The presence of antibiotic residue in foods inhibits the spores germination so that the medium remains purple.

TEST PROCEDURE

Food sample preparation.

Solid sample (e.g. meat)

Triturate about 50 g sample. Put 2 g finely ground sample into a tube containing 6 ml distilled water (sample/water ratio, 1:3). Homogenize for 15 seconds. Centrifuge at 4000 rpm for 15 minutes. Use the supernatant to perform the test.

Liquid sample (e.g. milk)

Use the original sample.

Performing the test.

1. Take a Rapid Antibiotic Test tube and add 0.5 ml sample.
2. Incubate the tube at $60 \pm 2^\circ\text{C}$ for 4-6 hours.
3. Watch for the color changing and interpret the results.

INTERPRETATION OF THE RESULTS

The color change of the medium from purple to yellow indicates absence of antibiotic residue in concentrations equal to, or greater than the detection limits (Table 1). Purple color indicates presence of antibiotic residue.

Milk sample with
absence of antibiotic residue, after
incubation at $60 \pm 2^\circ\text{C}$ for 4 hours

Milk sample with
presence of antibiotic residue, after
incubation at $60 \pm 2^\circ\text{C}$ for 4 hours

Milk samples can react after 3 hours up to a maximum of 4 hours.
Meat samples can react after 5 hours up to a maximum of 6 hours.

Table 1. Rapid Antibiotic Test Sensitivity

Antibiotic	Detection limit (µg/l)	Antibiotic	Detection limit (µg/l)		
Penicillins	Ampicillin	6	Aminoglycosides	Neomycin	300
	Penicillin G	2		Kanamycin	400
	Cloxacillin	40		Gentamicin	100
	Oxacillin	10	Sulfonamides	Sulphadiazine	300
	Amoxicillin	6	Cephalosporins	Cephalexin	100
	Dicloxacillin	10	Nitrofurans	Nitrofurantoin	50
Tetracyclines	Oxytetracycline	250	Benzil Pirimidine	Trimethoprim	300
	Tetracycline	250	Macrolides	Erythromycin	100
	Doxycycline	100		Tylosin	150
			Lincomycin	400	

QUALITY CONTROL FOR THE USER

Each lot of Rapid Antibiotic Test is submitted to quality control using milk samples containing antibiotic such as Gentamicin and Penicillin G and milk samples without antibiotic residue.

PRECAUTIONS

For professional use only. The product must be used by properly trained operators with certain experience in the laboratory methods. It is recommended to consult the safety data sheet for its correct use.

STORAGE

Store at 2-8°C in its original package. In such conditions the product is valid until the expiry date shown on the label. Do not use beyond that date. Eliminate if signs of deterioration are evident.

DISPOSAL OF USED MATERIAL

After the use, the tubes and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.

REFERENCES

1. P. Aureli, A.M. Ferrini, V. Mannoni. Presumptive identification of sulphonamide and antibiotic residues in milk by microbial inhibitor tests. Food Control, 1996 7;3:165-8.
2. D. Currie, L. Lynas, D.G. Kennedy, W.J. McCaughey. Evaluation of a modified EC Four Plate Method to detect antimicrobial drugs. Food Additives and Contaminants, 1998, 15;6:651-60.
3. DM 26/03/1992 GU n° 90 16/04/1992: Attuazione della decisione n. 91/180/CEE concernente la fissazione di metodi di analisi e prova relativi al latte crudo e al latte trattato termicamente.
4. CONALAnalisi e controllo qualità dei prodotti lattiero-caseari Il latte VIII/5, 356-357 (1983).
5. 91/180/EEC: Commission Decision of 14 February 1991 laying down certain methods of analysis and testing of raw milk and heat-treated milk. Official Journal L 093, 13/04/1991 P. 0001 - 0048.
6. Van Oss J.L., Lameris S.A. Doodewaard J., Oostendrop J.G., Diffusion Test for the determination of antibiotic residues in milk. Neth. Milk Dairy J. 29, 16 (1975).
7. International Dairy Federation International Standard FIL-IDF 57:1970 Detection of penicillin in milk by a disk assay technique (1970).
8. Deborah Currie, Lorraine Lynas, D. Glenn Kennedy and W. John McCaughey. Evaluation of a modified EC Four Plate.
9. Method to detect antimicrobial drugs. Food Additives and Contaminants, 1998, Vol. 15, No. 6, 651-660.
10. Regulation 37/2010 EC.

PRESENTATION

Product	Ref.	Tests
Rapid Antibiotic Test	80351	50 tubes

TABLE OF SYMBOLS

LOT Batch code	 Do not reuse	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Consult instruction for use	

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
BIBLIOGRAFIA

1. P. Aureli, A.M. Ferrini, V. Mannoni. Presumptive identification of sulphonamide and antibiotic residues in milk by microbial inhibitor tests. Food Control, 1996 7;3:165-8.
2. D. Currie, L. Lynas, D.G. Kennedy, W.J. McCaughey. Evaluation of a modified EC Four Plate Method to detect antimicrobial drugs. Food Additives and Contaminants, 1998, 15;6:651-60.
3. DM 26/03/1992 GU n° 90 16/04/1992: Attuazione della decisione n. 91/180/CEE concernente la fissazione di metodi di analisi e prova relativi al latte crudo e al latte trattato termicamente.
4. CONAL Analisi e controllo qualità dei prodotti lattiero-caseari II latte VIII/5, 356-357 (1983).
5. 91/180/EEC: Commission Decision of 14 February 1991 laying down certain methods of analysis and testing of raw milk and heat-treated milk. Official Journal L 093 , 13/04/1991 P. 0001 - 0048.
6. Van Oss J.L., Lameris S.A. Doodewaard J., Oostendrop J.G., Diffusion Test for the determination of antibiotic residues in milk. Neth. Milk Dairy J. 29, 16 (1975).
7. International Dairy Federation International Standard FIL-IDF 57:1970 Detection of penicillin in milk by a disk assay technique (1970).
8. Deborah Currie, Lorraine Lynas, D. Glenn Kennedy and W. John McCaughey. Evaluation of a modified EC Four Plate.
9. Method to detect antimicrobial drugs. Food Additives and Contaminants, 1998, Vol. 15, No. 6, 651-660.
10. Regulation 37/2010 EC.

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Prodotto	Ref.	Test
Rapid Antibiotic Test	80351	50 provette

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