H₂O₂ / Catalase Reagent

Reagent for rapid detection of catalase enzyme.

INTENDED PURPOSE

Reagent for the execution of catalase test on bacterial isolates. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

H₂O₂ / Catalase Reagent is a reagent for microbiological use made of a H₂O₂ solution, for the performance of the catalase test.

KIT CONTENT

- 1 vial of H₂O₂ solution (5 ml)

METHOD PRINCIPLE

Hydrogen peroxide, final product of aerobic metabolism of carbohydrates, is extremely toxic for bacteria. Catalase-positive microorganisms, containing the enzyme catalase, break down hydrogen peroxide into water and oxygen.

Except for *Streptococci*, most of aerobe and facultative anaerobe microorganisms have catalase activity. Catalase reagent shows the presence of enzyme catalase with gas bubbles production.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, test tube, culture media, Catalase / Oxy Test, physiological solution (0.85% saline), quality control organisms.

REAGENTS

- A vial of H₂O₂ / Catalase Reagent contain 0.375 ml di H₂O₂ dissolved in 4.625 mL distilled water.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory.

Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Take the number of H₂O₂ / Catalase Reagent needed from the fridge and allow them to reach room temperature.
2. Add 4 drops of physiological solution (0.85 % saline) to each tube.
3. With a sterile loop pick up one or more than one well isolated colony from a fresh culture plate (18-24 hours old) and emulsify in the tube.
4. Add 3 drops of H₂O₂ Reagent.
5. Observe the almost immediate development of gas bubbles and the foam formation in the tube.

INTERPRETING RESULTS

The production of catalase enzyme is evidenced by the immediate bubbles and foam formation after the addition of H₂O₂ / Catalase Reagent to a suspension of the microorganism.

The absence of bubbles and foam formation indicates a microorganism unable to resolve H₂O₂, therefore lacking in catalase enzyme.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

3 years.

QUALITY CONTROL**Appearance of H₂O₂ / Catalase Reagent:** Limpid and clear solution, free of precipitate.

Control strain		Characteristic reactions
<i>Staphylococcus aureus</i>	WDCM 00034 ATCC® 25923 NCTC 12981	Positive reaction: Formation of bubbles of oxygen
<i>Campylobacter jejuni</i>	WDCM 00005 ATCC® 33291 NCTC 13367	Positive reaction: Formation of bubbles of oxygen
<i>Listeria monocytogenes</i>	WDCM 00021 ATCC® 13932 NCTC 10527	Positive reaction: Formation of bubbles of oxygen
<i>Listeria innocua</i>	WDCM 00017 ATCC® 33090 NCTC 11288	Positive reaction: Formation of bubbles of oxygen
<i>Listeria ivanovii</i>	WDCM 00018 ATCC® 19119 NCTC 11846	Positive reaction: Formation of bubbles of oxygen
<i>Enterococcus faecalis</i>	WDCM 00009 ATCC® 19433 NCTC 775	Negative reaction: No formation of bubbles of oxygen
<i>Enterococcus faecium</i>	WDCM 00177 ATCC® 6057	Negative reaction: No formation of bubbles of oxygen

PERFORMANCE CHARACTERISTICS

Performance testing of H₂O₂ / Catalase Reagent was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

H₂O₂ / Catalase Reagent may be used as an aid in the identification of microorganisms. Additional biochemical testing using pure cultures is recommended for complete identification.

Culture media that contain blood are unacceptable for this test due to the presence of catalase in erythrocytes. When picking colonies from blood media, avoid carryover when performing the slide test.

If a blood medium is used, a control slide catalase test should be performed where a small loopful of the blood-containing agar is tested with the reagent on the same slide as the organism. If the catalase reaction from the colony is much stronger than that from the agar alone, the test can be considered positive. Dirty glassware can cause false-positive results. Use properly cleaned tubes and slides. Inoculating needles or loops containing iron may produce false-positive results. Avoid contact of the catalase reagent with iron-containing needles and loops.

Bacteria (e.g., *Lactobacilli*) grown on media with low levels or no glucose may yield confusing reactions from pseudocatalase, a non-iron enzyme. The pseudocatalase reaction can be prevented by the addition of 1% glucose to the medium. Anaerobic cultures grown on a blood-free medium must be exposed to air for a minimum of 30 min before testing.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.

- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
H2O2 / Catalase Reagent	1 vial x 5 ml	80057

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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