# ESC swab (ISO 18593, ISO 17604) Neutralizing Rinse Solution

## DESCRIPTION

ESC (Easy Surface Checking) swab - Neutralizing Rinse Solution is a device used for the microbiological monitoring of surfaces by swab sampling method. Neutralizing agents included in the medium allow to assess the effectiveness of sanitization procedures in areas where residual disinfectants are presumably present.

## CONTENTS OF THE PACKAGES

The product consists of a transparent plastic tube containing transport medium and a sterile swab inserted in the cap (Fig 1).

Each package contains:

- 100 tubes with swab and liquid medium;
- 1 instructions sheet.

## METHOD PRINCIPLE

This solution contains the following neutralizers: histidine, lecithin, polysorbate 80 (Tween) and sodium thiosulphate. L-histidine inactivates aldehydes, lecithin neutralizes quaternary ammonium compounds, polysorbate 80 is effective against phenolic compounds and mercurial derivatives and sodium thiosulfate neutralizes halogen compounds.

## FORMULA (per liter of solution)

Casein peptone 1 g, histidine 1 g, lecithin 3 g, Tween 80 30 g, sodium thiosulphate 5 g, sodium chloride 8.5 g, pH 7.0±0.2 (25°C).

## TEST PROCEDURE

1. Take a tube from the package and unscrew the cap with the swab.
2. Press the tip of the swab against the wall of the tube to remove excess liquid.
3. Streak the swab on a surface horizontally and vertically. A Sampling Template 10x10 (ref. 96762) can be used to delimit and sample an area of 100 cm<sup>2</sup> (Fig 2).
4. Put the swab into the tube with the medium.
5. Screw the cap and record date and the sampling point.
6. Transport the device to the laboratory and examine according to methods in use.

Fig. 1



Fig. 2

## QUALITY CONTROL

Each batch of ESC swab - Neutralizing Rinse Solution is submitted to quality control using the following strains: *E. coli* ATCC® 25922, *S. aureus* ATCC® 6538.

## STORAGE

10-25°C away from light, until the expiry date on the label. Eliminate if signs of deterioration or contamination are evident.

## WARNING AND PRECAUTIONS

ESC Swab - Neutralizing Rinse Solution does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product must be used only by properly trained operators.

## BIBLIOGRAPHY

- ISO 18593:2018. Microbiology of the food chain - Horizontal method for surface sampling.
- ISO 17604:2015. Microbiology of the food chain - Carcass sampling for microbiological analysis.

## PRESENTATION

	Ref.	Contents
ESC swab - Neutralizing Rinse Solution	85601	100 x 10 ml tubes
ESC swab - Neutralizing Rinse Solution	85602	100 x 5 ml tubes

## TABLE OF SYMBOLS

<b>LOT</b> Batch code	Fragile, handle with care	Manufacturer	Use by
<b>REF</b> Catalogue number	Do not reuse	Contains sufficient for <n> tests	Caution, consult Instruction For Use

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