Indole Test Stick

Rapid test on stick for detection of tryptophanase enzyme

INTENDED PURPOSE

Rapid test for the Indole test of bacterial isolates. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Indole Test Stick is a biochemical test performed on bacterial species to determine the ability of the organism to convert tryptophan into indole.

KIT CONTENT

- 30 Stick of Indole Test Stick
- 1 tube of Indole Test Reagent (3.0 ml)

METHOD PRINCIPLE

Organisms that produce enzyme tryptophanase degrades the amino acid tryptophan into pyruvic acid, ammonia and indole. Indole is underlined by adding an aldehydic indicator that develops a strong green coloration. The test can be used to differentiate *Proteus* species and for presumptive identification of *Escherichia coli*.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, culture media, quality control organisms.

REAGENTS

- Indole Test Stick contains desiccated biochemical substrates and nutritive components.
- Indole Test Reagent consist of p-dimethylaminocinnamaldehyde dissolved in a solution of hydrochloric acid.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Take the number of Indole Test Stick from the container.
2. With a pointed loop, touch the colony to subject to test and slide on the stick's edge.
3. Add 2 drops of Indole Test Reagent.
4. Development of an immediate green coloration indicates positive reaction.

INTERPRETING RESULTS

A positive reaction is indicated by a colour change to green.

If the organism is Indole negative it is possible *Salmonella typhimurium* (Colorless).

Further tests are required for confirmation.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Appearance of Indole Test Stick: White

Appearance of Indole Test Reagent: Yellow-green liquid, clear

Control strain	Inoculum	Characteristic reactions
<i>Escherichia coli</i>	1-2 x 10 ⁸ CFU/ml	Positive reaction: Colour change to green
<i>Proteus mirabilis</i>		Negative reaction: Colorless
Salmonella Tiphymurium		Negative reaction: Colorless

PERFORMANCE CHARACTERISTICS

Performance testing of Indole Test Stick was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Allow at least 30 seconds for the color to develop before considering the indole test negative. Pigmented organisms may yield ambiguous results.

Indole production by certain anaerobes, such as *Clostridium* species, can be rapidly degraded resulting in false negative reactions.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Indole Test Stick	30 tests	88032

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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References / Riferimenti

1. EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. Public Health England. (2019). Urease test. UK Standards for Microbiology Investigations. TP 36 Issue 4. <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
3. Public Health England. (2018). Indole test. UK Standards for Microbiology Investigations. TP 19 Issue 4
4. Isemberg H.D. and Sundheim L.H.(1980) Indole reactions in bacteria. J. Bacteriol. 75:682-690. Baltimore, Williams & Wilkins.
5. MacFaddin J.F. (1980) Biochemical tests for identification for medicalbacteria. 2nd ed. 173-183. Baltimore, Williams & Wilkins.
6. Blazevic D.J. and Ederer, G.M.(1975) Principles of biochemical tests in diagnostic microbiology. 63-67. New York, John Wiley & Sons.

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