MAIA Pesticide MultiTest[®]

ENGLISH

Colorimetric system for the detection of organophosphate, organochloride and carbamate pesticides residues in soil, water, food and drink

PRINCIPLE OF THE METHOD

The method of Acetylcholinesterase – *Microplate Acetylcholinesterase Inhibition Assay (MAIA)* – for the detection of organophosphate, organochloride and carbamate pesticides residues in hydro-acetonitrilic extracts of solid/liquid food matrices is based on testing in microtitre plates of Acetylcholinesterase (AChE) activity inhibition by pesticide molecules, that belong to 3 important families. It is a semiquantitative analytic method called *screening method*, or first level method. The AChE is preincubated with the desiccated extract to accentuate the pesticide's inhibition effect on AChE, if present in the food matrix in examination. Afterwards, the enzyme-extract system is integrated with an appropriate reaction substratum, the acetylthiocholine, and with a chromogenic detector for thiocoline that has developed. The reaction is stopped by a denaturant of the enzymatic protein. Proceed to **direct visual evaluation** by comparing the color in Sample wells with the Color Card included in the kit or with the color of the Control Columns (ref. 79703). The method makes a biological test in which the analyte pesticide is selectively intercepted because of its specific “noxious” action on a critical physiological event of the animal organism, the AChE activity in nervous and neuromuscular junctions.

KIT CONTENTS

- 1 microtiter plate (96 wells);
- 1 sealing film
- 10 vials of freeze-dried **MAIA MEDIUM**, containing mammal AChE and a pH buffer;
- 10 vials of **MAIA STARTER** containing ATCl¹ and DTNB²;
- 1 bottles containing salts for matrix extraction (**MAIA Salts**);
- 1 basin for multichannel pipette.

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

- Adjustable pipette 1-5 mL.
- Adjustable multichannel pipette 50-300 µL.
- Chronometer.
- Mixer.
- Dispersing homogenizer for 10 mL tubes.
- Tabletop Centrifuge.
- Dryer.
- Shaker for microplate.
- Microplate reader with filter $\lambda = 405$ nm.

Other reagents as:

- Ultrapure water HPLC degree.
- Acetonitrile ACS degree, used as extractant and then as reaction STOPPER.

¹ Acetylthiocholine iodide

² 5,5'-Dithiobis (2-nitrobenzoic acid)

TEST PROCEDURE

The method of AChE inhibition assay (MAIA) for detection of organophosphate, organochloride and carbamate pesticides in water, food and drinks is characterized by the sequence of activities described below.

1. Preparation of liquid matrix Extract ($\Delta t \sim 20'$).
 - Introduce **4 mL** of liquid sample in a test tube.
 - Add **4 mL** of acetonitrile.
 - Put the cup to the tube and shake strongly for 1' roughly (by hand or *Vortex*).
 - Remove the cup and add **2 g** of salts.
 - Shake strongly for 1' roughly.
 - Centrifuge at 3000 rpm per 10'. The postcentrifugation supernatant is the Extract.
2. Preparation of the solid matrix Extract ($\Delta t \sim 25'$).
 - Mince by electric mixer an aliquot of **50 g**, statistically representative of the animal or vegetal, fresh or frozen matrix.
 - Introduce **4 g** of minced foodstuff, **1 mL** of ultrapure water and **4 mL** of acetonitrile into a 10 mL tube.
 - Homogenize by dispersing homogenizer for 30" roughly.
 - Add **2 g** of salts and put the cup to the tube.
 - Shake strongly for 1' roughly.
 - Centrifuge at 3000 rpm per 10'. The postcentrifugation supernatant is the Extract.
3. Dispense **100 μ L** of each Extract into Samples wells.
4. Desiccate the Extracts in wells under a direct air flow, at room temperature.
The desiccation time varies depending on laboratory temperature (*read paragraph Recommendations point 1*).
5. Prepare the **MEDIUM solution** by adding **3.5 mL** of ultrapure water to one MEDIUM vial.
One MEDIUM vial is sufficient for the test of 4 columns. Shake gently the vial until dissolution.
Dispense the MEDIUM solution into the basin for multichannel pipette.
6. Add **100 μ L** of **MEDIUM solution** by Multichannel pipette into all wells under examination.
7. Cover the microplate through the supplied sealing film.
8. Preincubate on the microplate shaker for 50' roughly at room temperature.
9. Prepare the **STARTER solution** by adding **2.5 mL** of ultrapure water to one vial of STARTER and mix until dissolution at room temperature. Avoid direct light. One STARTER vial is sufficient for the test of 4 columns. Dispense the STARTER solution into the basin for multichannel pipette.
10. Add **50 μ L** of **STARTER solution** into wells under examination. Dispense the STARTER solution into the columns by multichannel pipette. In the case more that one column are tested in the same session, keep 1' interval between each column. After dispensing the STARTER solution cover the microplate so that the reaction occurs in the dark.


EXAMPLE: aspirate STARTER by multichannel pipette, dispense into the first column, start chronometer, cover the microplate, start the Shaker. At 45" on chronometer stop the Shaker, aspirate Starter from the basin by multichannel pipette, at 60" dispense into the second column, cover the microplate, start the Shaker.
11. At 8' (*read paragraph Recommendations point 2*) from the addition of STARTER solution to the first column, dispense **100 μ L** of **STOPPER** into every well in examination. Aspirate and dispense the system until complete mixture.

EXAMPLE: At 7'45" stop the Shaker, aspirate STOPPER from he basin by multichannel pipette. At 8' stop the reaction by adding the STOPPER into the first column; aspirate and dispense the system until complete mixture, cover the microplate, start the Shaker. At 8'45" stop the Shaker, aspirate the Stopper; at 9' dispense the STOPPER into the second column, start the Shaker. (*The indicated times refer to a laboratory temperature of 25 °C. Read paragraph Recommendations point 2*).
12. Let the STOPPER act on the Shaker for 3'.
13. Determine the results of reaction by comparison with the Color Card or the Control Columns.

INTERPRETATION OF RESULTS

1. Direct visual evaluation of pesticide titers in the examined Extracts by comparing with the Color Card.
Compare the color in Sample wells with the yellow to colorless band: the yellow color intensity is inversely proportional to the pesticide quantity contained in the sample.
2. Direct visual evaluation of pesticide titers in the examined Extracts by comparing with the Control Columns (ref. 79703):
Compare the color in Sample wells with those in Control wells: a pesticide concentration estimate in the extract is obtained by the yellow colour intensity developed in the sample well, compared to the colour in control wells. The yellow colour intensity is inversely proportional to the pesticide quantity.

In the photo: Control Columns containing Paraoxon

MAIA Pesticide MultiTest was tested with milk containing the following pesticides:

Pesticide	Detection limit
Paraoxon	< 10 µg/kg
Carbaryl	< 100 µg/kg
DDT	< 100 µg/kg

The detection limits of analyte in matrix are below Italian and European **MRLs**.

RECOMMENDATIONS

1. The desiccation time varies depending on laboratory temperature: 15' at 25 °C, 30' at 20 °C, 45' at 18 °C, roughly. Check the laboratory temperature is 18-25 °C during the desiccation. Do not extend the desiccation phase over necessary time.
2. The enzyme activity varies depending on laboratory temperature. Consequently, set incubation time of **STARTER** (chromogenic reaction) according to as follows: 8' at 25 °C, 10' at 20 °C, 12' at 18 °C.
3. Avoid foam while resuspending **MEDIUM** with water. Shake gently the vial by circular movements.
4. Check the sealing film adhesion on wells after the addition of the **MEDIUM solution**.
5. Avoid direct light on **STARTER solution** in basin and wells during the color development reaction.
6. The mixture of **STOPPER** with the incubated solution in wells must be completely carried through 2 successive aspirations and dispensations at least for each well. Change the multichannel pipette tips after the stop of each column.
7. Adjust the multichannel pipette volume before the phases: from **MEDIUM solution** to **STARTER solution**: 100 µL → 50 µL; from **STARTER solution** to **STOPPER**: 50 µL → 100 µL.

QUALITY CONTROL FOR THE USER

Every batch of **MAIA Pesticide MultiTest** is submitted to quality control. The user can test negative and positive control by using the Control Columns (ref. 79703) which contain pesticide-free dried milk, and dried milk fortified with Paraoxon (organophosphate) and Carbaryl (carbamate).

PRECAUTIONS

Reagent may contain some non-reactive and preservative components.
It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general Good Laboratory Practice guidelines (GLP).

STORAGE

Store at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label.
Do not use beyond that date. Eliminate without using if there are signs of deterioration.

WASTE DISPOSAL

Product is intended for professional laboratories.
Waste products must be handled as per relevant security cards and local regulations.

REFERENCES

- Anastassiades, M., S. J. Lehotay, D. Stajnbaher and F. J. Schenck. (2003). *Fast and easy multiresidue method employing acetonitrile extraction/partitioning and dispersive solid-phase extraction for the determination of pesticide residues in produce*. **Journal of AOAC International** 86(2), 412-31
- Ellman, G.L., K.D. Courtney, V. Andres, and R.M. Featherstone. (1961). *A new and rapid colorimetric determination of acetylcholinesterase activity*. **Biochem. Pharmacol.** 7, 88-95.
- Galgani, F., and G. Bocquenne, (1989). *A method for routine detection of organophosphates and carbamates in sea water*. **Environm. Tech. Lett.** 10, 311-322.
- Mishra, N.N., J.A. Pedersen, and K.R. Rogers. (2001). *Highly sensitive assay for anticholinesterase compounds using 96 well plate format*. pp. 289-305. *In* Lipnick, R.L., R.P. Mason, M.L. Phillips, C.U. Pittnam, Jr. (eds.) *Chemicals in the Environment: Fate, Impacts and Remediation*, **ACS Symposium Series** No. 806; Oxford University Press: New York.
- Wilson, B.W, J.N. Seiber, M.E. Stelljes, J.D. Henderson, T.E. Archer, G.A. Pollock, and J.B. Knaak. (1989). *Bioassays for detection of aldicarb in watermelon*. **Bull. Environ. Contam. Toxicol.** 42, 159-166.

MAIA Pesticide MultiTest® trademark and the relevant process are covered by Italian (no. 0001332238) and European (Request no. 03425778.2) patents.

PRESENTATION

Product	REF	Kit contents
MAIA Pesticide MultiTest ®	79700	1 MAIA microplate 96 wells 10 MAIA <i>STARTER</i> vials 10 MAIA <i>MEDIUM</i> vials 1 bottle MAIA Salts

TABLE OF SYMBOLS

 Store away from light	 Do not reuse	 Manufacturer	 Contents of the package	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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
BIBLIOGRAFIA

- Anastassiades, M., S. J. Lehotay, D. Stajnbaher and F. J. Schenck. (2003). *Fast and easy multiresidue method employing acetonitrile extraction/partitioning and dispersive solid-phase extraction for the determination of pesticide residues in produce*. **Journal of AOAC International** 86(2), 412-31
- Ellman, G.L., K.D. Courtney, V. Andres, and R.M. Featherstone. (1961). *A new and rapid colorimetric determination of acetylcholinesterase activity*. **Biochem. Pharmacol.** 7, 88-95.
- Galgani, F., and G. Bocquenne, (1989). *A method for routine detection of organophosphates and carbamates in sea water*. **Environm. Tech. Lett.** 10, 311-322.
- Mishra, N.N., J.A. Pedersen, and K.R. Rogers. (2001). *Highly sensitive assay for anticholinesterase compounds using 96 well plate format*. pp. 289-305. In Lipnick, R.L., R.P. Mason, M.L. Phillips, C.U. Pittnam, Jr. (eds.) *Chemicals in the Environment: Fate, Impacts and Remediation*, **ACS Symposium Series** No. 806; Oxford University Press: New York.
- Wilson, B.W, J.N. Seiber, M.E. Stelljes, J.D. Henderson, T.E. Archer, G.A. Pollock, and J.B..Knaak. (1989). *Bioassays for detection of aldicarb in watermelon*. **Bull. Environ. Contam. Toxicol.** 42, 159-166.

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