Kit May-Grünwald Giemsa

Staining solutions for the preparation of cyto-histological samples for optical microscopy.

DESCRIPTION

May-Grünwald and Giemsa solutions are used for staining of blood and bone marrow smears and clinical-cytological specimens for optical microscopy.

KIT CONTENTS

Each kit contains:

- 1 bottle of May-Grünwald stain;
- 1 bottle of Giemsa stain;
- instruction sheet.

PRINCIPLE OF THE METHOD

The staining solutions contain methylene blue (basic dye) azures (basic dyes) and eosin (acid dye). The basic dyes carry net positive charge, consequently they stain nuclei, granules of basophil granulocytes and RNA molecules of the cytoplasm of white blood cells. The eosin carries net negative charge and stains red blood cells and granules of eosinophil granulocytes.

COMPOSITION

May-Grünwald is a methanol solution of 0.25% (w/v) eosin methylene blue.

Giemsa is a buffered methanol solution of 0.4% (w/v) azure eosin methylene blue.

PROCEDURE FOR USE

May-Grünwald and Giemsa solutions are supplied ready-to-use, although the Giemsa solution may be diluted 1:20 before use in either deionized water or in phosphate buffer solutions.

1. Place slides in May-Grünwald solution for 5 minutes.
2. Wash in running tap water for 1 minutes or place slides in phosphate buffer for 1.5 minutes.
3. Place slides in diluted or undiluted Giemsa solution for 15-20 minutes or 2-4 minutes, respectively.
4. Wash slides briefly in deionized water, dry with filter paper then in the air for 5 minutes.
5. Observe slides under microscope.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

The nuclei of white blood cells appear in blue, whereas the cytoplasm in light blue color. Basophils demonstrate dark blue black granules in the cytoplasm. The red blood cells appear red. Eosinophils demonstrate bright orange granules in the cytoplasm.

Note: the pH of washing water and dilution buffer as well as fixation and staining time can influence the staining results.

QUALITY CONTROL

1. Visual inspection:
 - May-Grünwald stain, clear blue solution;
 - Giemsa stain, clear blue solution.

WARNING AND PRECAUTIONS

The product is not classified as dangerous. it is recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* Diagnostic Use only and must be used by properly trained operators.

STORAGE

10-25°C away from light, until the expiry date on the label. Eliminate if signs of deterioration or contamination are evident.

REFERENCES

- E. V. Cowdry (1943) Microscopy Technique in Biology and Medicine. The Williams and Wilkins Company.
- P. L. Cardozo (1954) Clinical Cytology Using the May-Grunwald-Giemsa Stained Smear, Volume 1. Layden, L. Stafleu.

PRESENTATION

| Product | Ref. | Content |
|-------------------------|-------|--|
| Kit May-Grünwald Giemsa | 80282 | 1 x 250 ml bottle of May-Grünwald stain 1 x 250 ml bottle of Giemsa stain |

TABLE OF SYMBOLS

| | | | | |
|---|---|--|--|--|
|  Batch code |  <i>In vitro</i> Diagnostic Medical Device |  Manufacturer |  Contains sufficient for <n> tests |  Temperature limitation |
|  Catalogue number |  Fragile, handle with care |  Use by |  Caution, consult accompanying documents |  Do not reuse |

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LIOFILCHEM® S.r.l.