MiRA Test

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Rev.4 / 06.03.2013

Microbiological test with spores of *Geobacillus stearothermophilus* for the detection of antibiotic and sulphamide residues in milk.

DESCRIPTION

Antimicrobial substances are given to cattle for therapeutic treatment of infections and so they can be present in milk as residual drugs. The presence of residual antimicrobial drugs in milk is a potential hazard for the consumers since those substances interfere with the intestinal flora, can cause allergic reactions and can contribute to the spreading of resistant microorganisms.

Furthermore, the presence of antimicrobial drugs in milk is a serious technological problem for dairy industry, since antimicrobial substances interfere with milk fermentation which is the basis of many dairy productions.

The detection of inhibiting substances in milk is ruled in Europe by the Regulation 37/2010 EC.

MiRA (Milk Rapid Antibiotic) Test is a microbiological test, set for the detection of residual antibiotics and sulphamides in milk, that uses spores of *Geobacillus stearothermophilus* and a culture medium designed for that purpose.

CONTENT OF THE PACKAGE

Each package contains:

- n. 50 tubes containing liquid culture medium
- n. 50 discs containing spores of *Geobacillus stearothermophilus*
- 1 instructions sheet

PRINCIPLE OF THE METHOD

The test detects residues of antibiotics and sulphamides through their inhibition action against the growth *in vitro*, at $64 \pm 0,5$ °C, of *Geobacillus stearothermophilus*. The inhibition of bacterial growth is made evident by the presence of a fermentable sugar and a pH indicator in the culture medium.

PROCEDURE FOR USE

- Remove the cup from the tube and add 1 disc of spores. Shake gently.
- Put the cup and introduce the tube in a waterbath at $64 \pm 0,5$ °C or in a Termoblock set at $64 \pm 0,5$ °C. Incubate for 20'.
- Remove the tube from the incubator and let it reach room temperature (roughly 5').
- Remove the cup and introduce 1 mL of the milk under test.
- Prepare one or more analyte-positive controls, by the addition of antibiotic/sulphamide into the tubes as indicated in the paragraph "Quality Control for the user".
- Put the cup to the tube and let the antibiotic/sulphamide, if present, act at room temperature for 20'.
- Reintroduce the tube in the waterbath or in the Termoblock at $64 \pm 0,5$ °C for the second incubation for 4h - 6h.

INTERPRETATION OF THE RESULTS

At the end of the second incubation observe colour of the liquid in the tubes that, if green, indicates antibiotic/sulphamide residues in the examined milk in concentrations above the detection limits; otherwise, a yellow colour indicates the absence of antibiotic/sulphamide residues in the examined milk or their presence in concentrations under the detection limits.

Photo 1.

Photo 1.

Results obtained by **MiRA Test** after 3h 30' of incubation at $64 \pm 0,5$ °C:

- in the left side, no colour change of the pH indicator: analyte-positive milk sample containing antibiotic/sulphamide residues in concentration above the detection limits;
- in the right side, colour change of the pH indicator: analyte-negative milk sample containing no antibiotic/sulphamide residues or containing antibiotic/sulphamide residues in concentration under the detection limits.

USER QUALITY CONTROL

Every batch of **MiRA Test** is submitted to quality control. The analyte-negative control is obtained using a certified milk sample free of antibiotics. The analyte-positive controls are obtained using milk samples with the addition of antibiotics/sulphamides such as gentamicin (100 µg/kg), penicillin G (4 µg/kg) and sulphadiazine (150 µg/kg).

PRECAUTIONS

The product **MiRA Test** is not classified as dangerous under current legislation. **MiRA Test** is a disposable device to be utilized only for *in vitro* use, it is intended for use in a professional environment and must be used in laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store **MiRA Test** at 2-8 °C in its original packaging. Keep away from heat sources and avoid excessive temperature changes.

In such conditions **MiRA Test** will remain valid until the expiry date indicated on the label. Do not use beyond that date.

Eliminate without using it if there are deterioration signs.

DISPOSAL

After use and before disposal, tubes must be properly decontaminated with chemical substances or autoclaved.

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Table 1: Sensitivity of **MiRA Test** to some antibiotics/sulphamides, commonly used in veterinary medicine and respective MRL (Maximum Residue Limits) values valid for Europe¹.

Table 1

ANTIBIOTIC/SULPHAMIDE	MRL for milk ¹ MAXIMUM RESIDUE LIMITS (µg/kg)	MiRA Test Sensitivity in 6h DETECTION LIMITS (µg/kg)
BETA-LACTAMS		
Penicillin G	4	2-4
Ampicillin	4	2-4
Oxacillin	30	<10
Cloxacillin	30	<10
Dicloxacillin	30	<10
Amoxicillin	4	2-4
Benzylpenicillin	4	2-4
Nafcillin	30	15-30
Penethamate	4	2-4
Cefalexin	100	50-100
Cefazolin	50	25-50
Ceftiofur	100	50-100
Cephapirin	60	5-10
Cefquinome	20	10-20
TETRACYCLINES		
Tetracycline	100	50-100
Clortetracycline	100	50-100
Oxytetracycline	100	50-100
MACROLIDES		
Erythromycin	40	<10
Tylosin	50	25-50
Tilmicosin	50	20-40
Spiramycin	200	100-200
LINCOSAMIDES		
Lincomycin	150	75-150
Pirlimycin	100	<50
AMINOGLYCOSIDES		
Gentamicin	100	50-100
Neomycin	1500	<100
Streptomycin	200	<100
Dihydrostreptomycin	200	<100
SULPHAMIDES		
Sulfadiazine	100	<150
SULFANILAMIDES		
Sulfadimidine	100	<200
BENZIL PIRIMIDINE		
Trimethoprim	50	25-50
QUINOLONES		
Flumequine	50	50-100
Enrofloxacin	100	50-100
NOVOBIOCIN		
Novobiocin	50	100-200

¹ Regulation 37/2010 EC.

TABLE OF SYMBOLS

LOT Batch code	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	

LIOFILCHEM® S.r.l.

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Phone +390858930745 Fax +390858930330 Website: www.liofilchem.net E-Mail: liofilchem@liofilchem.net

Rev.4 / 06.03.2013