DESCRIZIONE

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 F13213
Rev.0 / 17.02.2005

DESCRIPTION

USP (United States Pharmacopoeia), EP (European Pharmacopoeia) and DAB (Deutsches Arzneibuch) standards recommend to use bioindicators during sterilization processes.

Biological indicators **SUSPENSION SPORE GST VHP E8/E7/E6/E5** are used for regular control of H₂O₂ vapour superficial sterilization cycles and validation of sterilizers. These bioindicators consist of *Geobacillus stearothermophilus* (ATCC 7953) spore suspension in distilled water in predefined concentrations. Biological indicators **SUSPENSION SPORE GST VHP E8/E7/E6/E5** are produced under strictly controlled conditions in order to satisfy the requirements indicated in the USP current edition and in accordance with ISO 11138 and EN 866 standards.

PRINCIPLE

SUSPENSION SPORE GST VHP E8/E7/E6/E5 contain *Geobacillus stearothermophilus* (ATCC 7953) spores in predefined concentrations: E8=1-5x10⁸CFU/ml; E7=1-5x10⁷CFU/ml; E6=1-5x10⁶CFU/ml; E5=1-5x10⁵CFU/ml.

Spores are completely killed off if the sterilization cycle has been efficient. In this case, during the following incubation of biological indicator in a suitable medium, spores are not able to grow and to modify the medium's aspect. On the contrary if the sterilization cycle has not been efficient, spores partially survive and therefore during the incubation are able to grow and to modify the medium's aspect.

INSTRUCTIONS FOR USE

Take a suitable volume of the suspension and put it on the selected material. Protect appropriately the inoculated material in order to avoid spore dispersion.

Exposure of biological indicators

Put the biological indicators, prepared from the **SUSPENSION SPORE GST VHP E8/E7/E6/E5**, inside the sterilization environment in the points considered the most difficult to reach and to sterilize with the hydrogen peroxide vapour.

Follow the operative instructions indicated by sterilizer manufacturer for the sterilization process. After sterilization, during the biological indicators removal, it is important to work aseptically in order to avoid possible contaminations that could invalidate results.

Biological indicators cultures

At the end of sterilization/aeration process, remove the biological indicators from the sterilizers and deliver them to the laboratory for the culture together with a not exposed control (positive control), prepared with the others starting from the **SUSPENSION SPORE GST VHP E8/E7/E6/E5** belonging to the same batch. In optimal conditions all culture procedures have to be performed under a laminar flow cabinet following aseptic methods. Submit to culture the biological indicators inoculating them in a suitable medium container. Ensure that the biological indicators are completely dipped in the medium gently shaking the container.

Incubation

Incubate the containers with the biological indicators at 55-60 °C (131°-140 °F) for 7 days or for a shorter time validated by user.

EVALUATION

Geobacillus stearothermophilus (ATCC 7953) spores are killed off if the sterilization cycle has been efficient: in this case the medium remains unchanged even after incubation at 55-60 °C (131-140 °F) for the selected time.

If the sterilization cycle has not been efficient, spores partially survive and the medium's aspect change after incubation at 55-60 °C (131-140 °F) for the selected time. The medium in the container inoculated with the biological indicator not submitted to the sterilization cycle and used as spore growth control (positive control), has to change aspect after incubation. On the contrary, the test must be repeated after having investigated causes of the negative result.

EVALUATION TABLE		
MEDIUM'S ASPECT	SPORE	STERILIZATION
Unchanged	Killed off	Successful
Changed	Vital	Unsuccessful

MATERIAL TREATMENT AFTER USE

After use, sterilize the positive containers (changed medium's aspect) in autoclave at 121 °C for at least 30 minutes and eliminate them in accordance with the procedures of the laboratory.

STORAGE

Store the product at 2-8 °C: in these conditions it maintains its validity until the expiry date indicated on the label.

BIBLIOGRAPHY

- United States Pharmacopoeia latest edition.
- Deutsches Arzneibuch latest edition.
- European Pharmacopoeia latest edition.
- ISO 11138 and EN 866 latest edition.

PRESENTATION

PRODUCT	CODE	PACKAGING	Spore: CFU / ml	D _{VHP} (1,6±0,5mg/L, 30°C)
SUSPENSION SPORE GST VHP E8	82030	10 ml vial	1-5x10 ⁸	1,5-4,0 minutes
SUSPENSION SPORE GST VHP E7	82031	10 ml vial	1-5x10 ⁷	1,5-4,0 minutes
SUSPENSION SPORE GST VHP E6	82032	10 ml vial	1-5x10 ⁶	1,5-4,0 minutes
SUSPENSION SPORE GST VHP E5	82033	10 ml vial	1-5x10 ⁵	1,5-4,0 minutes

TABLE OF SYMBOLS

 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Caution, consult accompanying documents
 Use by	 Batch code	 Do not reuse

LIOFILCHEM Bacteriology Products

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745

Fax +390858930330

Website: www.liofilchem.net

E-Mail: liofilchem@liofilchem.net

F13213

Rev.0 / 17.02.2005