STERILtest BAT E6/E5 self-contained

Biological indicators of ethylene oxide sterilization processes containing *Bacillus atropheus* (ATCC 9372) spores self-contained.

DESCRIPTION

USP (United States Pharmacopoeia) and EP (European Pharmacopoeia) standards, recommend to use bioindicators during ethylene oxide sterilization processes.

STERILtest BAT E6/E5 self-contained are made of a thermoplastic tube, with a cup permeable to ethylene oxide, that contains an ampoule with a nutrient broth with a pH indicator and a disc impregnated with *Bacillus atropheus* (ATCC 9372) spores in predefined concentrations (E6=1-5x10⁶UFC/disc; E5=1-5x10⁵UFC/disc).

These bioindicators are used for validation, re-validation and process monitoring of ethylene oxide sterilizers.

Biological indicators **STERILtest BAT E6/E5 self-contained** are produced under strictly controlled conditions in order to satisfy the requirements indicated in the USP current edition and in accordance with ISO 11138 and EN 866 standards.

PRINCIPLE

STERILtest BAT E6/E5 self-contained contains, inside a thermoplastic tube, a paper disc impregnated with *Bacillus atropheus* (ATCC 9372) spores in predefined concentrations (E6=1-5x10⁶UFC/disc; E5=1-5x10⁵UFC/disc).

The tube contains also a small glass ampoule with nutrient broth with a pH indicator. Spores are completely killed off during the ethylene oxide sterilization process if the process has been efficient. Therefore during the following incubation in the culture medium, contained inside the ampoule, they are not able to grow and to modify the medium's colour.

On the contrary if the sterilization process has not been efficient, spores partially survive and therefore during the incubation are able to grow and to modify the medium's colour.

INSTRUCTIONS FOR USE

- Take one or more **STERILtest BAT E6/E5 self-contained** from the package.
- Put the **STERILtest BAT E6/E5 self-contained** on the bottom, in the centre, inner sides and on the critical points of the sterilizer.
- For sterilizer with capacity up to 250 litres put two **STERILtest BAT E6/E5 self-contained** for each selected point of the sterilizer. For sterilizers with capacity higher than 250 litres put six or more **STERILtest BAT E6/E5 self-contained** in each selected point.
- Remove the **STERILtest BAT E6/E5 self-contained** after sterilization/aeration cycle, and activate them by breaking the glass ampoule through pressure to permit contact between culture medium and spores.
- Incubate the activated **STERILtest BAT E6/E5 self-contained** at 30-35 °C (86-95 °F) for 7 days or for a different time validated by user.
- Incubate, at the same conditions of time and temperature, a **STERILtest BAT E6/E5 self-contained**, activated by breaking the glass ampoule through pressure, not submitted to the sterilization cycle and belonging to the same batch, as spore growth control (positive control).
- Examine **STERILtest BAT E6/E5 self-contained** colour and interpret results as per **EVALUTATION TABLE**: a change of medium's colour from violet/clear to yellow/turbid indicates a microbial growth and therefore an unsuccessful sterilization. On the contrary, the persistence of the medium's initial colour (violet/clear) indicates absence of microbial growth and therefore a successful sterilization.

EVALUATION

Bacillus atropheus (ATCC 9372) spores are killed off if the sterilization cycle has been efficient: in this case the broth contained in **STERILtest BAT E6/E5 self-contained** remains violet/clear even after incubation at 30-35 °C (86-95 °F) for the selected time. If the sterilization cycle has not been efficient, spores partially survive and the broth contained in **STERILtest BAT E6/E5 self-contained** turns yellow/turbid after incubation at 30-35°C (86-95°F) for the selected time. The broth contained in **STERILtest BAT E6/E5 self-contained** not submitted to the sterilization cycle and used as spore growth control has to turn yellow/turbid after incubation. On the contrary, the test must be repeated after having investigated the causes of the negative result.

EVALUATION TABLE		
MEDIUM COLOUR	SPORE	STERILIZATION
Violet / Clear	Killed off	Successful
Yellow / Turbid	Vital	Unsuccessful

STERILTEST TREATMENT AFTER USE

After use, sterilize the positive **STERILtest BAT E6/E5 self-contained** (yellow/turbid) in autoclave at 121 °C for at least 30 minutes and eliminate them in accordance with the procedures of the laboratory.

STORAGE

Store the product at 10-25 °C: in these conditions it maintains its validity until the expiry date indicated on the label.

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D value for ethylene oxide is calculated with the MPN and SC methods with 600±30 mg/L of ethylene oxide, at a temperature of 54±1 °C, with 50-70% humidity.

TABLE OF SYMBOLS

 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Caution, consult accompanying documents
 Use by	 Batch code	 Do not reuse

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LIOFILCHEM® S.r.l.