(Strips)

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Rev.0/ 22.10.2007

(Strips)

Biological indicators of dry heat or ethylene oxide sterilization processes containing *Bacillus atropheus* (ATCC 9372) spores inoculated on special filter paper strips

DESCRIPTION

USP (United States Pharmacopoeia), EP (European Pharmacopoeia) and DAB (Deutsches Arzneibuch) standards, recommend to use bioindicators during dry heat and ethylene oxide sterilization processes.

Biological indicators **STRIP CONTROL BAT E2/E3 (Strips)** are special filter paper strips, inoculated with *Bacillus atropheus* (ATCC 9372) spores in predefined concentrations and contained in a special envelope. These bioindicators are used for validation, re-validation and process monitoring of dry heat and ethylene oxide sterilizers.

Biological indicators **STRIP CONTROL BAT E2/E3 (Strips)** are produced under strictly controlled conditions in order to satisfy the requirements indicated in the USP current edition and in accordance with ISO 11138 and EN 866 standards.

PRINCIPLE

STRIP CONTROL BAT E2/E3 (Strips) contains *Bacillus atropheus* (ATCC 9372) spores in predefined concentrations: E2=1-5x10²CFU/strip; E3=1-5x10³CFU/strip.

Spores in the strips are completely killed off during dry heat or ethylene oxide sterilization process, if the process has been efficient. In this case the aspect of the medium inoculated with the strips and incubated for a suitable time will remain unchanged.

On the contrary, if the sterilization process has not been efficient, spores partially survive causing a change in the aspect of the nutrient medium used for evaluation of microbial growth.

INSTRUCTIONS FOR USE

- Take one or more strips from the kit leaving them inside their original envelope.
- Put envelopes with strips on the bottom, in the centre, inner sides and on the critical points of the sterilizer.
- For sterilizers with capacity up to 250 litres put two envelopes for each selected point of the sterilizer. For sterilizers with capacity higher than 250 litres put six or more envelopes in each selected point.
- Remove envelopes after sterilization/aeration cycle and open them aseptically with sterile scissors or by tearing the edges.
- Transfer aseptically each strip from its envelope to a Tryptic Soy Broth tube (cod. 21113) or a Steri-Test Medium tube (cod. 20199).
- Incubate tubes containing strips at 30-35 °C (86-95 °F) for 7 days or for a shorter time validated by user.
- Incubate, at the same conditions of time and temperature, a strip contained in the envelope not submitted to the sterilization cycle, belonging to the same batch, as spore growth control (positive control).
- Examine tube medium's aspect and interpret results as per **EVALUATION TABLE**: the presence of turbidity or a sediment in TSB or a change of medium's colour from violet/clear to yellow/turbid in Steri Test Medium indicates a microbial growth and therefore an unsuccessful sterilization. On the contrary, the persistence of the medium's clearness or initial colour indicates absence of microbial growth and therefore a successful sterilization.

EVALUATION

Bacillus atropheus (ATCC 9372) spores are killed off if the sterilization cycle has been efficient: in this case Tryptic Soy Broth remains clear while the colour of the Steri Test Medium remains violet/clear, even after incubation at 30-35 °C (86-95 °F) for the selected time.

If the sterilization cycle has not been efficient, spores partially survive and the content of the TSB tubes becomes turbid while the Steri Test Medium turns yellow/turbid after incubation at 30-35 °C (86-95 °F) for the selected time.

The tube inoculated with the strip contained in the envelope, not submitted to the sterilization, has to show the growth of the spores. On the contrary, the test must be repeated after having investigated causes of the negative result.

| EVALUATION TABLE | | | |
|-------------------------|-------------------|------------|---------------|
| TRYPTIC SOY BROTH | STERI TEST MEDIUM | SPORES | STERILIZATION |
| Clear | Violet | Killed off | Successful |
| Turbid or with sediment | Yellow | Vital | Unsuccessful |

TREATMENT OF TEST TUBES AFTER USE

After use, sterilize the positive tubes (turbid/ yellow) in autoclave at 121 °C for at least 30 minutes and eliminate them in accordance with the procedures of the laboratory.

STORAGE

Store the product at 2°C -8 °C: in these conditions it maintains its validity until the expiry date indicated on the label.

BIBLIOGRAPHY

1. United States Pharmacopoeia latest edition.
2. Deutsches Arzneibuch latest edition.
3. European Pharmacopoeia latest edition.
4. ISO 11138 and EN 866 latest edition.

PRESENTATION

| PRODUCT | CODE | PACKAGING | Spores: CFU/strip | D _{EO} (600 ± 30 mg/L, 60% ± 10% RH, 54 ± 1 °C) | D _{DH} (160 ± 1 °C) |
|-------------------------------|-------|-----------|----------------------|---|---------------------------------|
| STRIP CONTROL BAT E2 (Strips) | 91169 | 100 buste | 1-5x10 ⁶ | 2,6-4,5 minutes | 1,0-3,0 minutes |
| STRIP CONTROL BAT E3 (Strips) | 91166 | 100 buste | 1-5x10 ⁶ | 2,6-4,5 minutes | 1,0-3,0 minutes |

D value for ethylene oxide (D_{EO}) is calculated with the MPN and SC methods with 600±30 mg/L of ethylene oxide, at a temperature of 54±1°C, with 50-70% humidity.
D value for dry heat (D_{DH}) is calculated with MPN method at a temperature of 160±1°C.

TABLE OF SYMBOLS

| | | |
|---|---|---|
|  Manufacturer |  Contains sufficient for <n> tests |  Temperature limitation |
|  Catalogue number |  Fragile, handle with care |  Caution, consult accompanying documents |
|  Use by |  Batch code |  Do not reuse |

LIOFILCHEM Bacteriology Products

Via Scozia- Zona Ind.le - 64026 ROSETO D.A. (TE) ITALY

Tel.+39 085 8930745 - Fax +39 085 8930330 Web site: <http://www.liofilchem.net> E-Mail: liofilchem@liofilchem.net

Rev.0/ 22.10.2007