(Strips)

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Rev.0 / 22.10.2007

(Strips)

Biological indicators of irradiation sterilization processes containing
Bacillus pumilus (ATCC 27142) spores inoculated on special filter paper strips

DESCRIPTION

USP (United States Pharmacopoeia), EP (European Pharmacopoeia) and DAB (Deutsches Arzneibuch) standards, recommend to use bioindicators during irradiation sterilization processes.

Biological indicators **STRIP CONTROL PUM E6 (Strips)** are special filter paper strips, inoculated with *Bacillus pumilus* (ATCC 27142) spores in predefined concentrations and contained in a special envelope. These bioindicators are used for validation, re-validation and process monitoring of irradiation sterilizers.

Biological indicators **STRIP CONTROL PUM E6 (Strips)** are produced under strictly controlled conditions in order to satisfy the requirements indicated in the USP current edition and in accordance with ISO 11138 and EN 866 standards.

PRINCIPLE

STRIP CONTROL PUM E6 (Strips) contain *Bacillus pumilus* (ATCC 27142) spores in predefined concentrations: E6=1-5x10⁶CFU/strip.

Spores in the strips are completely killed off during irradiation sterilization process, if the process has been efficient. In this case the aspect of the medium inoculated with the strips and incubated for a suitable time will remain unchanged. On the contrary, if the sterilization process has not been efficient, spores partially survive causing a change in the aspect of the nutrient medium used for evaluation of microbial growth.

INSTRUCTIONS FOR USE

- Take one or more strips, from the package leaving them inside their original envelope.
- Put a suitable number of envelopes with the strips on the critical points of material to sterilize.
- Remove envelopes after sterilization cycle and open them aseptically with a steril scissors or by tearing the edges.
- Transfer aseptically each strip from its envelope to a Tryptic Soy Broth tube (cod. 21113) or a Steri-Test Medium tube (cod. 20199).
- Incubate tubes containing strips at 30-35 °C (86-95 °F) for 7 days or for a shorter time validated by user.
- Incubate, at the same conditions of time and temperature, a strip contained in the envelope not submitted to sterilization cycle, belonging to the same batch, as spore control growth (positive control).
- Examine tube medium's aspect and interpret results as per **EVALUATION TABLE**: the presence of turbidity or a sediment in TSB or a change of medium's colour from violet/clear to yellow/turbid in Steri Test Medium indicates a microbial growth and therefore an unsuccessful sterilization. On the contrary, the persistence of the medium's clearness or initial colour indicates absence of microbial growth and therefore a successful sterilization.

EVALUATION

Bacillus pumilus (ATCC 27142) spores are killed off if the sterilization cycle has been efficient : in this case Tryptic Soy Broth remains clear while the colour of the Steri Test Medium remains violet/clear, even after incubation at 30-35 °C (86-95 °F) for the selected time.

If the sterilization cycle has not been efficient, spores partially survive and the content of the TSB tubes becomes turbid while the Steri Test Medium turns yellow/turbid after incubation at 30-35 °C (86-95 °F) for the selected time.

The tube inoculated with the strip contained in the envelope, not submitted to the sterilization, has to show the growth of the spores. On the contrary, the test must to be repeated after having investigated causes of the negative result.

EVALUATION TABLE			
TRYPTIC SOY BROTH	STERI TEST MEDIUM	SPORES	STERILIZATION
Clear	Violet / Clear	Killed off	Successful
Turbid or with sediment	Yellow / Turbid	Vital	Unsuccessful

TREATMENT OF STERI-TEST MEDIUM TEST TUBES AFTER USE

After use, sterilize the positive tubes (turbid/yellow) in autoclave at 121 °C for at least 30 minutes and eliminate them in accordance with procedures of the laboratory.

STORAGE

Store the product at 2-8 °C: in these conditions it maintains its validity until the expiry date indicated on the label.

BIBLIOGRAPHY

- United States Pharmacopoeia latest edition.
- Deutsches Arzneibuch latest edition.
- European Pharmacopoeia latest edition.
- ISO 11138 and EN 866 latest edition.

PRESENTATION

PRODUCT	CODE	PACKAGING	Spores: CFU / strip	D _{kGy}
STRIP CONTROL PUM E6 (Strips)	91172	100 envelopes	1-5x10 ⁶	1,0-2,0 kGy

D value for irradiation is calculated in accordance with USP method

TABLE OF SYMBOLS

 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Caution, consult accompanying documents
 Use by	 Batch code	 Do not reuse

LIOFILCHEM Bacteriology Products

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745

Fax +390858930330

Website: www.liofilchem.net

E-Mail: liofilchem@liofilchem.net

Rev.0 / 22.10.2007