# DISC DISPENSER 8

Antibiotic discs dispenser

Ref. 91200

Contents	page
Italiano	1
English	4

Rev.1 / 03.10.2012

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**ENGLISH**

### DESCRIPTION

DISC DISPENSER 8 is a device that allows to place up to 8 antibiotic discs on the agar surface of Petri dishes.

### FEATURES


1. Airtight container with desiccant.  
Label for the identification of antibiotics positioned on top of the container.

2. Block-cartridges system positioned at the base of the plunger.


3. Height adjustable base (1 cm  $\updownarrow$ ).

### Discs placing


- The discs are positioned at equidistance along a circumference of about 16 cm.

- Each disk is pressed on the agar surface with a pressure sufficient to ensure a homogeneous diffusion of the antibiotic.

## HOW TO USE

1. Adjust the height of the base according to the thickness of the agar.

2. Insert the cartridge into the guides by sliding them until you hear a "click".

3. Turn the block-cartridges system positioned on the plunger.

4. Place the dispenser over the Petri dish.

5. Push the plunger down until it stops with a firm, quick motion. Then release the plunger.

6. Before proceeding to a further distribution of discs on a new plate ensure the presence of other cartridges inside the dispenser.

**STORAGE**

Store DISC DISPENSER 8 in its container with desiccant provide to protect from moisture. With cartridges inside store the dispenser at 2-8°C.

The desiccant is active whilst pale yellow or white in color.

When green regenerate as indicated below:

1. Dry in a hot air oven at 50-55°C for 8-12 hours (or in a vacuum oven for 2-3 hours).
2. Ensure that the desiccant is pale yellow after drying before re-use.
3. Discard desiccant when a pale yellow color is not achieved after drying.

Before cleaning, always remove all cartridges from the dispenser.

Use only a solution of 2% quaternary ammonium salts and a solution of 1% sodium hypochlorite.

Clean the bottom of the dispenser and the ejection holes of the discs using a swab.

**DECONTAMINATION**

Using properly DISC DISPENSER 8 according to the instructions in this manual no component should become contaminated.

In case of accidental contamination, follow the instructions below:

1. Remove all cartridges from the dispenser.
2. Immerse the dispenser until the upper end of the base ring in a solution of 2% quaternary ammonium salts for 8-12 hours.
3. Rinse the dispenser with distilled water and dry thoroughly.

**TROUBLESHOOTING GUIDE**

If the dispenser fails to dispense, not push the plunger further and proceed to check in the order specified below:

1. Check for empty cartridges (replace as necessary).
2. Check that the cartridge is fitted correctly.
3. Remove all cartridges.
4. Check for discs trapped in dispenser (avoid to use metal implements)
5. Check for faulty discs remaining in cartridges.
6. If the problem persists, contact your local dealer, or [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

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Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.