EASY RID MULTIPLATE

DESCRIZIONE

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EASY RID MULTIPLATE

DESCRIPTION

EASY RID MULTIPLATE are plates for quantitative determination in radial immunodiffusion of human plasma proteins in serum and plasma. One plate allows to detect more plasma proteins. The simultaneous determination is permitted by the presence, in the same plate, of different gels, distinguished by the color and specific for the protein to search. EASY RID MULTIPLATE made by Liofilchem s.r.l. are in 4 configurations.

CONTENT OF THE KIT

Each kit contains a plate and an instructions sheet.

MATERIAL REQUIRED BUT NOT PROVIDED

- Precision pipette - Suitable measuring device (mm)

CONFIGURATION

EASY RID MULTIPLATE proposed by Liofilchem S.r.l. are available in the following configurations:

Table n°1

DESCRIPTION - CODE	PLASMA PROTEIN	GEL COLOR
h-IgG / IgA / IgM Code 93104	IgG IgA IgM	green red blue
h-C3c / C4 Code 93106	C3c C4	yellow red
h-ApoA1 / ApoB Code 93110	Apolipoprotein A1 Apolipoprotein B	red green
h-cat.k / cat.λ Code 93115	Human Ig light chain k Human Ig light chain λ	green blue

PRINCIPLE OF THE METHOD

EASY RID MULTIPLATE allows the determination of human plasma proteins in radial immunodiffusion. The antigen (protein), shown in the well of the plate, diffuses radially in the agarose gel reacting with specific antibodies incorporated in agarose gel and forming immune complexes visible as precipitin rings. Diameter of precipitin ring is directly proportional to the concentration of the relevant protein in the sample.

COMPOSITION

EASY RID MULTIPLATE contains, in different agarose gels, monospecific antisera produced immunizing rabbits or goats.

TAKING AND STORAGE OF SAMPLES

Human serum or plasma samples, fresh or stored correctly at 2-8 °C, can be used. Turbid samples must be clarified by centrifugation prior to the assay.

Fresh samples can be deep-frozen immediately after collection and used afterwards.

PREPARATION OF SAMPLES

According to the kind of plasma protein you are looking for, use a diluted sample in PBS or in isotonic salty solution, or an undiluted sample as indicated in table n°2.

Table n°2

PRODUCT	PLASMAPROTEIN	SAMPLE
h-IgG / IgA / IgM	IgG IgA IgM	Undiluted Undiluted Undiluted
h-C3c / C4	C3c C4	Undiluted Undiluted
h-ApoA1 / ApoB	Apolipoprotein A1 Apolipoprotein B	Undiluted Undiluted
h-cat.k / cat.λ	Human Ig light chain k Human Ig light chain λ	1:6 (1+5) 1:6 (1+5)

TEST PROCEDURE

1. Remove EASY RID MULTIPLATE from the envelope, open the plate and leave it for about 5' at room temperature so that any condensed water in the wells can evaporate.
2. Fill the well relevant to the plasma protein to search with 5 µL of sample as indicated in table n°2.
3. After the samples have diffused into the gel for about 20', close the plate with the lid and leave to stand, overturned, in the envelope, at room temperature for 48 hours.

INTERPRETATION OF THE RESULTS

After 48 h, measure the precipitin rings diameters using a suitable measuring device (mm).

The protein concentration for the precipitin ring diameters can be read using the attached tables.

CLINICAL INTERPRETATION

Refer to the data of bibliography.

QUALITY CONTROL

Quality Control is carried out performing a control of accuracy and precision inoculating a control serum on the plate in examination.

The expected results for the described method are: Precision (CV) ≤ 5%; Accuracy = ± 15% reference value.

LIMITS AND WARNINGS

The inoculation of sample with a protein concentration inferior to the minimum measurable value leads to a ring formation difficult to measure; in this case it is recommended to use more sensitive reading tools.

Otherwise, the inoculation of sample with a protein concentration higher than the maximum reference value brings to a non measurable ring formation; in such case it is recommended to dilute the sample.

PRECAUTIONS

All EASY RID MULTIPLATE configurations contain sodium azide. When diagnostics *in vitro* containing sodium azide are used, observe the following precautions: do not swallow and avoid contact with skin or mucous membranes. Sodium azide can form explosive azides when it is in contact with heavy metals like copper or lead.

All EASY RID MULTIPLATE configurations contain sodium azide with a concentration < 1g/L; they are classified as not dangerous according to current law therefore safety data sheet is not necessary.

STORAGE

Store EASY RID MULTIPLATE at 2-8 °C in their original packagings. Avoid to store near heat sources and avoid change of temperature. In these conditions EASY RID MULTIPLATE is valid until expire date indicated on the label. Do not use over this date.

A plate partially used can be used for 3 weeks starting from the package opening. After each use put plate overturned into the envelope and store at 2-8 °C.

Eliminate if there are deteriorations before use (dehydration or media separation from support, presence of bacteria or fungi). Do not freeze.

ELIMINATION OF THE USED MATERIAL

After use, the EASY RID MULTIPLATE and material come in contact with samples must be decontaminated and potentially digested in accord with the techniques in use in the laboratory for decontamination and disposal of infected material.

PRESENTATION

PRODUCT	CODE	TEST
h-IgG / IgA / IgM	93104	4
h-C3c / C4	93106	6
h-ApoA1 / ApoB	93110	6
h-cat.k / cat.λ	93115	6

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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