Code 93201

20 Test

BENCE JONES TEST
Piastra a doppia immunodiffusione
Double immunodiffusion plate

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COMPOSITION:

The BENCE JONES plate consists of a layer of agarose gel. The plate contains sodium azide (1 g/l) as preservative.

Warnings and Precautions

For In Vitro Diagnostic Use only.

Reagents containing sodium azide must be handled with due caution. Do not ingest or allow to contact with skin or mucous membranes.

STORAGE AND STABILITY:

Stored at +2°/+8°C in the original unopened containers, BENCE JONES plates are stable up to labelled expiry date.

PRINCIPLE OF THE METHOD:

Determination of K and λ chains by double immunodiffusion method.

The protein under analysis (antigen) forms immune complexes with specific antibody in the agarose gel of the plate. These complexes become visible as precipitation line.

PROCEDURE:

Remove BENCE JONES plate from the envelope, open the plate and before loading samples, leave to stand for about five minutes at room temperature (RT) so that any condensed water in the wells can evaporate.

Fill the central well with 40 μ l of patient sample (urine), the right well with 10 μ l of anti-K chains, the left well with 10 μ l of anti- λ chains.

Leave the plate to stand for about 20 min., until the samples are adsorbed in the agarose. Close the lid, overturn the plate into the envelope and allow to complete the diffusion for 9-15 hours at 37°C.

Results evaluation:

The lack of immunoprecipitation lines excludes presence of K and λ chains in urine. The presence of precipitation lines indicates presence of K and λ chains in urine.

It is possible to roughly evaluate the quantity of protein in the sample by the position of the precipitation lines. If it is near the antiserum well, sample's concentration is high. If it is near the sample well, sample's concentration is low.

The detection limit is about 2,8 mg/dl for K chains and 2,1 mg/dl for λ chains.

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CONTENT OF KIT

- 20 Plates with agarose gel
- 1 Tube of Anti h-Ig/ K chains x 400 μ l
- 1 Tube of Anti h-Ig/ λ chains x 400 μ l
- 1 Tube of positive Control K and λ x 100 μ l

Bibliografia/ Bibliografy:

1. Bulletin of the World Health Organization, 59,717, 1981
2. Sommer R., Hohenwalner W. Laboratoriumsblätter 29, 27, 1979
3. Dammacco F. Immunologia in medicina; Edi Ermes, 1989.

Rev.2/02.02.2010