# RID CONTROL SERUM

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# RID CONTROL SERUM

## COMPOSITION

The serum (human) used for quality control in **RID CONTROL SERUM** is constituted by a pool of healthy adult's sera. The preparation is liquid and it is stabilized by removing weak components (lipids).

## USE OF REAGENT

The vial's content is supplied ready for use. Bring at room temperature, then invert gently to mix, avoiding vigorous shaking (foam formation).

**RID CONTROL SERUM** has to be used according to the procedure and the dilutions requested in radial immunodiffusion.

## WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic use only.
2. Reagents containing sodium azide must be handled with due caution: do not ingest or allow to contact with skin or mucous membranes !
4. Each individual blood donation used to manufacture of rid control serum, is tested for hepatitis B surface antigen (HbsAg), anti-HCV, anti-HIV1 and anti-HIV2, according to FDA requirements. Only donations with negative findings are used for manufacture.

Nevertheless, as complete absence of infectious agents can never be assured, all materials derived from human blood should be handled with due care following the precautions recommended for bio hazardous material.

## STORAGE AND STABILITY

Store the product at +2/+8°C : in this conditions it remains stable up to the labelled expiry date indicated on the label. Once opened, vials can be used within a week. Avoid contamination.

## BIBLIOGRAPHY

- Whicher J T, Ritchie R, Johnson M A, Baudner S, Bienvenu J, Blirup-Jensen S, et al: New international reference preparation for proteins in human serum. 1994 *Clin Chem*; **40**:934-938.
- Baudner S, Haupt H, Hubner R.: Manufacture and characterization of a new reference preparation for 14 plasmaproteins. 1994. *Clin Chem*; 8:177-190

## PRESENTATION

PRODUCT	CODE	PACKAGE
RID CONTROL SERUM	940010	4 vials x 0.5 ml

## TABLE OF SYMBOLS

SYMBOL	MEANINGS	SYMBOL	MEANINGS
	In Vitro Diagnostic Medical Device		Do not reuse
	Manufacturer		Contains sufficient for <n> tests
	Catalogue number		Fragile, handle with care
	Use by		Caution, consult accompanying documents
	Temperature limitation		Batch code

Rev.5/31.03.2005

**LIOFILCHEM Bacteriology Products**

64026 ROSETO D.A. (TE) ITALY- Via Scozia- Zona Ind.le

Tel.+39 085 8930745 - Fax +39 085 8930330

Sito Web: <http://www.liofilchem.net> E-Mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)