DESCRIZIONE

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EASY RID

DESCRIPTION

EASY RID are plates with 12 wells for quantitative determination in radial immunodiffusion of human plasmaproteins in serum and plasma. EASY RID made by Liofilchem S.r.l. is available in 17 configurations.

CONTENT OF THE KIT

Each kit contains:- 1 plate - 1 instructions sheet.

MATERIAL REQUIRED BUT NOT PROVIDED:- Precision pipette - Suitable measuring device (mm)

CONFIGURATION

EASY RID proposed by Liofilchem S.r.l. are available in the following configurations:

Table no.1

PRODUCT	CODE	PLASMAPROTEIN	GEL COLOUR	PRODUCT	CODE	PLASMAPROTEIN	GEL COLOUR
h-IgG	93001	IgG	green	h-Alfa 1 Glicoproteina acida	93010	Alpha 1 Acid Glycoprotein	red
h-IgA	93002	IgA	red	h-Fibrinogeno	93011	Fibrinogen	violet
h-IgM	93003	IgM	blue	h-Antitrombina III	93012	Antithrombin III	yellow
h-C3c	93004	C3c	yellow	h-Ig/Catena kappa	93013	Human Ig light chain K	green
h-C4	93005	C4	red	h-Ig/Catena Lambda	93014	Human Ig light chain λ	blue
h-Transferrina	93006	Transferrin	green	h-Alfa 1 Antitripsina	93015	Alpha 1 Antitrypsin	yellow
h-Albumina	93007	Albumin	orange	h-Ceruloplasmina	93016	Ceruloplasmin	red
h-Apolipoproteina A1	93008	Apolipoprotein A1	red	h-Aptoglobina	93018	Haptoglobin	green
h-Apolipoproteina B	93009	Apolipoprotein B	green				

PRINCIPLE OF THE METHOD

EASY RID allows the determination of human plasmaproteins in radial immunodiffusion. The antigen (protein) inoculated in the well of the plate, diffuses radially in the agarose gel, reacting with specific antibodies incorporated in agarose gel and forming immune complexes visible as precipitin rings. Diameter of precipitin ring is directly proportional to the concentration of the relevant protein in the sample.

COMPOSITION

EASY RID contains, in a layer of agarose gel the monospecific antiserum to human plasmaprotein, produced by immunizing rabbits or goats.

PICKING AND STORAGE OF SAMPLES

It can be used human serum or plasma samples, fresh or stored correctly at 2-8 °C. Use plasma for Fibrinogen and Antithrombin III and serum for the other proteins.

Turbid samples must be clarified by centrifugation prior to the assay.

Fresh samples can be deep-frozen immediately after collection and used after.

PREPARATION OF SAMPLES

According to the kind of plasmaprotein you are looking for, use a diluted sample in PBS or in isotonic salty solution, or an undiluted sample as indicated in table no.2.

Table no.2

PRODUCT	SAMPLE	PRODUCT	SAMPLE	PRODUCT	SAMPLE	PRODUCT	SAMPLE
h-IgG	(Undiluted)	h-C4	(Undiluted)	h-Apolipoproteina B	(Undiluted)	h-Ig/Catena kappa	1:6 (1+5)
h-IgA	(Undiluted)	h-Transferrina	(Undiluted)	h-Alfa 1 Glicoproteina acida	1:2 (1+1)	h-Ig/Catena Lambda	1:6 (1+5)
h-IgM	(Undiluted)	h-Albumina	1:11 (1+10)	h-Fibrinogeno	(Undiluted)	h-Alfa 1 Antitripsina	(Undiluted)
h-C3c	(Undiluted)	h-Apolipoproteina A1	(Undiluted)	h-Antitrombina III	(Undiluted)	h-Ceruloplasmina	(Undiluted)
						h-Aptoglobina	(Undiluted)

TEST PROCEDURE

1. Remove EASY RID from the envelope, open the plate and leave to stand for about 5 minutes at room temperature so that any condensed water in the wells can evaporate.

2. Fill the wells with 5 µl of undiluted or diluted patient samples as indicated in table no.2.

3. Close the plate with the lid, after the samples have diffused into the gel for about 20', leave to stand, overturned into the envelope, at room temperature for 48 hours.

INTERPRETATION OF THE RESULTS

After 48h, measure the diameters of the precipitin rings using a suitable measuring device (mm). The protein concentration for the precipitin ring diameters can be read using the attached tables.

CLINICAL INTERPRETATION

Make reference to the data of bibliography.

QUALITY CONTROL

Quality Control is realized performing a control of accuracy and precision sowing on the plate in examination a serum of control.

The results attended with the described method are: Precision CV ≤ 5% Accuracy = ± 15% reference value.

LIMITS AND WARNINGS

The filling of sample with protein concentration inferior to the value least measurable lead to the formation of a ring difficult to measure; in this case it is recommended to use tools of reading with greater sensibility.

In the contrary case, the filling of sample with protein concentration superior to the maximum value of reference brings to the formation of a non measurable ring; in such case it is recommended to dilute the sample.

PRECAUTIONS

All EASY RID configurations contain sodium azide. When are employed diagnostics *in vitro* containing sodium azide, observe the following precautions: do not swallow and avoid contact with skin or mucous membranes. Sodium azide can form explosive azides when in contact with heavy metals like copper or lead. All EASY RID configurations contain sodium azide with a concentration < 1g/l; they are classified as not dangerous according to in force law so availability of safety data sheet is not necessary.

STORAGE

Store EASY RID at 2-8 °C in their original packagings. Avoid to store near heat sources and avoid change of temperature. In these conditions EASY RID is valid until expire date indicated on the label. Don't use over this date. A plate partially used can be used for 3 weeks starting from opening of the packaging. After each use put plate overturned into the envelope and store at 2-8 °C. Eliminate if there are deteriorations before use (great change of media colour, dehydration or media separation from support, presence of bacteria or fungi). Do not freeze.

ELIMINATION OF THE USED MATERIAL

After use, the EASY RID and material come in contact with samples must be decontaminated and potentially digested in accord with the techniques in use in the laboratory for decontamination and disposal of infected material.

BIBLIOGRAPHY

- Bulletin of the World Health Organization. 1981. **59**: 717.
- Dammacco, F. 1989. *Immunologia in Medicina*; Edi Ermes.
- Fundenberg, H.H., et al. *Basic & Clinical Immunology*, Lange Medical Pub.,California (1980).
- Marcovina S., G.Di Cola and C.Rapetto. 1985. *Cli.Chim.Acta* **147**:117-25
- Marcovina S., G.Di Cola and A.L.Catapano.1986.*Clin.Chem.* **32**:12, 2155-2159.
- Parker, C.W. ed.,*Clinical Immunology*, Saunders Company, Philadelphia, Vol.I, II (1980).
- Putnam, F.W.ed., *The plasma proteins*, Academic Press, New York, Vol. 1,2 (1975), Vol. 3 (1977).
- Schultze, H.E. and J.F.Heremans. 1966. *Molecular Biology of Human proteins*;Elsevier

PRESENTATION

PRODUCT	CODE	TEST	PRODUCT	CODE	TEST	PRODUCT	CODE	TEST	PRODUCT	CODE	TEST
h-IgG	93001	12	h-C4	93005	12	h-Apolipoproteina B	93009	12	h-Ig/Catena kappa	93013	12
h-IgA	93002	12	h-Transferrina	93006	12	h-Alfa 1 Glicoproteina acida	93010	12	h-Ig/Catena Lambda	93014	12
h-IgM	93003	12	h-Albumina	93007	12	h-Fibrinogeno	93011	12	h-Alfa1 Antitripsina	93015	12
h-C3c	93004	12	h-Apolipoproteina A1	93008	12	h-Antitrombina III	93012	12	h-Ceruloplasmina	93016	12
									h-Aptoglobina	93018	12

TABLE OF SYMBOLS

 Manufacturer	 Content sufficient for <n> tests	 Temperature limitations	 Use by	 In Vitro Diagnostic Medical Device
 Catalogue number	 Fragile, handle with care	 Caution, consult accompanying documents	 Batch code	 Do not reuse

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