OXI CONTROL GST E6/E5/E4/E3/E2 SS

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F12611
Rev.2 / 10.07.2019

OXI CONTROL GST E6/E5/E4/E3/E2 SS

Biological indicators of H₂O₂ vapour superficial sterilization processes containing *Geobacillus stearothermophilus* (ATCC 7953) spores inoculated on stainless-steel coupons.

DESCRIPTION

USP (United States Pharmacopoeia), EP (European Pharmacopoeia) and DAB (Deutsches Arzneibuch) standards recommend to use bioindicators during sterilization processes. Biological indicators **OXI CONTROL GST E6/E5/E4/E3/E2 SS** are used for regular control of H₂O₂ vapour superficial sterilization cycles and validation of sterilizers. These biological indicators use a stainless-steel support (coupon), that is not permeable to hydrogen peroxide. Stainless-steel coupon has been inoculated with *Geobacillus stearothermophilus* (ATCC 7953) spores and has been introduced into a Tyvek®/Polytene envelope. Tyvek® side of the envelope is permeable to hydrogen peroxide vapours, while the Polytene side of the envelope is not permeable to these vapours. A small circle on the coupon visibly ensures that it has been inoculated and that is correctly oriented in the envelope. Biological indicators **OXI CONTROL GST E6/E5/E4/E3/E2** are produced under strictly controlled conditions in order to satisfy the requirements indicated in the USP and EP current editions and in accordance with ISO 11138 and EN 866 standards.

PRINCIPLE

OXI CONTROL GST E6/E5/E4/E3/E2 contain *Geobacillus stearothermophilus* (ATCC 7953) spores in predefined concentrations: E6 = 1-5 x 10⁶ CFU/coupon; E5 = 1-5 x 10⁵ CFU/coupon; E4 = 1-5 x 10⁴ CFU/coupon, E3 = 1-5 x 10³ CFC/coupon, E2 = 1-5 x 10² CFU/coupon. Spores are completely killed if the sterilization cycle has been efficient. In this case, during the following incubation in the Steri-Test Medium tubes, included in the package, spores are not able to grow and to modify medium's aspect. On the contrary, in the case the sterilization process has not been efficient, spores partially survive and therefore during the following incubation will be able to grow and to modify medium's aspect.

INSTRUCTIONS FOR USE

Before using biological indicators **OXI CONTROL GST E6/E5/E4/E3/E2 SS**, allow them to reach room temperature (about 1 hour). Coupons are exposed to vapours, during sterilization, inside their original envelope or aseptically removed from their envelope and hanged utilizing the special hole.

Exposure of indicators

Put envelopes containing coupons or coupons without their envelopes inside the sterilization environment in the points considered the most difficult to reach and to sterilize by hydrogen peroxide vapour.

Follow the operative instructions indicated by sterilizer manufacturer for the sterilization process.

NOTE 1 (exposure inside the envelope): do not cover the Tyvek® side (white side) of the biological indicator.

NOTE 2 (exposure without envelope): during transfer of coupons without envelope it is essential to operate in asepsis to avoid accidental contaminations which could invalidate results.

Coupons culture

At the end of sterilization/aeration process, remove the biological indicators, with or without envelope, and deliver them to the laboratory for the culture together with a not exposed control (positive control), belonging to the same batch

In optimal conditions all culture procedures have to be performed under a laminar flow cabinet following aseptic methods.

Submit to culture the biological indicators, opening the Tyvek®/Polytene envelope with a cooled sterile scissors and carefully transferring both coupons, the exposed ones with envelope and those exposed without envelope, inside a Steri-Test Medium tube, taking them by cooled sterile tweezers and letting them fall into the medium. Ensure that the coupons are completely dipped in the medium gently shaking the tube.

Incubation

Incubate the Steri-Test Medium tubes, that contain coupons, at 55-60°C (131-140°F) for 7 days or for a shorter time validated by user.

EVALUATION

Geobacillus stearothermophilus (ATCC 7953) spores are killed off if the sterilization cycle has been efficient: in this case the medium's colour remains violet/clear even after incubation at 55-60°C (131-140°F) for the selected time.

If the sterilization cycle has not been efficient, spores partially survive and medium's colour turns yellow/turbid after incubation at 55-60°C (131-140°F) for the selected time.

The tube inoculated with the coupon not submitted to the sterilization cycle and used as spore growth control (positive control), has to turn yellow/turbid after incubation. On the contrary, the test must be repeated after having investigated the causes of the negative result.

EVALUATION TABLE		
MEDIUM COLOUR	SPORE	STERILIZATION
Violet / Clear	Killed off	Successful
Yellow / Turbid	Vital	Unsuccessful

TREATMENT OF STERILTEST MEDIUM TEST TUBES AFTER USE

After use, sterilize the positive tubes (yellow/turbid) in autoclave at 121°C for at least 30 minutes and eliminate them in accordance with the procedures of the laboratory.

STORAGE

Store the product at 2-8°C: in these conditions it maintains its validity until the expiry date indicated on the label.

BIBLIOGRAPHY

1. United States Pharmacopoeia latest edition.
2. Deutsches Arzneibuch latest edition.
3. European Pharmacopoeia latest edition.
4. ISO 11138 and EN 866 latest edition.

PRESENTATION

PRODUCT	REF	PACKAGING	Spore: UFC / coupon	D _{VHP} (1,6 ± 0,5 mg/l, 30 °C)
OXI CONTROL GST E6 SS	91082	20 coupon+20 Steri-Test Medium tubes	1-5x10 ⁶	1,5-4,0 minutes
OXI CONTROL GST E6 SS	91182	100 coupon	1-5x10 ⁶	1,5-4,0 minutes
OXI CONTROL GST E5 SS	91083	20 coupon+20 Steri-Test Medium tubes	1-5x10 ⁵	1,5-4,0 minutes
OXI CONTROL GST E4 SS	91081	20 coupon+20 Steri-Test Medium tubes	1-5x10 ⁴	1,5-4,0 minutes
OXI CONTROL GST E4 SS	91181	100 coupon	1-5x10 ⁴	1,5-4,0 minutes
OXI CONTROL GST E3 SS	91080	20 coupon+20 Steri-Test Medium tubes	1-5x10 ³	1,5-4,0 minutes
OXI CONTROL GST E3 SS	91180	100 coupon	1-5x10 ³	1,5-4,0 minutes
OXI CONTROL GST E2 SS	91079	20 coupon+20 Steri-Test Medium tubes	1-5x10 ²	1,5-4,0 minutes

TABLE OF SYMBOLS

 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Caution, consult accompanying documents
 Use by	 Batch code	 Do not reuse

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 www.liofilchem.com

F12611
Rev.2 / 10.07.2019