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Biological indicators of dry heat or ethylene oxide sterilization processes containing *Bacillus atropheus* (ATCC 9372) spores in suspension

DESCRIPTION

USP (United States Pharmacopoeia), EP (European Pharmacopoeia) and DAB (Deutsches Arzneibuch) standards recommend to use bioindicators during dry heat and ethylene oxide sterilization processes.

Biological indicators **SUSPENSION SPORE BAT E10/E9/E8/E7/E6/E5** are used for validation, re-validation and monitoring process of dry heat and ethylene oxide sterilizers.

Biological indicators **SUSPENSION SPORE BAT E10/E9/E8/E7/E6/E5** are produced under strictly controlled conditions in order to satisfy the requirements indicated in the USP current edition and in accordance with ISO 11138 and EN 866 standards.

PRINCIPLE

SUSPENSION SPORE BAT E10/E9/E8/E7/E6/E5 contain *Bacillus atropheus* (ATCC 9372) spores in distilled water in predefined concentrations: E10=1-5x10¹⁰CFU/ml; E9=1-5x10⁹CFU/ml; E8=1-5x10⁸CFU/ml; E7=1-5x10⁷CFU/ml; E6=1-5x10⁶CFU/ml; E5=1-5x10⁵CFU/ml.

Spores are completely killed off during dry heat or ethylene oxide sterilization processes, if the process has been efficient. In this case the suitable medium's aspect, inoculated with spores and incubated at a suitable time, will remain unchanged.

On the contrary if the sterilization process has not been efficient, spores partially survive and medium's aspect will change.

INSTRUCTIONS FOR USE

- Take a suitable volume of the suspension and put it on the selected material. Protect appropriately the inoculated material in order to avoid spore dispersion.
- Put the inoculated and protected material on the bottom, in the centre, inner sides and on the critical points of the sterilizer.
- For sterilizers with capacity up to 250 litres put two units for each selected point of the sterilizer. For sterilizers with capacity higher than 250 litres put six or more units for each selected point.
- Remove the inoculated material from the sterilizer after the sterilization cycle.
- Transfer aseptically the material in a container with a culture medium suitable for microorganism growth as for example Steri-Test Medium (Ref. 20199).
- Incubate the containers with the medium and the material at 30-35°C (86-95°F) for 7 days or for a shorter time validated by user.
- Incubate, at the same conditions of time and temperature, an unit of material not submitted to the sterilization cycle, inoculated with a spore suspension belonging to the same batch, as spore growth control (positive control).
- Examine medium's aspect in the container and interpret results as per **EVALUATION TABLE**: a change of medium's aspect indicates microbial growth and therefore an unsuccessful sterilization. On the contrary, the persistence of the medium's initial aspect indicates absence of microbial growth and therefore a successful sterilization.

EVALUATION

Bacillus atropheus (ATCC 9372) spores are killed off if the sterilization cycle has been efficient: in this case the medium remains unchanged even after incubation at 30-35°C (86-95°F) for the selected time. If the sterilization cycle has not been efficient, spores partially survive and the medium's aspect change after incubation at 30-35°C (86-95°F) for the selected time. The medium in the container inoculated with the material not submitted to the sterilization cycle and used as spore growth control, has to change aspect after incubation. On the contrary, the test must be repeated after having investigated the causes of the negative result.

EVALUATION TABLE		
APPAREANCE OF THE MEDIUM	SPORE	STERILIZATION
Unchanged	Killed off	Successful
Changed	Vital	Unsuccessful

MATERIAL TREATMENT AFTER USE

After use, sterilize the positive containers (changed medium's aspect) in autoclave at 121°C for at least 30 minutes and eliminate them in accordance with the procedures of the laboratory.

STORAGE

Store the product at 2-8°C: in these conditions it maintains its validity until the expiry date indicated on the label.

BIBLIOGRAPHY


- United States Pharmacopoeia last edition.
- Deutsches Arzneibuch last edition.
- European Pharmacopoeia last edition.
- ISO 11138 ed EN 866 last edition.

PRESENTATION

PRODUCT	REF	PACKAGING	Spores: CFU/ml	D _{EO} (600 ± 30 mg/l, 60% ± 10% RH, 54 ± 1°C)	D _{DH} (160 ± 1°C)
SUSPENSION SPORE BAT E10	82006	10 ml vial	1-5x10 ¹⁰	2.6-4.5 minutes	1.0-3.0 minutes
SUSPENSION SPORE BAT E9	82007	10 ml vial	1-5x10 ⁹	2.6-4.5 minutes	1.0-3.0 minutes
SUSPENSION SPORE BAT E8	82008	10 ml vial	1-5x10 ⁸	2.6-4.5 minutes	1.0-3.0 minutes
SUSPENSION SPORE BAT E7	82009	10 ml vial	1-5x10 ⁷	2.6-4.5 minutes	1.0-3.0 minutes
SUSPENSION SPORE BAT E7	82010	10 ml vial	1-5x10 ⁶	2.6-4.5 minutes	1.0-3.0 minutes
SUSPENSION SPORE BAT E5	82011	10 ml vial	1-5x10 ⁵	2.6-4.5 minutes	1.0-3.0 minutes

D value for ethylene oxide (DEO) is calculated with the MPN and SC methods with 600±30 mg/l of ethylene oxide, at a temperature of 54±1°C, with 50-70% humidity. D value for dry heat (DDH) is calculated with MPN method at a temperature of 160±1°C.

TABLE OF SYMBOLS

LOT	Batch code	 Do not reuse	 Manufacturer	 Use by	 Fragile, handle with care
REF	Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult accompanying documents	

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LIOFILCHEM® S.r.l.