

## Steam Control GST E6 SS (coupons)

Biological indicators of steam sterilization processes with spores of *Geobacillus stearothermophilus* (ATCC 7953) inoculated on stainless steel coupons.

### DESCRIPTION

Steam Control GST E6 SS are used for regular control of steam sterilization cycles and validation of sterilizers. These biological indicators use a stainless steel support (coupon) previously inoculated with *Geobacillus stearothermophilus* (ATCC 7953) spores and enclosed in a steam permeable envelope.

The coupons are produced under strictly controlled conditions in order to satisfy the requirements indicated in the current editions of USP and EP and in accordance with ISO 11138 and EN 866 standards.

### COMPOSITION

The coupons are inoculated with spores of *Geobacillus stearothermophilus* (ATCC 7953) in a predefined concentration: E6 =  $1.5 \times 10^6$  CFU/coupon.

A certificate of analysis is included in the package providing the following data:

- Certified population
- Species
- D-value (121°C)
- Z-value
- Survival time
- Kill time
- Lot number
- Expiration date

### PRINCIPLE

The spore containing coupons are exposed to the sterilization cycle and then incubated in a culture medium to promote germination (either Liofilchem® Steri Test Medium or Tryptic Soy Broth can be used for spores germination). Sterilization is achieved if no spores survive: there will be any growth after incubation with no change of the medium. If the sterilization process is not performed correctly, spores will grow during the incubation period and cause turbidity of the medium or even a colour change if Steri Test Medium is used.

### TECHNIQUE

1. Allow the coupons to reach room temperature before use. Coupons are exposed to steam, during sterilization, inside their original envelope or aseptically removed from their envelope and hanged utilizing the special hole.
2. Put the coupons, with or without the envelopes, inside the autoclave or plant steam treated with SIP, in the points considered the most difficult to sterilize by steam.
3. Follow the operative instructions indicated by the sterilizer manufacturer for the sterilization process. **NOTE:** when coupons are moved without envelope, it is essential to operate in asepsis to avoid accidental contaminations that could invalidate the results.
4. At the end of sterilization cycle, remove the coupons and deliver them with the unexposed control (positive control) from the same batch to the laboratory for culturing. In optimal conditions all culture procedures have to be performed under a laminar flow cabinet following aseptic methods.
5. Cut the envelope with a cold sterile scissors and by means of a cold sterile tweezers carefully transfer the coupon to a tube of Steri Test Medium or Tryptic Soy Broth (TSB). Gently shake the tube to ensure the coupon is completely dipped in the medium.
6. Incubate the tubes containing the coupons at 55-60°C (131-140°F) for at least 7 days.

### INTERPRETATION OF RESULTS

*Geobacillus stearothermophilus* (ATCC 7953) spores are killed off if the sterilization cycle is efficient: no changes in the medium after the incubation period. If the sterilization cycle is not performed correctly, spores partially survive: the Steri Test Medium changes from violet/clear to yellow/turbid after incubation; in TSB spores survival is indicated by turbidity only. The medium inoculated with the unexposed coupon (positive control) should demonstrate microbial growth after incubation. If not, the test must be repeated after investigating the causes of such result.

EVALUATION TABLE		
MEDIUM CHARACTERISTICS	SPORE	STERILIZATION
Violet (Steri Test Medium) / Clear	Killed off	Successful
Yellow (Steri Test Medium) / Turbid	Vital	Unsuccessful

### STORAGE

Store at 2-8°C: in these conditions the product maintains its validity until the expiry date indicated on the label. **NOTE:** The product may be held at room temperature (25°C) for up to 72 hours prior to use without any impact to the performance.

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for professional use only and must be used by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.



**LIOFILCHEM® S.r.l.**

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY  
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.com

#### REFERENCES

1. United States Pharmacopoeia current edition.
2. Deutsches Arzneibuch current edition.
3. European Pharmacopoeia current edition.
4. ISO 11138 and EN 866 current edition.



**LIOFILCHEM® S.r.l.**

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY  
Tel +39 0858930745 Fax +39 0858930330 Website: [www.liofilchem.com](http://www.liofilchem.com)

## PRODUCT SPECIFICATIONS

### NAME

Steam Control GST E6 SS (coupons)

### PRESENTATION

Coupons inoculated with a predefined concentration of spores: *G. stearothermophilus* ATCC 7953

### STORAGE

2-8°C

### PACKAGING

Ref.	Content	Packaging
91097	50 coupons	Coupons individually packed

### TECHNICAL PROPERTIES

**Primary packaging:** Steam-permeable envelope

**Spore carrier:** AISI 316L stainless steel coupon (approximately 20 x 7 mm)

**Species:** *Geobacillus stearothermophilus* ATCC 7953

**Mean Population Recovery:**  $1 \times 10^6$ - $5 \times 10^6$  spores per coupon

**Purity:** Bacterial contaminates less than 1 percent of the labeled population; this is the detection limit using the pour plate method with an aliquot that yields at least 100 colony forming units (CFU)

**Resistance data:** D-value, Z-value, survival time and kill time

### USE

Steam Control GST E6 SS (coupons) are biological indicators used for regular control of steam sterilization cycles and validation of sterilizers

### TECHNIQUE

Refer to technical sheet of the product

### APPEARANCE

Coupons are grey in colour










### SHELFLIFE

3 years

### QUALITY CONTROL

1. Control of general characteristics, label and print
2. Purity: <1% contamination. No moulds
3. Heat shocked population:  $1-5 \times 10^6$  Spores/coupons
4.  $D_{121^\circ\text{C}}$  : 1.5-4.5 minutes
5. Z-value  $\geq 6^\circ\text{C}$

### TABLE OF SYMBOLS

 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Caution, consult accompanying documents
 Use by	 Batch code	 Do not reuse



**LIOFILCHEM® S.r.l.**

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY  
 Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.com