Vaseline Oil Droppers

Vaseline oil for microbiological use, in dropper bottles.

INTENDED PURPOSE

Vaseline Oil for microbiological use, used to keep an anaerobic environment for the microbial growth and for the performance of some biochemical tests suitable for microbial identification. This reagent is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Vaseline Oil Droppers is a material used to keep an anaerobic environment for the microbial growth and for the performance of some biochemical tests suitable for microbial identification. The product use is recommended specially for Liofilchem microbial identification systems, when specified in test procedure of the system.

KIT CONTENT

- 10 bottles of Vaseline Oil Droppers (10.0 ml)

METHOD PRINCIPLE

Vaseline Oil Droppers, due to its physical-chemical properties, prevents environmental oxygen from coming into contact with medium and microorganism, assuring therefore the keeping of an anaerobiosis level sufficient to the aim the product is set for.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: culture media, quality control organisms.

REAGENTS

- Vaseline Oil

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

For the product use refer to the instruction sheet of the relevant identification systems and to the laboratory procedures.

INTERPRETING RESULTS

For the interpretation of results refer to the instruction sheet of the relevant identification systems and to the laboratory procedures.

STORAGE

10-25°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

3 years.

QUALITY CONTROL

Control of the appearance: limpid colourless solution.

WARNING AND PRECAUTIONS

- 1) For *in vitro* diagnostic use (IVD).
- 2) For laboratory professional use only.

- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Vaseline Oil	10 x 10 ml	87006

Revision History

Revision	Release Date	Change Summary
3	2024-09-16	Added notice to report any malfunction, defect or incident.
2	2023-05-15	Updated layout and content in compliance with IVDR 2017/746

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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