Oxi Control GST E6/E5/E4 PLX

Geobacillus stearothermophilus ATCC® 7953 spores inoculated on poliflex coupons for vaporized hydrogen peroxide sterilization control.

DESCRIPTION

USP (United States Pharmacopoeia), EP (European Pharmacopoeia) and DAB (Deutsches Arzneibuch) standards recommend to use biological indicators (BIs) during sterilization processes. BIs Oxi Control GST E6/E5/E4 PLX are used for regular control of H₂O₂ vapour superficial sterilization cycles and validation of sterilizers. These BIs use a poliflex support (coupon), that is not permeable to hydrogen peroxide. Poliflex coupon is inoculated with *Geobacillus stearothermophilus* ATCC® 7953 spores and introduced into a Tyvek®/Polytene envelope. The Tyvek® side of the envelope is permeable to hydrogen peroxide vapours, while the Polytene side is not permeable to these vapours. A small circle on the coupon visibly ensures that it has been inoculated and that is correctly oriented in the envelope.

Oxi Control GST E6/E5/E4 PLX are produced under strictly controlled conditions in order to meet the current USP/EP requirements and in accordance with ISO 11138 and EN 866 standards.

PRINCIPLE

Oxi Control GST E6/E5/E4 PLX contain *Geobacillus stearothermophilus* ATCC® 7953 spores in predefined concentrations: 1-5 x 10⁶ (E6), 1-5 x 10⁵ (E5), and 1-5 x 10⁴ (E4) CFU per coupon. Spores are completely killed if the sterilization cycle has been efficient. In that case, during the following incubation in the Steri-Test Medium, which is included in the package, spores are not able to grow and the medium remains violet/clear. On the contrary, in case of an inefficient sterilization, spores partially survive and will grow causing a change of the medium to yellow/turbid.

INSTRUCTIONS FOR USE

Before using Oxi Control GST E6/E5/E4 PLX, allow them to reach room temperature (about 1 hour). Coupons are exposed to vapours, during sterilization, inside their original envelope.

Exposure of indicators

Put envelopes containing coupons inside the sterilization environment in the points considered the most difficult to reach and to sterilize by hydrogen peroxide vapour. Follow the operative instructions indicated by sterilizer manufacturer for the sterilization process.

NOTE: Do NOT cover the Tyvek® side (white side) of the BI.

Coupons culture

At the end of sterilization/aeration process, remove the BIs and deliver them to the laboratory for the culture together with a not exposed control (positive control), belonging to the same batch.

In optimal conditions all culture procedures have to be performed under a laminar flow cabinet following aseptic techniques to prevent contamination. Open the Tyvek®/Polytene envelope with sterile scissors and carefully transfer the coupons into a Steri-Test Medium tube by means of sterile tweezers. Ensure that the coupons are completely dipped in the medium by gently shaking the tube.

Incubation

Incubate the Steri-Test Medium tubes, with the coupons, at 55-60°C (131-140°F) for 7 days or for a shorter time validated by the user.

EVALUATION

Geobacillus stearothermophilus ATCC® 7953 spores are killed off if the sterilization cycle has been efficient: After incubation at 55-60°C (131-140°F) for the selected time, the medium remains violet/clear.

If the sterilization cycle has not been efficient, the spores partially survive and the medium turns yellow/turbid after incubation at 55-60°C (131-140°F) for the selected time.

The tube inoculated with the coupon but not submitted to the sterilization cycle (positive control), has to turn yellow/turbid after incubation. On the contrary, the test must be repeated after having investigated the causes of the negative result.

EVALUATION TABLE

MEDIUM COLOUR	SPORE	STERILIZATION
Violet / clear	Killed off	Successful
Yellow / turbid	Vital	Unsuccessful

TREATMENT OF STERILTEST MEDIUM TEST TUBES AFTER USE

After use, sterilize the positive tubes (yellow/turbid) in autoclave at 121°C for at least 30 minutes and eliminate them in accordance with the procedures of the laboratory.

STORAGE

Store the product at 2-8°C: in these conditions it maintains its validity until the expiry date indicated on the label.

REFERENCES

1. United States Pharmacopoeia latest edition.
2. European Pharmacopoeia latest edition.

PRESENTATION

PRODUCT	REF	PACKAGING	Spores: CFU / coupon	D _{VHP} (1,6 ± 0,5 mg/l, 30 °C)
Oxi Control GST E6 PLX	91085	20 coupons + 20 tubes	1-5 x 10 ⁶	1.5-4.0 minutes
Oxi Control GST E5 PLX	91086	20 coupons + 20 tubes	1-5 x 10 ⁵	1.5-4.0 minutes
Oxi Control GST E4 PLX	91048	20 coupons + 20 tubes	1-5 x 10 ⁴	1.5-4.0 minutes

TABLE OF SYMBOLS

LOT Batch code	 Do not reuse	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult accompanying documents	

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LIOFILCHEM® S.r.l.