




Integral System Yeasts Plus

System for the simultaneous identification and susceptibility testing of most clinically important yeasts



**NEW CHROMOGENIC
WELL (13-CHR)**

**ID and AST
IN ONE SINGLE TEST
RESULTS in 24/48 h.**



**EASY RESULT
INTERPRETATION BY
CLEAR COLOR CHANGE**

PERFORMANCE

- The susceptibility testing is according to the M27-A2 method of CLSI for evaluation of yeasts antimicrotics.
- Easy procedure of use
- Simultaneous identification and susceptibility testing in a single test.

NECESSARY ITEMS NOT CONTAINED IN THE PACKAGE

- VASELINE OIL (ref. 80278)
- Mc Farland 0.5 Barium Sulphate Standard (ref. 80400)

EXAMPLE OF RESULTS



IDENTIFICATION:

13-CHR: GREEN = *Candida albicans*

SUSCEPTIBILITY TESTING:

- **Sensibility:** all the antimicrotics.



IDENTIFICATION:

13-CHR: PINK = *Candida krusei*

SUSCEPTIBILITY TESTING:

- **Resistance:** *Nystatin* 13-NY.
- **Intermediate Sensibility:** *Amphotericin* 14-AMB; *Econazole* 16-ECN; *Clotrimoxazole* 18-CLO.
- **Sensibility:** all the other antimicrotics.



Take with a sterile loop one or more microbial colonies from a selective medium for the growth and isolation of yeasts.



Suspend the microbial colonies in *Physiological Solution* vial contained in the kit (**Suspension A**).



Suspension A

The **Suspension A** must have a 0,5 Mc Farland turbidity.



Suspension A

Suspension B

Dispense 0,02 mL (20 μ L) of **Suspension A** into another *Physiological Solution* vial contained in the kit (**Suspension B**).



Insert one xylose disc in the well **9-XYL** of the system.



Transfer 0,2 mL (4 drops) of **Suspension A** into wells from 1 to 13 for the **Identification test**.



Transfer 0,2 mL (4 drops) of **Suspension B** into wells from 14 to 24 for the **Susceptibility testing**.



Cover all wells of the system with two drops of vaseline oil and incubate at 36 ± 1 °C for 48 hours.



Write down the obtained results on the **Test Result Form**. Obtain the 5 digit code and identify the microorganism using the **Table of codes**.