Oxi/FermPluri-Test

System for the identification of Gram negative, oxidase positive bacteria.

DESCRIPTION

Oxi/FermPluri-Test is a 12-sector system containing special culture media for the identification of the Gram negative, oxidase positive bacteria.

The system allows the simultaneous inoculation of all media present in the sectors and the execution of 14 biochemical reactions.

Microorganism is identified evaluating the colour change of the different culture media after 48 hours of incubation at $36 \pm 1^\circ\text{C}$ and by a code number obtained from biochemical reaction interpretation.

CONTENT OF THE PACKAGE

Each package contains 10 or 25 **Oxi/FermPluri-Test**, 1 Instructions sheet and 1 data chart for biochemical reaction results.

ITEMS NECESSARY BUT NOT INCLUDED IN THE PACKAGE

KOVAC'S Reagent	Ref. 80271
Oxi/FermPluri-Test Codebook	Ref. 71708
OXIDASE TEST STICK / SWABS / DISC	Ref. 88029 / 88003 / 88004
Sundry microbiology laboratory materials	

CONFIGURATION

The configuration of the system is shown in Table no.1.

Table no.1

Sector	BIOCHEMICAL REACTIONS
Anaerobic Glucose	Glucose fermentation
Arginine	Arginine decarboxylation in anaerobiosis
Lysine	Lysine decarboxylation in anaerobiosis
Lactose / N₂	Lactose fermentation and nitrogen production in anaerobiosis
Sucrose / Indole	Sucrose oxidation and indole production
Xylose	Xylose oxidation
Aerobic Glucose	Glucose oxidation
Maltose	Maltose oxidation
Mannitol	Mannitol oxidation
Phenyl-alanine	Phenylalanine deamination
Urea	Urea hydrolysis
Citrate	Citrate utilization

PRINCIPLE OF THE METHOD

Oxi/FermPluri-Test makes possible the identification of Gram negative, oxidase positive bacteria isolated from clinical and environmental samples.

The identification is based on biochemical tests performed on culture media containing specific substrates. The combination of positive and negative reactions allows to build up a code number that permits to identify bacteria by using the **Codebook**.


SAMPLE COLLECTION

For the inoculation of **Oxi/FermPluri-Test** you must use a growth from non-selective agar media such as: Columbia Agar (Sheep Blood 5%), ref. 11025 or Tryptic Soy Agar (Sheep Blood 5%), ref. 11037. The isolated microorganism to be identified must be a pure culture of a Gram negative, oxidase positive bacteria.

TEST PROCEDURE

The microorganism to be identified should be recently isolated (18-24 hours): bacteria from cultures older than 48 hours can provide unreliable results.

Before inoculating the microorganism to be identified, it's compulsory to perform a gram staining and oxidase test on the microorganism. Only Gram negative, oxidase positive bacteria should be inoculated on **Oxi/FermPluri-Test**. For the correct performance of both tests please consult appropriate bacteriology manuals.

- Pick up an **Oxi/FermPluri-Test** system from the package and note identification name of bacterial strain to submit to identification, date of test and other useful information.
- Remove both caps of the system. Using the tip of inoculating needle, placed under the white cap, and without flaming, pick up a well isolated colony from a selective or non selective agar medium, without penetrating into the agar.
- Inoculate **Oxi/FermPluri-Test** turning and withdrawing the needle throughout the sectors of the system.
- Reinsert the needle with a turning movement until the breakage notch; break the inoculating needle folding it in correspondence with the notch. The portion of the needle remaining inside the system keeps anaerobic conditions necessary for reactions of the sectors **Anaerobic Glucose, Arginine, Lysine and Lactose/N₂**.
- Use the broken portion of the needle, remained in the user hands, to punch the plastic film in correspondence of the holes of the sectors **Sucrose/Indole, Xylose, Aerobic Glucose, Maltose, Phenyl-alanine, Urea, Citrate** in order to support aerobic growth.
- Screw again both caps and incubate **Oxi/FermPluri-Test** at $36 \pm 1^\circ\text{C}$ for 48 hours, putting it on its flat surface or vertically in a test-tube holder with the sector **Anaerobic Glucose** pointing upward.

INTERPRETATION OF RESULTS

After 24 hours incubation only read the urea reaction and note the results.

After 48 hours incubation:

- Watch for the color change of the culture media in the different sectors, interpret the results using the table no. 2 and, possibly, an **Oxi/FermPluri-Test** not inoculated and kept at room temperature.
- Write down the results, including the results oxidase, in the enclosed data chart, with the exception of indole test (**Sucrose/Indole** sector) that is executed later.
- **Indole test**
Place **Oxi/FermPluri-Test** with its flat surface pointing upward and, by punching the plastic film, inject with a syringe 3 or 4 drops of KOVAC'S Reagent in the sector **Sucrose/Indole**.
The reaction is positive if the added reagent turns pink-red in about 10-15 seconds. The color development is particularly visible at the point of syringe insertion.
- Form the 5-digit code following the instructions provided in the paragraph **CODE NUMBER FORMING**. Identify the bacterium using the **Codebook**.

Table no.2

Sector	BIOCHEMICAL REACTION	Sector color	
		Positive reaction	Negative reaction
Anaerobic Glucose	Glucose fermentation	Yellow	Green-blue
Arginine	Arginine decarboxylation in anaerobiosis	Purple	Yellow-grey
Lysine	Lysine decarboxylation in anaerobiosis	Purple	Yellow
Lactose	Lactose fermentation	Yellow	Red
- N₂	Nitrogen production in anaerobiosis	Wax detached	Wax bonded
Sucrose	Sucrose oxidation	Yellow	Green-blue
- Indole	Indole production	Red	Colorless
Xylose	Xylose oxidation	Yellow	Green-blue
Aerobic Glucose	Glucose oxidation	Yellow	Green-blue
Maltose	Maltose oxidation	Yellow	Green-blue
Mannitol	Mannitol oxidation	Yellow	Green-blue
Phenyl-alanine	Phenylalanine deamination	Light brown	Beige
Urea	Urea hydrolysis	Purple	Beige
Citrate	Citrate utilization	Blue	Green

CODE NUMBER FORMING

1) The 15 biochemical tests are divided into 5 groups each containing 3 tests and each one is indicated with a positivity value of 4, 2, 1.

- Value 4 : first test positive in each group (**Anaerobic Glucose, Lactose, Indole, Maltose, Urea**)
- Value 2 : second test positive in each group (**Arginine, N₂, Xylose, Mannitol, Citrate**)
- Value 1 : third test positive in each group (**Lysine, Sucrose, Aerobic Glucose, Phenyl-alanine, Oxidase**)
- Value 0 : every negative test

2) Adding the number of positive reactions in each group, it is obtainable a 5 digit code which, by the use of the **Codebook**, allows the identification of the microorganism under examination as in the following example.

Test	Gruppo 1			Gruppo 2			Gruppo 3			Gruppo 4			Gruppo 5		
	Anaerobic Glucose	Arginine	Lysine	Lactose	N ₂	Sucrose	Indole	Xylose	Aerobic Glucose	Maltose	Mannitol	Phenyl-alanine	Urea	Citrate	Oxidase
Positivity code	4	2	1	4	2	1	4	2	1	4	2	1	4	2	1
Results	-	+	-	-	+	-	-	+	+	-	-	-	0	+	+
Code	0+2+0=2			0+2+0=2			0+2+1=3			0+0+0=0			0+2+1=3		
CODE: 22303	IDENTIFICATION: <i>Pseudomonas aeruginosa</i>														

USER QUALITY CONTROL

Inoculate **Oxi/FermPluri-Test** using the reference bacterial strains indicated in the table no. 3.

For inoculation, incubation and reading please follow the instructions indicated in the paragraph **TEST PROCEDURE**.

Table no.3

Microorganism	Anaerobic Glucose	Arginine	Lysine	Lactose	N ₂	Sucrose	Indole	Xylose	Aerobic Glucose	Maltose	Mannitol	Phenyl-alanine	Urea	Citrate	Oxidase	Typical Biocodes
<i>Acinetobacter Iwoffii</i> ATCC 15309	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	00000
<i>Aeromonas hydrophila</i> ATCC 7966	+	+	-	+	-	+	+	-	+	+	+	-	-	±	+	65563 / 65561
<i>Bordetella bronchiseptica</i> ATCC 19395	-	-	-	-	-	-	-	-	-	-	-	-	+	+	+	00007
<i>Brevundimonas diminuta</i> ATCC 11568	-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	00011
<i>Plesiomonas shigelloides</i> ATCC 14029	±	+	+	+	-	-	+	-	+	+	-	-	-	-	+	74541 / 34541
<i>Pseudomonas aeruginosa</i> ATCC 27853	-	+	-	-	+	-	-	+	+	-	-	-	+	±	+	22307 / 22303
<i>Pseudomonas aeruginosa</i> ATCC 10145	-	+	-	-	+	-	-	+	+	-	-	-	+	+	+	22307
<i>Stenotrophomonas maltophilia</i> ATCC 13637	-	-	-	±	-	-	-	-	-	+	-	-	±	+	-	00042 / 04042 / 00046

TABLE OF BIOCHEMICAL REACTIONS**Table no.4 (Percentage of strains giving positive reactions with 18-24 h incubation at 36 ± 1°C)**

Microorganismo		Anaerobic Glucose	Arginine	Lysine	Lactose	N ₂	Sucrose	Indole	Xylose	Aerobic Glucose	Maltose	Mannitol	Phenyl-alanine	Urea	Citrate	Oxidase
<i>Acinetobacter</i>	<i>calcoaceticus</i>	18	0	6	0	0	0	0	100	100	0	0	0	7	10	0
	<i>iwoffii</i>	4	0	0	4	0	0	0	0	0	0	0	0	50	50	0
<i>Achromobacter</i>	<i>xylooxidans</i>	2	0	100	100	0	100	0	0	0	0	0	0	100	100	100
<i>Aeromonas</i>	<i>hydrophila</i>	100	100	67	100	89	100	100	89	100	89	11	0	0	89	100
<i>Alcaligenes</i>	<i>denitrificans</i>	0	0	100	0	0	0	0	0	0	0	0	0	100	100	100
<i>Bordatella</i>	<i>brochiseptica</i>	0	0	0	0	0	0	0	0	0	0	0	0	100	100	100
<i>Flavobacterium</i>	species	20	0	0	40	0	0	80	0	40	20	20	0	20	40	100
<i>Moraxella</i>	species	0	0	0	0	0	0	0	0	0	0	0	0	60	40	100
<i>Pasteurella</i>	<i>multocida</i>	0	0	0	83	0	50	33	0	0	0	0	0	0	0	100
<i>Plsesiomonas</i>	<i>shigelloides</i>	100	100	100	100	0	0	100	0	100	100	0	0	0	0	100
<i>Delftia</i>	<i>acidovorans</i>	0	0	100	0	0	0	0	0	0	0	0	0	22	100	100
<i>Burkholderia</i>	<i>cepacia</i>	0	8	100	0	0	0	0	0	42	58	0	0	8	83	100
<i>Pseudomonas</i>	<i>aeruginosa</i>	0	100	40	0	50	0	10	90	100	0	0	0	90	100	100
	<i>fluorescens</i>	0	0	0	0	0	0	0	80	100	0	0	0	0	100	100
	<i>mendocina</i>	0	100	100	0	0	0	0	0	100	0	0	0	50	50	100
	<i>putida</i>	4	92	88	0	0	0	0	80	88	0	0	0	56	100	100
	species	0	0	23	0	6	0	0	0	0	0	0	0	0	24	88
	<i>stutzeri</i>	0	0	0	0	67	0	0	33	33	0	0	0	0	100	100
<i>Ralstonia</i>	<i>pickettii</i>	0	0	0	0	0	0	0	100	100	0	0	0	100	100	100
<i>Stenotrophomonas</i>	<i>maltophilia</i>	0	0	88	0	0	0	0	0	0	75	0	13	88	100	75
<i>Vibrio</i>	<i>alginolyticus</i>	100	0	100	100	0	0	100	0	100	100	100	0	0	0	100
	<i>parahaemolyticus</i>	100	0	100	100	0	0	100	0	100	100	100	0	0	0	100

FACTORS THAT MAY INVALIDATE THE RESULTS

- Use of mixed cultures.
- Application of the method to other bacteria than Gram negative, oxidase positive bacteria.
- Use of expired systems.
- Test procedure different from the one suggested.

PRECAUTION

The product, **Oxi/FermPluri-Test**, cannot be classified as hazardous under current legislation and does not contain harmful substances in concentrations $\geq 1\%$. It therefore does not require a Safety Data Sheet to be available. **Oxi/FermPluri-Test** is a disposable device to be used only for *in vitro* diagnostic use; it is intended for use in a professional environment and should be used in laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8°C away from light. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using them if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **Oxi/FermPluri-Test** should be decontaminated and disposed off in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

BIBLIOGRAPHY

- Abbott, S.L. (2003) Klebsiella, Enterobacter, Citrobacter, Serratia, Plesiomonas, and other Enterobacteriaceae. In: Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenenbaum. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Murray P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenenbaum: (ed. 2003) Manual of clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Edwin, H.Lenette.: Manual of clinical Microbiology (1985), 4th Edition (ASM).
- Dito W. R., Bulmash J., Campbell J., Roberts E. : A numerical Coding and Identification System for Enterobacteriaceae. Chicago: American Society of Clinical Pathologists, Commission on Continuing Education, Illinois, 1972.
- Ewing, W. H.: Differentiation of Enterobacteriaceae by Biochemical Reactions. Washington, DC: U.S. dept. of Health, Education and Welfare, Public Health Service, National center for Disease Control, 1973.
- Oxi/FermPluri-Test Archivio Liofilchem, 2012.

PRESENTATION

Product	Ref.	Package
Oxi/FermPluri-Test	78620	10 test
	78621	25 test

TABLE OF SYMBOLS

 IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 LOT Batch code
 Store away from light				

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BIBLIOGRAFIA

- Abbott, S.L. (2003) Klebsiella, Enterobacter, Citrobacter, Serratia, Plesiomonas, and other Enterobacteriaceae. In: Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenenbaum. Manual of clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Murray P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenenbaum: (ed. 2003) Manual of clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Edwin, H.Lenette.: Manual of clinical Microbiology (1985), 4th Edition (ASM).
- Dito W. R., Bulmash J., Campbell J., Roberts E. : A numerical Coding and Identification System for Enterobacteriaceae. Chicago: American Society of Clinical Pathologists, Commission on Continuing Education, Illinois, 1972.
- Ewing, W. H.: Differentiation of Enterobacteriaceae by Biochemical Reactions. Washington, DC: U.S. dept. of Health, Education and Welfare, Public Health Service, National center for Disease Control, 1973.
- Oxi/FermPluri-Test Archivio Liofilchem, 2012.

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**LIOFILCHEM® S.r.l.**