Liofilchem®

Microbiology products

URITEST REF 500152/51015	URITEST N REF 500232/51023/51123	URITEST EC REF 500242/51024
URITEST EC REF 500412/51041	URITEST 2 REF 500302/51030/51130	URITEST MALTO REF 500402/51040/51140
URITEST EF REF 500702/51070/51170	URITEST M REF 500182/51018/51118	URITEST PENTA REF REF 500142/51014

Instructions For Use

Table of Contents

Italiano
English
Español
Français
Ελληνικά
Srpski
Dansk

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URITEST

DESCRIPTION

URITEST are systems that are ready for use with solidified culture media on a plastic support for creating the urine culture. These systems are universally known as dip slide. The results of the urine culture can be used to assist with diagnoses of infections of the urinary tract and must be used together with the patient's clinical history. URITEST by Liofilchem s.r.l. for urine culture are available in nine configurations: URITEST, URITEST N, URITEST C, URITEST EC, URITEST 2, URITEST MALTO, URITEST EF, URITEST M, URITEST PENTA. Each package is available in the 20 tests and 120 tests versions, while URITEST N, URITEST 2, URITEST MALTO, URITEST EF and URITEST M are also available in the package of 500 tests.

CONTENTS OF THE PACKAGES

Each package contains 20 cultural systems ready for use (dip slides) in the 20 tests version, 120 cultural systems ready for use (dip slide) in the 120 tests version and 500 cultural systems ready for use (dip slide) in the 500 tests version.

CONFIGURATIONS

URITEST for urine culture proposed by Liofilchem S.r.l. are available in the following configurations:

Table 1

URITEST ref 500152- ref 51015	URITEST N ref 500232 - ref 51023 - ref 51123	URITEST C ref 500242 - ref 51024
1. Mac Conkey Agar (red-purple)	1. Mac Conkey Agar (red-purple)	1. Mac Conkey MUG Agar (red-purple)
2. Cetrimide Agar (light amber)	2. Cetrimide Agar (light amber)	2. Cetrimide Agar (light amber)
3. Cled Andrade Agar (grey-green)	3. Cled Agar (green)	3. Cled Andrade Agar (grey-green)
URITEST EC ref 500412- ref 51041	URITEST 2 ref 500302 - ref 51030 - ref 51130	URITEST MALTO ref 500402 - ref 51040 - ref 51140
1. Mac Conkey Agar (red-purple)	1. Cled Agar (green)	1. Mac Conkey Agar (red-purple)
2. E.Coli Agar (grey)	2. Mac Conkey Agar (red-purple)	2. Malt Selective Agar (amber)
3. Cled Agar (green)		3. Cled Agar (green)
URITEST EF ref 500702 - ref 51070 - ref 51170	URITEST M ref 500182 - ref 51018 - ref 51118	URITEST PENTA ref 500142 - ref 51014
1. Cled Agar (green)	1. Mac Conkey Agar (red-purple)	1. Cled Agar (green)
2. Mac Conkey Agar (red-purple)	2. Malt Selective Agar (amber)	2. Mac Conkey Agar (red-purple)
3. Bile Aesculin Agar (dark amber)	3. Cled Andrade Agar (grey-green)	3. Malt Selective Agar (amber)
		4. Cetrimide Agar (light amber)
		5. Bile Aesculin Agar (dark amber)

METHOD PRINCIPLE

URITEST for urine culture are available in the variations with 3 culture media (URITEST, URITEST N, URITEST C, URITEST MALTO, URITEST EC, URITEST EF, URITEST M), with 2 culture media (URITEST 2) and 5 culture media (URITEST PENTA). URITEST N, inoculated through immersion in the urine enables the determination of the total microbe population on the Cled Agar medium, the gram-negative bacteria count on the Mac Conkey Agar and the count of the bacteria of *Pseudomonas* spp. on the Cetrimide Agar medium. URITEST is different from URITEST N because it substitutes for the Cled Agar medium the Cled Andrade Agar medium, which enables the Andrade indicator to be introduced to better differentiate the colonies. URITEST C differs from URITEST N by substituting Mac Conkey Agar MUG for Mac Conkey Agar; on the former medium *Escherichia coli* can be differentiated and the negative-gram bacteria can be counted. This differentiation is made possible by hydrolysis of the MUG by *E.coli* with the formation of fluorescent methylumbelliferone when the colonies are observed under the Wood lamp. URITEST EC differs from URITEST N because it has the medium E. Coli Agar instead of Cetrimide Agar. This enables *E. coli* to be distinguished from the other negative-gram bacteria. URITEST 2 is prepared with two cultural media, Cled Agar and Mac Conkey Agar and enables the total bacterial count (on Cled Agar) and the gram-negative bacteria (on Mac Conkey Agar) to be determined. URITEST MALTO differs from URITEST N as it replaces the medium Cetrimide Agar with the Malt Selective Agar medium, on which the fungi are determined. URITEST M differs from URITEST MALTO as it replaces the Cled Agar medium with the medium Cled Andrade Agar. URITEST PENTA carries 5 culture media and delivers total bacterial count (Cled Agar), gram negative bacteria (Mac Conkey Agar), yeasts and moulds (Malt Selective Agar), *Pseudomonas* spp. (Cetrimide Agar) and enterococci (Bile Aesculin Agar) growth.

TYPICAL FORMULATIONS OF THE MEDIA USED IN THE URITEST FOR URINE CULTURE (grams/litre)

CLED AGAR		MAC CONKEY AGAR		MAC CONKEY MUG AGAR	
Beef extract	3.0	Peptone	20.0	Peptone	20.0
Peptone	4.0	Lactose	10.0	Lactose	10.0
Tryptone	4.0	Bile salts 3	1.5	Bile salts n.3	1.5
Lactose	10.0	Sodium chloride	5.0	Sodium chloride	5.0
L-cystine	0.0128	Neutral red	0.03	Neutral red	0.03
Bromothymol blue	0.020	Crystal purple	0.001	Crystal purple	0.001
Agar	15.0	Agar	15.0	MUG	0.1
pH 7.3 ± 0.2		pH 7.1 ± 0.2		Agar	15.0
				pH 7.1± 0.2	
CLED ANDRADE AGAR		CETRIMIDE AGAR		BILE AESCULIN AGAR	
Beef extract	3.0	Peptogelatin	20.0	Peptone	17.0
Peptone	4.0	Magnesium chloride	1.4	Tryptone	3.0
Tryptone	4.0	Sulphate of potassium	10.0	Yeast extract	5.0
Lactose	10.0	Cetrimide	0.3	Ox bile	10.0
L-cystine	0.128	Agar	15.0	Sodium chloride	5.0
Bromothymol blue	0.020	Glycerol	10.0 ml	Aesculine	1.0
Andrade	0.1	pH 7.0 ± 0.2		Ferric ammonium citrate	0.5
Agar	15.0			Thallium acetate	0.9
pH 7.5 ± 0.2				Agar	15.0
				pH 7.3 ± 0.2	
MALT SELECTIVE AGAR		E. COLI AGAR			
Malt extract	30.0	Peptone	10.0		
Agar	17.0	Peptospecial	5.0		
pH 4.5 ± 0.2		Yeast extract	3.0		
		Lactose	12.5		
		Bile salts number 3	1.5		
		Sodium chloride	5.0		
		Aniline blue	0.1		
		Agar	15.0		
		pH 7.4 ± 0.2			

COLLECTION AND CONSERVATION OF SAMPLES

Collect the urine at least 6 hours after the last urysis (it is preferable to take a sample from the first urine passed that day). The recommended collection methods are *mitto intermedio* for continent adults and adhesive plastic bags for infants. Samples must be obtained from bladder catheterism and from suprapubic injection. The samples must be collected in sterile containers of suitable volume. The urine samples should be inoculated on the dip slide immediately after being taken. If this is not possible, keep the samples in a refrigerator at 2-8°C for no longer than 24 hours. Sterile containers containing preservatives are available for collecting urine. If these containers are used, follow the manufacturer's instructions. The urine samples should be taken before any antibacterial treatment or at least 48 hours after the last dose of antimicrobial pharmaceutical is administered.


TEST PROCEDURE

1. Choose the most suitable URITEST for the required use, taking into account the instructions in the previous chapters.
2. Mark the sample with an adhesive label or write the information directly on the container of the dip slide.
3. Unscrew and extract the slide from the cylindrical container without touching the surfaces of the culture media.
4. Hold the slide by the cap and immerse it in the urine sample that has just been taken in such a way that the surfaces of the culture media are completely immersed.
5. Let the excess urine drip off from the slide.
6. Use absorbent paper to mop up the last drops of urine.
7. Reinsert the slide in the container and screw it tight.
8. Incubate in a thermostat set at 36 ± 1°C for 18-24 hours.

INTERPRETING RESULTS

After incubation, the presence of bacteria is indicated by the colonies appearing on the surfaces of the culture media. The number of colonies indicates the concentration of Colony Forming Units per ml of urine (CFU/ml). The total microbial load of the urine sample is measured on Cled Agar and Cled Andrade Agar. There is a correlation between microbial load and the density of the colonies developed on the surface of the culture media. This correlation is set out in Figure 1, which can be taken as a reference for assessing the number of CFU/ml.

Figure 1: Correlation between density of microbial growth and CFU/ml on Cled Agar and Cled Andrade Agar media.

The culture media combined in the different formats with Cled Agar and Cled Andrade Agar are used to assess the quality and quantity of the different microbial strains that sustain the urinary infection. Table 2 shows an interpretation guide of the colonies on the different cultures present in the URITEST in relation to the growth characteristics.

Table 2: Reading and interpretation of the colonies grown on different culture media

Medium	<i>Escherichia coli</i>	<i>Proteus spp.</i>	<i>Klebsiella spp.</i>	<i>Candida spp.</i>	<i>Pseudomonas spp.</i>	<i>Enterococcus spp.</i>	<i>Stafilococcus spp.</i>
Cled Agar	Yellow	Green	Yellow	White	Green	Yellow	Yellow
Cled Andrade Agar	Red	Green	Pink	White	Green	Pink	Pink
Mac Conkey Agar	Red	Colourless	Pink	Inhibited	Brown	Inhibited	Inhibited
Mac Conkey MUG Agar	Red, fluorescent ⁽¹⁾	Colourless	Pink	Inhibited	Brown	Inhibited	Inhibited
Cetrimide Agar	Inhibited	Inhibited	Inhibited	Inhibited	Yellow-green/ Colourless	Inhibited	Inhibited
Bile Aesculine Agar	Inhibited	Inhibited	Inhibited	Inhibited	Inhibited	Black	inhibited
Malt Selective Agar	Inhibited	Inhibited	Inhibited	White	Inhibited	Inhibited	Inhibited
E. Coli Agar	Blue	Colourless	Grey	Inhibited	Colourless	Inhibited	Inhibited

(1) Fluorescence blue-green is observed in a dark room under a Wood lamp

CLINICAL INTERPRETATION

As the urethra is not a sterile site but is colonised by different types of microorganisms, the presence of colonies on the URITEST does not in itself indicate an infection of the urinary tract unless the sample has been taken in sterile conditions by suprapubic injection. To interpret the quantitative results of the urine culture, the following guidelines can be considered.

Table 3

Collection Method	CFU/ml microbial load	Meaning
MITTO INTERMEDIO	> 10 ⁵	Significant bacteruria: sign of infection
	from 10 ⁴ to 10 ⁵	Borderline value: request second sample. In women of childbearing age the lower urinary tract may be infected (dysuria-polysuria syndrome).
	< 10 ⁴	Generally considered to be a sign of contamination.
SUPRAPUBIC ASPIRATION	Any microbial load	Any bacteria or fungus even if small in number, must be reported because of their potential pathogenic action.
CATHETERISM	> 10 ⁵	Significant bacteruria: sign of infection
	< 10 ⁵	Any bacteria or fungus, even if small in number, must be reported because of their potential pathogenic action.

QUALITY CONTROL

For quality control of the URITEST, overnight broth cultures are used that are diluted in a sterile physiological solution so as to obtain approximately 10⁵-10⁶ cells/ml of the following strains: *Staphylococcus aureus* ATCC® 25923, *Escherichia coli* ATCC® 25922, *Pseudomonas aeruginosa* ATCC® 27853, *Enterococcus faecalis* ATCC® 19433, *Candida albicans* ATCC® 10231. The slides are inoculated through immersion in the diluted broth culture and are incubated at 36 ± 1°C for 18-24 hours. Read and interpret the results as described in Figure 1 and in Table 2. The appearance of the media must also be checked: for the colour of the media see Table 1 of the configuration chapter.

LIMITS AND NOTICES

1. Instruct the patient properly on collecting urine as sample contamination values between 5% and 40% have been reported (Valenstein P. et al., 1998).
2. Although most microorganisms grow within 24 hours of incubation the sample may contain slow-growing microorganisms or substances with antimicrobial activity that prevents growth during overnight incubation. Urine cultures that are negative after 24 hours should be incubated for another 24 hours before the final result is supplied.
3. A large quantity of epithelial cells in the sample indicates possible contamination with periurethral and/or vaginal flora. In such cases the urine sampling and analysis must be repeated.

PERFORMANCE

The approximate determination of the microbial load of the urine using the dip-slide system correlates well with the standard microbial plate count (Rosenberg, M, 1992).

PRECAUTIONS

The Bile Aesculine Agar medium contains Thallium acetate that is classed as harmful by current legislation; see the safety datasheet on use. The other media described are not classified as hazardous by current legislation and do not contain harmful substances in concentrations of $\geq 1\%$. No safety datasheets are therefore required for them. URITEST are disposable devices. URITEST are used only for *in vitro* diagnosis. They are designed for professional use and must be used in laboratories by properly trained operators using approved aseptic and safety methods to deal with pathogenic agents. Sterilise the inoculated matter after use and before eliminating it as waste.

DIP- SLIDE CONSERVATION

Conserve at 10-25°C in their original packaging that is protected against the light. Do not keep close to sources of light and avoid excessive temperature variations. In these conditions the products described here will be valid until the expiry date shown on the label. Do not use beyond this date. Eliminate if there signs of deterioration (noticeable changes in the colour of media or detachment from the sport, bacterial colonies or fungus before use). DO NOT FREEZE.

DISPOSING OF USED MATERIAL

After use, URITEST and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with laboratory techniques for decontamination and disposal of potentially infected material.

BIBLIOGRAPHY

- Manual of Clinical Microbiology. Ed. In Chief Patrick R. Murray. 2003. 8 th Ed., Vol 1, pp 320-330. ASM Press.
- Rosenberg M., S. A. Berger, M.Barki, S.Goldberg, A.Fink, and A. Cunha. 1992. Initial testing of a novel urine culture device. J.Clin. Microbiol. 30:2686-2691.
- Valenstein P., and F.Meier. 1998. Urine culture contamination: a College of American Pathologists Q-Probes study of contaminated urine culture in 906 institutions. Arch. Pathol. Lab. Med. 122: 123-129.

PRESENTATION

PRODUCT	REF	PACKAGE	PRODUCT	REF	PACKAGE
URITEST	500152	20 Dip Slide	URITEST MALTO	500402	20 Dip Slide
	51015	120 Dip Slide		51040	120 Dip Slide
				51140	500 Dip Slide
URITEST N	500232	20 Dip Slide	URITEST EF	500702	20 Dip Slide
	51023	120 Dip Slide		51070	120 Dip Slide
	51123	500 Dip Slide		51170	500 Dip Slide
URITEST C	500242	20 Dip Slide	URITEST M	500182	20 Dip Slide
	51024	120 Dip Slide		51018	120 Dip Slide
URITEST EC	500412	20 Dip Slide		51118	500 Dip Slide
	51041	120 Dip Slide	URITEST PENTA	500142	20 Dip Slide
URITEST 2	500302	20 Dip Slide		51014	120 Dip Slide
	51030	120 Dip Slide			
	51130	500 Dip Slide			

TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Consult instructions for use	 Do not reuse

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
- Manual of Clinical Microbiology. Ed. In Chief Patrick R. Murray. 2003. 8 th Ed., Vol 1, pp 320-330. ASM Press.
- Rosenberg M., S. A. Berger, M. Barki, S. Goldberg, A. Fink, and A. Cunha. 1992. Initial testing of a novel urine culture device. J.Clin. Microbiol. 30:2686-2691.
- Valenstein P., and F. Meier. 1998. Urine culture contamination: a College of American Pathologists Q-Probes study of contaminated urine culture in 906 institutions. Arch. Pathol. Lab. Med. 122: 123-129.

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PRODUCTO	REF	ESTUCHE	PRODUCTO	REF	ESTUCHE
URITEST	500152	20 Dip Slide	URITEST MALTO	500402	20 Dip Slide
	51015	120 Dip Slide		51040	120 Dip Slide
				51140	500 Dip Slide
URITEST N	500232	20 Dip Slide	URITEST EF	500702	20 Dip Slide
	51023	120 Dip Slide		51070	120 Dip Slide
	51123	500 Dip Slide		51170	500 Dip Slide
URITEST C	500242	20 Dip Slide	URITEST M	500182	20 Dip Slide
	51024	120 Dip Slide		51018	120 Dip Slide
URITEST EC	500412	20 Dip Slide	URITEST PENTA	51118	500 Dip Slide
	51041	120 Dip Slide		500142	20 Dip Slide
URITEST 2	500302	20 Dip Slide	51014	120 Dip Slide	
	51030	120 Dip Slide			
	51130	500 Dip Slide			

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BIBLIOGRAFI

- Manual of Clinical Microbiology. Ed. In Chief Patrick R.Murray. 2003. 8th Ed., Vol 1, pp 320-330. ASM Press.
- Rosenberg M., S.A. Berger, M. Barki, S. Goldberg, A. Fink, and A. Cunha. 1992. Initial testing of a novel urine culture device. J.Clin. Microbiol. 30:2686-2691.
- Valenstein P., and F. Meier. 1998. Urine culture contamination: a College of American Pathologists Q-Probes study of contaminated urine culture in 906 institutions. Arch.Pathol. Lab. Med. 122: 123-129.

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LIOFILCHEM® S.r.l.