Czapek Dox Agar

Medium for cultivating fungi and bacteria
from clinical specimens and environmental samples.

DESCRIPTION

Czapek Dox Agar is a medium used for the cultivation of fungi and bacteria capable of using inorganic nitrogen.

Czapek Dox Agar is recommended by APHA for isolating *Aspergillus*, *Penicillium*, *Paecilomyces* and related fungi, yielding characteristic mycelia and conidia. The only bacteria which are able to grow in this medium are non-fastidious soil bacteria.

TYPICAL FORMULA

| | (g/l) |
|----------------------------|-------|
| Saccharose | 30.0 |
| Sodium Nitrate | 2.0 |
| Dipotassium Phosphate | 1.0 |
| Magnesium Sulfate | 0.5 |
| Potassium Chloride | 0.5 |
| Ferrous Sulfate | 0.01 |
| Agar | 15.0 |
| Final pH 7.3 ± 0.2 at 25°C | |

METHOD PRINCIPLE

Saccharose is the sole source of carbon and sodium nitrate is the sole nitrogen source. Dipotassium phosphate is the buffering agent. Magnesium sulfate, potassium chloride and ferrous sulfate provide essential ions. Agar is the solidifying agent.

PREPARATION

Dehydrated medium Suspend 49 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.

Notice that the acidity of the medium may be increased to pH 3.5 ± 0.2 for the cultivation of acidophilic organisms, such as yeasts, by adding 10 ml of 10% lactic acid per liter after sterilization.

Medium in bottles Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

TEST PROCEDURE

Inoculate the medium by spreading the sample over the agar surface. Time and temperature of incubation vary considerably according to the species being cultivated. As a general guide, incubate for 1-2 weeks at 25°C. Most *Penicillium* spp have an optimum growth temperature between 20°C and 25°C, whilst many *Aspergillus* spp grow best at about 30°C. However, different fungi grow over a wide range of temperatures.

INTERPRETING RESULTS

Examine for microbial growth and chlamydospore production.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, very light beige.

Prepared medium: slightly opalescent, light amber, may have a slight precipitate.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.
 Medium in bottles: 2 years.
 Ready-to-use plates: 6 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.
 Inoculum for productivity: 50-100 CFU
 Incubation conditions: aerobically at 25-30°C for 48-72 hours..

QC Table.

| Microorganism | | Growth |
|---------------------------------|-------------|--------|
| <i>Aspergillus brasiliensis</i> | ATCC® 16404 | Good |
| <i>Candida albicans</i> | ATCC® 10231 | Good |
| <i>Saccharomyces cerevisiae</i> | ATCC® 9763 | Good |

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. Eaton, Rice and Baird (2005) Standard methods for the microbiological examination of water and wastewater, 21st ed. American Public Health Association, Washington D.C.
2. Dox (1910) U.S. Dept. Agr. Bur. Anim. Ind. Bull. 120:70
3. Czapek (1902-1903) Beitr. Chem. Physiol. Pathol. 1:540.

| PRESENTATION | | Contents | Ref. |
|-----------------|---------------------------|--------------------|--------|
| Czapek Dox Agar | 90 mm ready-to-use plates | 20 plates | 10017 |
| Czapek Dox Agar | Bottles | 6 x 500 ml bottles | 470220 |
| Czapek Dox Agar | Dehydrated medium | 500 g of powder | 610095 |
| Czapek Dox Agar | Dehydrated medium | 100 g of powder | 620095 |

TABLE OF SYMBOLS

| | | | | |
|-----------------------------|--|---|---|---|
| LOT Batch code | IVD <i>In vitro</i> Diagnostic Medical Device |  Manufacturer |  Use by |  Fragile, handle with care |
| REF Catalogue number |  Temperature limitation |  Contains sufficient for <n> tests |  Caution, consult Instruction For Use |  Do not reuse |

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