## CLED Agar

Medium for the isolation, enumeration and differentiation of pathogenic bacteria in urine specimens.

### DESCRIPTION

CLED Agar is a medium used for microbiological urine analysis since it supports the growth of all urinary potential pathogens providing good colony differentiation.

### TYPICAL FORMULA

	(g/l)
Enzymatic Digest of Gelatin	4.0
Enzymatic Digest of Casein	4.0
Beef Extract	3.0
Lactose	10.0
L-Cystine	0.128
Bromothymol Blue	0.02
Agar	15.0
Final pH 7.3 ± 0.2 at 25°C	

### METHOD PRINCIPLE

Enzymatic digest of gelatin, enzymatic digest of casein and beef extract provide amino acids, nitrogen, carbon, vitamins and minerals required for organisms growth. Lactose is the fermentable carbohydrate. L-cystine is a growth supplement for cystine-dependent organisms. Bromothymol blue is the pH indicator changing color from green to yellow when lactose fermentation lowers the pH. Agar is the solidifying agent. Lack of electrolytes suppresses the swarming of *Proteus* and *Shigella* species.

### PREPARATION

<u>Dehydrated medium</u>	Suspend 36.1 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.
<u>Medium in bottles</u>	Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

### TEST PROCEDURE

Urine must be directly streaked over the agar surface no later than 2 h after collection or must be kept at 2-8°C (not longer than 24 h) to avoid microbial overgrowth. Use a calibrated loop (0.01 or 0.001 ml) to inoculate the medium with the undiluted, well-mixed urine sample. Incubate aerobically at 35 ± 2°C for 18-48 h.

### RESULTS INTERPRETATION

Count the number of colonies on the plate. Each colony correspond to 100 or 1000 CFU/ml of urine, using a 0.01 ml or 0.001 ml loop, respectively. Observe the color and the morphology of the colonies for presumptive identification according to the ID table. Further tests should be performed for confirmation.

#### ID Table.

Microorganism	Appearance of the colonies
<i>Escherichia coli</i>	Yellow, opaque
<i>Klebsiella, Enterobacter</i> spp	Yellow to whitish-blue, mucoid
<i>Proteus</i> spp	Translucent blue
<i>Salmonella</i> spp	Flat blue
<i>Pseudomonas aeruginosa</i>	Green, matt surface and rough periphery
<i>Enterococcus</i> spp	Yellow, about 0.5 mm in diameter
<i>Staphylococcus aureus</i>	Deep yellow
Coagulase negative staphylococci	Pale yellow

### APPEARANCE

Dehydrated medium: free-flowing, homogeneous, light beige.  
Prepared medium: slightly opalescent, green to blue-green.

### STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

**SHELF LIFE**

Dehydrated medium: 4 years.  
 Medium in bottles: 2 years.  
 Ready-to-use plates: 6 months.

**QUALITY CONTROL**

Slides are inoculated with the microbial strains indicated in the QC table.  
 Inoculum for productivity: 50-100 CFU  
 Inoculum for selectivity: 10<sup>4</sup>-10<sup>6</sup> CFU  
 Incubation conditions: aerobically at 35 ± 2°C for 18-24 hours.

**QC Table.**

Microorganism		Growth	Specification
<i>Escherichia coli</i>	ATCC® 25922	Good	Yellow colonies
<i>Klebsiella pneumoniae</i>	ATCC® 13883	Good	Yellow to whitish-blue colonies
<i>Proteus mirabilis</i>	ATCC® 12453	Good	Blue colonies, no swarming
<i>Salmonella</i> Typhimurium	ATCC® 14028	Good	Blue colonies
<i>Pseudomonas aeruginosa</i>	ATCC® 27853	Good	Green colonies
<i>Enterococcus faecalis</i>	ATCC® 29212	Good	Yellow colonies
<i>Staphylococcus aureus</i>	ATCC® 25923	Good	Small yellow colonies

**WARNING AND PRECAUTIONS**

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use only and must be used by properly trained operators only.

**DISPOSAL OF WASTE**

Disposal of waste must be carried out according to national and local regulation in force.

**BIBLIOGRAPHY**

- Sandys G.H. (1960) A new method for preventing swarming of *Proteus* spp with a description of a new medium suitable for use in routine laboratory practice. J. Med. Lab. Technol. 17:224.
- Mackey J.P. and G.H. Sandys (1966) Diagnosis of urinary tract infections. Br. Med. J. 1:1173.

**PRESENTATION**

		Contents	Ref.
CLED Agar	90 mm ready-to-use plates	20 plates	10026
CLED Agar	90 mm ready-to-use plates	100 plates	10026*
CLED Agar	140 mm ready-to-use plates	10 plates	10236
CLED Agar	Bottles	6 x 100 ml bottles	402180
CLED Agar	Bottles	6 x 200 ml bottles	412180
CLED Agar	Bottles	6 x 500 ml bottles	470110
CLED Agar	Dehydrated medium	500 g of powder	610012
CLED Agar	Dehydrated medium	100 g of powder	620012
CLED Agar	Dehydrated medium	5 kg of powder	6100125

**TABLE OF SYMBOLS**

<b>LOT</b> Batch code	<b>IVD</b> <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
<b>REF</b> Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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