SS Agar

Selective medium for isolation of *Salmonella* spp and some strains of *Shigella* spp, from clinical samples and other materials.

DESCRIPTION

SS Agar is a selective medium used for the primary isolation of *Salmonella* spp from very contaminated samples. The medium complies with the recommendations of the APHA for the examination of food.

TYPICAL FORMULA	(g/l)
Pancreatic Digest of Casein	2.5
Peptic Digest of Animal Tissues	2.5
Beef Extract	5.0
Lactose	10.0
Sodium Thiosulfate	8.5
Sodium Citrate	8.5
Bile Salts	8.5
Ferric Ammonium Citrate	1.0
Brilliant Green	0.00033
Neutral Red	0.025
Agar	15.0
Final pH 7.0 ± 0.2 at 25°C	

METHOD PRINCIPLE

Pancreatic digest of casein, peptic digest of animal tissues and beef extract provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Lactose is the fermentable carbohydrate. Bile salts, sodium citrate and brilliant green inhibit Gram-positive-bacteria and several Enterobacteriaceae other than *Salmonella* and some species of *Shigella*. Sodium thiosulfate and ferric ammonium citrate enable the detection of hydrogen sulfide production as evidenced by colonies with black center. Neutral red is the pH indicator. Agar is the solidifying agent.

PREPARATION

Medium in bottles Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

TEST PROCEDURE

Inoculate the plates by directly streaking the specimen on the agar surface or spread the sample from an enrichment culture. Incubate aerobically at 35 ± 2°C for 18-24 h

INTERPRETING RESULTS

Lactose-nonfermenting organisms, such as *Salmonella* and *Shigella* spp, form colorless colonies and the strains of *Salmonella* which produce hydrogen sulfide are visualized as colonies with black center.

Also *Proteus* spp ferment lactose and produce H₂S developing translucent colonies, smaller than *Salmonella*, with black center.

Lactose-fermenting organisms, such as *Klebsiella* spp, grow as pink to red colonies with or without a zone of precipitated bile.

Enterococci, Staphylococci and other Gram-positive bacteria are partially or completely inhibited.

APPEARANCE OF THE MEDIUM

Slightly opalescent, red-orange.

STORAGE

Store at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Medium in bottles: 2 years.

Ready-to-use plates: 6 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Inoculum for selectivity: 10⁴-10⁶ CFU.

Incubation conditions: aerobically at 35 ± 2°C for 18-24 hours.

QC Table.

Microorganism		Growth	Specification
<i>Salmonella</i> Typhimurium	ATCC® 14028	Good	Colorless colonies with black center
<i>Shigella flexneri</i>	ATCC® 12022	Good	Colorless colonies
<i>Proteus mirabilis</i>	ATCC® 12453	Good	Colorless colonies with black center
<i>Klebsiella pneumoniae</i>	ATCC® 13883	Good	Pink colonies
<i>Escherichia coli</i>	ATCC® 25922	Poor or inhibited	Pink to red colonies
<i>Enterococcus faecalis</i>	ATCC® 29212	Inhibited	---

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. Murray, Baron, Jorgensen, Landry and Pfaller ed. (2007) Manual of clinical microbiology, 9th ed. American Society for Microbiology, Washington, D.C.
2. American Public Health Association (1992) Compendium of methods for the microbiological examination of food - 3rd ed.
3. Rose, H.M., and M.H. Kolodny (1942). J. Lab. Clin. Med. 27: 1081-1083.
4. Leifson E. (1935). J. Pathol. Bacteriol. 40: 581.

PRESENTATION		Contents	Ref.
SS Agar	90 mm ready-to-use plates	20 plates	10036
SS Agar	90 mm ready-to-use plates	100 plates	10036*
SS Agar	Bottles	6 x 100 ml bottles	402300
SS Agar	Bottles	6 x 200 ml bottles	412300

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Medical Diagnostic Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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