Legionella Agar (GVPC)

Selective medium for detection and enumeration of *Legionella* spp.

DESCRIPTION

Legionella Agar (GVPC) is a selective medium used for primary isolation of *Legionella* species from water samples and water related matrices according to ISO 11731.

This medium is also used for detection of *Legionella* in clinical specimens.

TYPICAL FORMULA* (Per litre of Purified Water)

Yeast Extract	10.0 g
Charcoal Activated	2.0 g
ACES Buffer	10.0 g
Potassium Hydroxyde	2.8 g
α-Ketoglutarate	1.0 g
Ferric Pyrophosphate	0.25 g
L-Cysteine	0.4 g
Agar	13.0 g
Glycine	3.0 g
Vancomycin	0.001 g
Polymyxin B	80000 IU
Cycloheximide	0.08 g

Final pH 6.9 ± 0.2 at 25°C

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Yeast extract provides amino acids, nitrogen, carbon, vitamins and minerals. Activated charcoal decomposes hydrogen peroxide, a metabolic product toxic to *Legionella* spp, and may also collect carbon dioxide and modify surface tension. ACES buffer (N-2- acetamido-2-aminoethane sulfonic acid) and potassium hydroxide maintain the proper pH for optimal growth. Alpha-ketoglutarate, cysteine and ferric pirophosphate are incorporated to satisfy the specific nutritional requirements of Legionella species. Agar is the solidifying agent. Glycine, vancomycin and polymyxin B inhibit or suppress most non-target bacterial species, both Gram-positive and Gram-negative, including common contaminants such as enterococci, coliform, and *Pseudomonas* spp, while cycloheximide suppresses the growth of yeasts and moulds.

TEST PROCEDURE

For use in clinical analysis, inoculate respiratory specimens directly by rolling the swab over the agar surface (90 mm plate) to obtain isolated colonies or following centrifugation. If contamination flora is abundant, it is recommended to treat specimens by acidification (HCl 0.2M), followed by neutralization (NaOH 0.01M) and/or heating (30 minutes at 50°C).

Contact plates are suitable for environmental monitoring, e.g. surfaces and aerosols generated from air conditioning systems. Selected surfaces are sampled by firmly pressing the agar medium against the test area for about 10 seconds. Do not move laterally the plate. Residues of the medium should be subsequently removed from the area tested. For active air sampling, place the plate, lid up, in the air sampler, remove the lid of the plate and sample a specific volume of air according to laboratory procedure.

For water testing according to ISO 11731, the choice of the method is down to the individual laboratory, as it depends on the type of the sample (e.g., whether the water has low, high, or extremely high bacterial background flora) and the reason of the investigation. To ensure detection, water samples may be concentrated by membrane filtration or, alternatively, by centrifugation (when the number of legionellae in any given sample is not known, concentration technique are usually performed). Dilution is necessary when high concentrations of *Legionella* and/or other bacteria are expected. Heat treatment, acid treatment, or a combination of both may be also required before culturing on selective media. For more details refer to the ISO standard.

Incubate inoculated plates inverted at 36 ± 2°C for 7 to 10 days in humidified atmosphere (air with 2.5% CO₂ can be beneficial for growth of some *Legionella* but is not essential).

Regarding sample transport and storage, analysis should begin as soon as possible after the sample has been taken, preferably on the same day. Samples should be protected from heat and sunlight, kept between 6-20°C and transported to the laboratory within 24 hours of collection.

INTERPRETING RESULTS

Examine for growth and fluorescence under long-wave UV light at regular intervals during the incubation period starting from the second day.

Colonies of *Legionella* are white-grey in general but can also appear in other colours. They are smooth with an entire edge and exhibit a characteristic ground-glass appearance. Under UV light, colonies usually exhibit brilliant white fluorescence.

Select colonies characteristic of *Legionella* on each plate for subculture onto plates of Legionella BCYE Agar (90 mm, ref. 10051) and Legionella BCYE Agar w/o cysteine (90 mm, ref. 10412). Incubate at $36 \pm 2^\circ\text{C}$ for at least 2 days in humidified atmosphere. Regard as *Legionella* those colonies which grow on Legionella BCYE Agar but fail to grow on Legionella BCYE Agar w/o cysteine.

APPEARANCE

Opaque, black.

STORAGE

Store at $2-8^\circ\text{C}$ away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

Avoid quick temperature shifts to prevent condensation.

SHELF LIFE

6 months.

QUALITY CONTROL

To check the performance of the medium the following microbial strains are used.

Strain		Inoculum	Incubation	Criteria	Specification
<i>Legionella pneumophila</i>	WDCM 00107	50-100 CFU	2-5 d / $36 \pm 2^\circ\text{C}$	$P_R \geq 0.5$	White-grey-blue-purple colonies with entire edge and characteristic ground-glass appearance
<i>Legionella anisa</i>	WDCM 00106		5-10 d / $36 \pm 2^\circ\text{C}$		
<i>Enterococcus faecalis</i>	WDCM 00087	10^4 - 10^6 CFU	3 d / $36 \pm 2^\circ\text{C}$	Total inhibition	—
<i>Escherichia coli</i>	WDCM 00013			Total or partial inhibition	—

A productivity ratio (P_R) of 0.5 is equivalent to a recovery rate of 50%
BCYE agar used as reference medium

LIMITATIONS

- Some strains of *Pseudomonas*, Enterobacteria or yeasts may develop on this medium but do not produce typical colonies.
- Some strains of *Legionella* may not develop on this medium due to their growth requirements and/or antimicrobial sensitivity.
- For use in medical microbiology, this medium is intended as an aid to the diagnosis of the infection: the interpretation of the results must be made considering the patient's clinical history, the origin of the sample and the results of the microscopic and/or other diagnostic tests.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the various configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Legionella Agar (GVPC)	90 mm Plate	20 plates	10128 ^a
Legionella Agar (GVPC)	90 mm Plate	100 plates	10128* ^a
Legionella Agar (GVPC)	60 mm Plate (membrane placement)	20 plates	163392 ^b
Legionella Agar (GVPC)	Contact Plate	20 plates	15376 ^b

^a For *In Vitro* Diagnostic Use (CE marked)

^b Not intended for the diagnosis of disease or other conditions in humans (non-CE marked)

This document is available from the online Support Center:

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liofilchem.com/ifu-sds

BIBLIOGRAPHY / BIBLIOGRAFIA

1. EN ISO 11133:2014+Amd1:2018. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. ISO 11731:2017. Water quality – Enumeration of *Legionella*.
3. Montagna M.T. et al. (2014) Legionella spp contamination in indoor air: preliminary results of an italian multicenter study. Epidemiol Prev; 38(6) Suppl 2: 62-65. <https://core.ac.uk/download/pdf/53299093.pdf>
4. Edelstein P.H. (1981) Improved semiselective medium for the isolation of *Legionella pneumoniae* from contaminated clinical and environmental specimens. J. Clin. Microbiol. 14(3):298.

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