Mueller Hinton II Agar (Sheep Blood 5%)

Medium for antimicrobial susceptibility testing of certain fastidious organisms by the disc diffusion technique, as standardized by the CLSI.

DESCRIPTION

Mueller Hinton II Agar (Sheep Blood 5%) is a medium recommended for antimicrobial susceptibility testing of *Streptococcus* spp. and *Neisseria meningitidis* by the disc diffusion technique (Kirby-Bauer method), as standardized by the Clinical and Laboratory Standard Institute (CLSI).

TYPICAL FORMULA

	(g/l)
Beef Extract	2.0
Acid Hydrolysate of Casein	17.5
Starch	1.5
Agar	17.0
Sheep Blood, defibrinated	50 ml
Final pH 7.3 ± 0.2 at 25°C	

METHOD PRINCIPLE

Acid hydrolysate of casein and beef extract supply amino acids, nitrogen, minerals, vitamins, carbon and other nutrients which support the growth of microorganisms. Starch acts as a protective colloid against toxic molecules which can be present in the medium. Hydrolysis of starch during autoclaving supplies a little amount of glucose, which is a source of energy. Agar is the solidifying agent. Sheep blood is added to support the growth of fastidious organisms.

The Kirby-Bauer method is based on the diffusion, through the agar, of the antimicrobial substance which soaks the paper disc. Each disc has a single concentration of the antimicrobial agent that inhibits the microorganism growth showing an halo around the disc. The diameter of the inhibiting halo is correlated with the Minimal Inhibitory Concentration (MIC).

Mueller Hinton II Agar (Sheep Blood 5%) is manufactured to contain low levels of thymine and thymidine and controlled levels of calcium and magnesium ions.

TEST PROCEDURE

1. Prepare a standardized suspension of the test organism using either the direct colony suspension or growth method.
2. Dip a sterile cotton swab into the adjusted suspension.
3. Inoculate the surface of the plate by streaking the swab over the entire agar surface.
4. Dispense antimicrobial discs onto the surface of the inoculated agar plate.
5. Incubate in an atmosphere of 5% CO₂ at 35 ± 2°C for 20-24 hours.

INTERPRETING RESULTS

After incubation measure the diameter of the zone of complete inhibition including the diameter of the disc. Interpret zone sizes with updated tables for the CLSI antimicrobial susceptibility testing standards and report the organism as susceptible, intermediate or resistant to the agents that have been tested. For streptococci, measure the diameter of the zone of growth inhibition, not of the zone of inhibition of hemolysis.

APPEARANCE

Opaque, cherry red.

STORAGE

2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

75 days.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table. Antibiotic discs are then promptly placed onto the agar surface.

Inoculum: 0.5 Mc Farland.

Incubation conditions: 5% CO₂ at 35 ± 2°C for 20-24 h.

QC Table.

Microorganism	Specification
<i>Streptococcus pneumoniae</i> ATCC® 49619	Zone diameters within the acceptable QC ranges indicated by the CLSI
<i>Escherichia coli</i> ATCC® 25922	Zone diameters within the acceptable QC ranges indicated by the CLSI

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. CLSI M02-A12 (2015) Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard - Twelfth Edition. Wayne, PA: Clinical and Laboratory Standards Institute.
2. CLSI M31-A3 (2008) Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated from Animals; Approved Standard - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute.
3. CLSI M6-A2 (2006) Protocols for Evaluating Dehydrated Mueller-Hinton Agar; Approved Standard - Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute.
4. Barry, Garcia and Thrupp (1970) Am. J. Clin. Pathol. 53:149.
5. Ryan, Schoenknecht and Kirby (1970) Hospital Practice 5:91.
6. Bauer, Kirby, Sherris and Turck (1966) Am. J. Clin. Pathol. 45:493.

PRESENTATION	Contents	Ref.
Mueller Hinton II Agar (Sheep Blood 5%) 90 mm ready-to-use plates	20 plates	10131
Mueller Hinton II Agar (Sheep Blood 5%) 90 mm ready-to-use plates	100 plates	10131*
Mueller Hinton II Agar (Sheep Blood 5%) 140 mm ready-to-use plates	10 plates	11231
Mueller Hinton II Agar (Sheep Blood 5%) 120 mm square ready-to-use plates	5 plates	12032

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care	 Keep away from sunlight
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse	

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