Brucella Blood Agar w Hemin and Vitamin K1

For isolation of obligate anaerobes and MIC determination.

INTENDED PURPOSE

Medium for the isolation and subculture of anaerobes from clinical specimens. This medium is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Brucella Blood Agar w Hemin and Vitamin K1 is a modification of Brucella Agar that has been supplemented with hemin and vitamin K1 to support the growth of obligate anaerobes, especially *Bacteroides*, *Prevotella*, and *Porphyromonas* when incubated anaerobically. It is also used with the gradient test method for the determination of the MIC (Minimum Inhibitory Concentration) for anaerobic microorganisms.

TYPICAL FORMULA* Per Litre of Purified Water

Pancreatic Digest of Casein	10.0 g
Peptic Digest of Animal Tissue	10.0 g
Glucose	1.0 g
Yeast Extract	2.0 g
Sodium chloride	5.0 g
Sodium Bisulfite	0.1 g
Vitamin K1	0.001 g
Hemin	0.005 g
Sheep Blood, defibrinated	50.0 ml
Agar	15.0 g
Final pH 7.2 ± 0.2 at 25°C	

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Pancreatic digest of casein and peptic digest of animal tissue provide amino acids, nitrogen, carbon, minerals and vitamins essential for bacterial growth. Glucose is the fermentable carbohydrate. Yeast extract is a source of vitamins, particularly of group B. Sodium chloride maintains the osmotic balance of the medium. Sodium bisulfite lowers the redox potential to a range suitable for strict anaerobes. Hemin and vitamin K1 support the growth of certain fastidious anaerobes. Sheep blood provides additional nutrients and allows to detect hemolytic reactions. Agar is the solidifying agent.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: Inoculating loop, quality control organisms, incubator, controlled atmosphere generators, jars. **Other Reagents (Optional):** Physiological solution (0.85% saline), MTS™ (MIC Test Strip).

SPECIMENS

Brucella Blood Agar w Hemin and Vitamin K1 can be used for the isolation and cultivation of strict anaerobes from all types of specimens.

For antimicrobial susceptibility testing (AST) with gradient tests, the method requires the use of pure cultures. Refer to specific guidelines for information about collection, transport and preparation of anaerobic specimens.

TEST PROCEDURE

Isolation

Inoculate the medium by directly streaking the specimen on the agar surface to obtain well-isolated colonies. Incubate the plates in an inverted position at 35 ± 2°C for 24-72 hours in an anaerobic atmosphere.

Note: Alternative incubation temperature and time may be required depending on the method used, type of specimen and the microorganism being tested. The user is responsible for choosing the appropriate incubation conditions.

Once the primary cultures are set up, it is recommended to also run the test using a selective medium (such as Schaedler Kanamycin-Vancomycin Agar with 5% Sheep Blood) and an aerobic medium (such as Columbia Agar with 5% Sheep Blood). This last plate is incubated aerobically enriched with carbon dioxide together with the anaerobic cultures.

For MIC testing using the MTS™ method, ensure the Brucella Blood Agar w Hemin and Vitamin K1 plates are dried prior to inoculation to avoid excess moisture, which may lead to incorrect results. Also, note that the medium must be inoculated with pure calibrated strains obtained by culture on a Petri plate.

For more detailed information, see the Application Guide as well as other supporting material available from liofilchem.com/MTS.

INTERPRETING RESULTS

After incubation, the plates are inspected for growth, also in comparison to the growth on the other media. Colonies which appear on this medium are suspected to be strict anaerobes if they do not grow on aerobically incubated blood agar plates.

Note: further testing should be conducted to confirm the presumptive identification of organisms isolated on this medium.

For how to read and interpret the MTS results consult appropriate guidance.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

Avoid quick temperature shifts to prevent condensation.

SHELF LIFE

2 months.

QUALITY CONTROL

Appearance of Medium: Cherry red opaque.

Control strain		Inoculum	Incubation	Specification
<i>Bacteroides fragilis</i>	ATCC® 25285	50-100 CFU	24-72 h/ 35 ± 2°C	Good growth, grey colonies
<i>Fusobacterium nucleatum</i>	ATCC® 25586			Good growth, grey-white colonies surrounded by dark-grey zones
<i>Peptostreptococcus anaerobius</i>	ATCC® 27337			Good growth, whitish colonies
<i>Clostridium perfringens</i>	ATCC® 13124			Good growth, grey-white colonies, beta hemolysis

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of Brucella Blood Agar w Hemin and Vitamin K1 was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Invalid results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology. The operator should understand the principles of the procedures, including its performance limitations, in advance of operation to avoid potential mistakes.

Growth depends on the requirements of each individual microorganism. It is therefore possible that certain strains which have specific requirements (substrate, temperature, etc.) may not develop.

This medium is not specifically selective for strict anaerobes; facultative organisms will also grow. Therefore, it is important to compare the result of the anaerobic culture with that of an aerobically incubated plate if mixed cultures are obtained.

All presumptive anaerobic organisms must be identified by confirmatory test.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

See ordering info below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Brucella Blood Agar w Hemin and Vitamin K1	Plate 90 mm	20 (2 x 10) plates	10245

Revision History

Revision	Release Date	Change Summary
0	2023-07-06	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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liofilchem.com/ifu-sds

References / Bibliografia / Bibliografia

1. CLSI. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th ed. CLSI standard M11. Wayne PA: Clinical and Laboratory Standards Institute 2018.
2. Performance Standards for Antimicrobial Susceptibility Testing; Twenty First Informational Supplement - CLSI M100-S21 Third edition Jan. 2011 ISBN 1-56238-742-1.
3. Citron DM et al 1991. Evaluation of the Etest for susceptibility testing of anaerobic bacteria. J. Clin. Microbiol. 29: 2197–2203.
4. BARON EJ. and FINEGOLD SM. - Bailey and Scott's Diagnostic Microbiology - 8th ed - Mosby Co, 1990.
5. MacFaddin JF 1985. Media for isolation cultivation-identification-maintenance of medical bacteria, vol. I. Williams & Wilkins, Baltimore, Maryland USA.
6. Dowell 1975 Am. J. Med. Technol. 41:402
7. Finegold Miller and Posnick 1965 Hernahrungsforschung 10:517.
8. Gibbons RJ and JB MacDonald 1960. Hemin and vitamin K compounds as required factors for the cultivation of certain strains of Bacteroides melaninogenicus J. Bacteriol. 80:164–170.

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