## Middlebrook 7H11 Agar

Selective medium for the isolation and cultivation of mycobacteria.

### DESCRIPTION

Middlebrook 7H11 Agar is a medium used for the selective isolation of *Mycobacterium* spp, especially *Mycobacterium tuberculosis*, from clinical specimens as well as for the cultivation of pure cultures of mycobacteria.

This agar-based medium provides improved susceptibility testing compared to media based on coagulated egg (Lowenstein formulations) where drugs may be more easily inactivated.

| TYPICAL FORMULA             | (g/l)   |
|-----------------------------|---------|
| Pancreatic Digest of Casein | 1.0     |
| L-Glutamic Acid             | 0.5     |
| Sodium Citrate              | 0.4     |
| Pyridoxine Hydrochloride    | 0.001   |
| Biotin                      | 0.0005  |
| Ferric Ammonium Citrate     | 0.04    |
| Ammonium Sulfate            | 0.5     |
| Disodium Phosphate          | 1.5     |
| Monopotassium Phosphate     | 1.5     |
| Magnesium Sulfate           | 0.05    |
| Malachite Green             | 0.001   |
| Agar                        | 15.0    |
| Glycerol                    | 5 ml    |
| Sodium Chloride             | 0.85    |
| Glucose                     | 2.0     |
| Bovine Albumin (Fraction V) | 5.0     |
| Catalase                    | 0.003   |
| Oleic Acid                  | 0.06 ml |
| Final pH 6.6 ± 0.2 at 25°C  |         |

### METHOD PRINCIPLE

Pancreatic digest of casein is the ingredient added to Middlebrook 7H10 Agar (ref. 30368) to create Middlebrook 7H11 Agar providing nitrogen, vitamins and amino acids to enhance the growth of fastidious strains of *M. tuberculosis*. Glutamic acid, sodium citrate, pyridoxine, biotin and ammonium sulfate supply growth factors. Ferric ammonium citrate and magnesium sulfate are sources of trace ions. Phosphates help maintaining the pH of the medium. Malachite green is the selective agent inhibiting the contaminant microbial flora. Malachite green serves as pH indicator as well. Agar is the solidifying agent. Glycerol and glucose are energy sources. Sodium chloride maintains the osmotic equilibrium. Albumin protects the tubercle bacilli against toxic agents. Catalase destroys toxic peroxides that may be present in the medium. Oleic acid is a fatty acid utilized in the mechanism of mycobacteria.

### TEST PROCEDURE

Inoculate the medium according to test procedures recommended by the Centers for Disease Control and Prevention (CDC) or consult appropriate references. Keep plates and tubes shielded from light and incubate at 35 ± 2°C for up to 8 weeks in aerobic atmosphere enriched with 5-10% carbon dioxide.

### INTERPRETING RESULTS

Examine weekly for growth, pigment production and colony morphology. Carry out identification tests according to established laboratory procedures.

N.B. Negative culture results do not rule-out active infection by mycobacteria. Some factors that may lead to unsuccessful cultures are:

- The specimens was not representative of the infectious materials; i.e. saliva instead of sputum.
- The mycobacteria were destroyed during digestion and decontamination of the specimen.
- Gross contamination interfered with the growth of mycobacteria.
- Proper aerobic conditions and increased CO<sub>2</sub> tension were not provided during incubation.

**APPEARANCE**

Slightly opalescent, light yellowish green.

**STORAGE**

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

**SHELF LIFE**

Ready-to-use plates: 6 months.

Medium in tubes: 1 year.

**QUALITY CONTROL**

The medium is inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 100-300 CFU.

Inoculum for selectivity: 10<sup>3</sup>-10<sup>4</sup> CFU.

Incubation conditions: aerobic atmosphere with 5-10%CO<sub>2</sub> at 35 ± 2°C for up to 21 days.

**QC Table.**

| Microorganism                       |             | Growth                            |
|-------------------------------------|-------------|-----------------------------------|
| <i>Mycobacterium intracellulare</i> | ATCC® 13950 | Good                              |
| <i>Mycobacterium scrofulaceum</i>   | ATCC® 19981 | Good                              |
| <i>Escherichia coli</i>             | ATCC® 25922 | Partially to completely inhibited |

**WARNING AND PRECAUTIONS**

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

**DISPOSAL OF WASTE**

Disposal of waste must be carried out according to national and local regulations in force.


**BIBLIOGRAPHY**

- Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfaller (2007) Manual of Clinical Microbiology. 9<sup>th</sup> ed. ASM Press, Washington, D.C.
- Clinical and Laboratory Standards Institute (2004) Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard, 3<sup>rd</sup> ed. M22-A3. CLSI, Wayne, PA
- Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. (1997) Color Atlas and Textbook of Diagnostic Microbiology. 5<sup>th</sup> ed. Lippincott Williams and Wilkins, Philadelphia, PA.
- Kent, P.T. and G.P. Kubica (1985) Public Health Mycobacteriology, A Guide for the Level III Laboratory. Dept. of H.H.S. and Centers for Disease Control (CDC), Atlanta, GA.
- Middlebrook, G., M.L. Cohn, W.E. Dye, Russell, Jr., and D. Levy (1960) Microbiologic procedures of value in tuberculosis. Acta. Tuberc. Scand. 38:66-81.
- Middlebrook, G., and M.L. Cohn (1958) Bacteriology of tuberculosis: laboratory methods. Am. J. Public Health. 48:844-853.
- Middlebrook, G., M.L. Cohn, and W.B. Scheffer (1954) Studies on isoniazid and tubercle bacilli. III. The isolation, drug-susceptibility and catalase testing of tubercle bacilli from isoniazid-treated patients. Am. Rev. Tuberc. 70:852-872.

**PRESENTATION**

|                      |                           | Contents          | Ref.  |
|----------------------|---------------------------|-------------------|-------|
| Middlebrook 7H11Agar | 90 mm Ready-to-use plates | 20 plates         | 10416 |
| Middlebrook 7H11Agar | Slant tubes               | 10 x 8.5 ml tubes | 37000 |

**TABLE OF SYMBOLS**

|                             |  |   |  |   |
|-----------------------------|--|---|--|---|
| <b>LOT</b> Batch code       | <b>IVD</b> <i>In vitro</i> Diagnostic Medical Device   |  Manufacturer                      |  Use by                               |  Fragile, handle with care |
| <b>REF</b> Catalogue number |  Temperature limitation |  Contains sufficient for <n> tests |  Caution, consult Instruction For Use |  Do not reuse              |

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