Gelatin Peptone Agar

Medium for the cultivation of non-fastidious organisms.

DESCRIPTION

Gelatin Peptone Agar is a general purpose medium used for supporting the growth of a wide range of organisms providing no advantage to any specific group.

TYPICAL FORMULA

	(g/l)
Gelatin Peptone	5.0
Agar	15.0

Final pH 7.0 ± 0.2 at 25°C

METHOD PRINCIPLE

Gelatin peptone (pancreatic digest of gelatin) provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Gelatin as the only source of nutrients having also low cystine and tryptophan content makes the medium absolutely nonselective. Agar is the solidifying agent.

PREPARATION

Dehydrated medium Suspend 20 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.

Medium in tubes Melt the content of the tube in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the tube upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into a Petri dish.

TEST PROCEDURE

Inoculate the medium by pour plating or spread plating technique.

Incubation conditions may vary depending on the organisms under study. For a general aerobic count, incubate aerobically at 30°C for 72 hours.

INTERPRETING RESULTS

Observe colony growth.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, light beige.

Prepared medium: clear to slightly opalescent, light amber.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store tubes and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Medium in tubes: 2 years.

Ready-to-use plates: 6 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Incubation conditions: aerobically at $30 \pm 1^\circ\text{C}$ for 72 ± 3 hours.

QC Table.

Microorganism		Growth
<i>Bacillus subtilis</i>	ATCC® 6633	Good
<i>Escherichia coli</i>	ATCC® 8739	Good
<i>Staphylococcus aureus</i>	ATCC® 6538	Good

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

- Eaton, A.D., L.S. Clesceri, and A.E. Greenberg (1995) Standards methods for the examination of water and wastewater, 19th ed.
- Association of Official Analytic Chemist (1995) Bacteriological analytic manual, 8th ed.

PRESENTATION		Contents	Ref.
Gelatin Peptone Agar	90 mm ready-to-use plates	20 plates	10442
Gelatin Peptone Agar	Tubes	20 x 22 ml tubes	31094
Gelatin Peptone Agar	Tubes	10 x 22 ml tubes	34094
Gelatin Peptone Agar	Dehydrated medium	500 g of powder	610054
Gelatin Peptone Agar	Dehydrated medium	100 g of powder	620054
Gelatin Peptone Agar	Dehydrated medium	5 kg of powder	6100545

TABLE OF SYMBOLS

LOT Batch code	 Keep away from sunlight	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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