Middlebrook 7H10 Agar

Selective medium for the isolation and cultivation of mycobacteria.

DESCRIPTION

Middlebrook 7H10 Agar is a medium used for the selective isolation of *Mycobacterium* spp, especially *Mycobacterium tuberculosis*, from clinical specimens as well as for the cultivation of pure cultures of mycobacteria.

This agar-based medium provides improved susceptibility testing compared to media based on coagulated egg (Lowenstein formulations) where drugs may be more easily inactivated.

TYPICAL FORMULA	(g/l)
L-Glutamic Acid	0.5
Sodium Citrate	0.4
Pyridoxine Hydrochloride	0.001
Biotin	0.0005
Ferric Ammonium Citrate	0.04
Ammonium Sulfate	0.5
Disodium Phosphate	1.5
Monopotassium Phosphate	1.5
Magnesium Sulfate	0.025
Malachite Green	0.00025
Agar	15.0
Calcium Chloride	0.0005
Zinc Sulfate	0.001
Copper Sulfate	0.001
Glycerol	5 ml
Sodium Chloride	0.85
Glucose	2.0
Bovine Albumin (Fraction V)	5.0
Catalase	0.003
Oleic Acid	0.06 ml

Final pH 6.6 ± 0.2 at 25°C

METHOD PRINCIPLE

Glutamic acid, sodium citrate, pyridoxine, biotin and ammonium sulfate supply growth factors. Ferric ammonium citrate, magnesium sulfate, calcium chloride, zinc sulfate and copper sulfate are sources of trace ions. Phosphates help maintaining the pH of the medium. Malachite green is the selective agent inhibiting the contaminant microbial flora. Malachite green serves as pH indicator as well. Agar is the solidifying agent. Glycerol and glucose are energy sources. Sodium chloride maintains the osmotic equilibrium. Albumin protects the tubercle bacilli against toxic agents. Catalase destroys toxic peroxides that may be present in the medium. Oleic acid is a fatty acid utilized in the mechanism of mycobacteria.

TEST PROCEDURE

Inoculate the medium according to test procedures recommended by the Centers for Disease Control and Prevention (CDC) or consult appropriate references. Keep plates and tubes shielded from light and incubate at 35 ± 2°C for up to 8 weeks in aerobic atmosphere enriched with 5-10% carbon dioxide.

INTERPRETING RESULTS

Examine weekly for growth, pigment production and colony morphology. Carry out identification tests according to established laboratory procedures.

N.B. Negative culture results do not rule-out active infection by mycobacteria. Some factors that may lead to unsuccessful cultures are:

- The specimens was not representative of the infectious materials; i.e. saliva instead of sputum.
- The mycobacteria were destroyed during digestion and decontamination of the specimen.
- Gross contamination interfered with the growth of mycobacteria.
- Proper aerobic conditions and increased CO₂ tension were not provided during incubation.

APPEARANCE

Slightly opalescent, light yellowish green.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Ready-to-use plates: 6 months.

Medium in tubes: 1 year.

QUALITY CONTROL

The medium is inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 100-300 CFU.

Inoculum for selectivity: 10³-10⁴ CFU.

Incubation conditions: aerobic atmosphere with 5-10%CO₂ at 35 ± 2°C for up to 21 days.

QC Table.

Microorganism		Growth
<i>Mycobacterium intracellulare</i>	ATCC® 13950	Good
<i>Mycobacterium scrofulaceum</i>	ATCC® 19981	Good
<i>Escherichia coli</i>	ATCC® 25922	Partially to completely inhibited

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

- Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfaller (2007) Manual of Clinical Microbiology. 9th ed. ASM Press, Washington, D.C.
- Clinical and Laboratory Standards Institute (2004) Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard, 3rd ed. M22-A3. CLSI, Wayne, PA
- Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. (1997) Color Atlas and Textbook of Diagnostic Microbiology. 5th ed. Lippincott Williams and Wilkins, Philadelphia, PA.
- Kent, P.T. and G.P. Kubica (1985) Public Health Mycobacteriology, A Guide for the Level III Laboratory. Dept. of H.H.S. and Centers for Disease Control (CDC), Atlanta, GA.
- Middlebrook, G., M.L. Cohn, W.E. Dye, Russell, Jr., and D. Levy (1960) Microbiologic procedures of value in tuberculosis. Acta. Tuberc. Scand. 38:66-81.
- Middlebrook, G., and M.L. Cohn (1958) Bacteriology of tuberculosis: laboratory methods. Am. J. Public Health. 48:844-853.
- Middlebrook, G., M.L. Cohn, and W.B. Scheffer (1954) Studies on isoniazid and tubercle bacilli. III. The isolation, drug-susceptibility and catalase testing of tubercle bacilli from isoniazid-treated patients. Am. Rev. Tuberc. 70:852-872.

PRESENTATION

		Contents	Ref.
Middlebrook 7H10 Agar	90 mm Ready-to-use plates	20 plates	10453
Middlebrook 7H10 Agar	Slant tubes	10 x 8.5 ml tubes	30368

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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