Chromatic™ MRSA

Chromogenic selective medium for the isolation of methicillin-resistant *Staphylococcus aureus*.

DESCRIPTION

Chromatic™ MRSA is a chromogenic selective medium used for the isolation of methicillin/oxacillin resistant *S. aureus*.

MRSA infections are endemic in both hospitals and the community in all parts of the world. The mortality of MRSA bloodstream infections is double that of similar infections caused by methicillin susceptible strains due to delayed adequate treatment and inferior alternative treatment regimens.

TYPICAL FORMULA

| | (g/l) |
|--------------------------------|-------|
| Peptone and Yeast Extract | 30.0 |
| Sodium Chloride | 10.0 |
| Sodium Phosphate Dibasic | 2.5 |
| Selective and Matting Agents | 16.5 |
| Chromogenic and Antibiotic Mix | 0.8 |
| Agar | 15.0 |

Final pH 6.9 ± 0.2 at 25°C

METHOD PRINCIPLE

Peptone and yeast extract supply amino acids, nitrogen, carbon, minerals, vitamins and other nutrients which support the growth of microorganism. Sodium chloride enhances growth of *S. aureus*. Sodium phosphate is the buffer. Selective agents inhibit the growth of yeast and the most of Gram-negative and Gram-positive bacteria other than methicillin-resistant staphylococci. Matting agents enhance colonies contrast on the medium. Chromogenic and antibiotic mix allows the optimal recovery of MRSA and the identification based on a mauve or orange-mauve coloration of the colonies. Agar is the solidifying agent.

TEST PROCEDURE

Inoculate the plates by streaking directly the specimen onto the agar surface. Incubate aerobically at 37°C for 18-24 hours.

INTERPRETING RESULTS

After incubation observe the color of the colonies and interpret the results as indicated in the ID table.

ID Table.

| Microorganism | Typical colony color |
|--|-----------------------|
| Methicillin/oxacillin resistant <i>S. aureus</i> | Mauve to orange-mauve |
| Other Gram-positive bacteria (if not inhibited) | White or blue |

See pictures in Appendix I.

APPEARANCE

Opaque, whitish.

STORAGE

2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

4 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC tables.

Inoculum for productivity: 10-100 CFU

Inoculum for selectivity: 10^4 - 10^5 CFU.

Inoculum for specificity: $\leq 10^4$ CFU.

Incubation conditions: aerobically at $35 \pm 2^\circ\text{C}$ for 18-24 h.

QC Table.

| Microorganism | | Growth | Specification |
|-------------------------------------|-------------|-----------|----------------|
| <i>Staphylococcus aureus</i> (MRSA) | ATCC® 43300 | Good | Mauve colonies |
| <i>Staphylococcus aureus</i> (MSSA) | ATCC® 25923 | Inhibited | --- |
| <i>Staphylococcus aureus</i> (MSSA) | ATCC® 6538 | Inhibited | --- |
| <i>Escherichia coli</i> | ATCC® 25922 | Inhibited | --- |
| <i>Proteus mirabilis</i> | ATCC® 25933 | Inhibited | --- |
| <i>Pseudomonas aeruginosa</i> | ATCC® 27853 | Inhibited | --- |

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. Evaluation of CHROMagar Staph aureus, a new chromogenic medium, for isolation and presumptive identification of *Staphylococcus aureus* from human clinical specimens. Gaillot O. et al. 2001. Journal of Clinical Microbiology, 38 : 1587-1591.
2. Dépistage nasal de *Staphylococcus aureus*. Nécessité de standardiser les protocoles. Laudat P. et al. 2000 Poster 343/P2 presented at RICAI in Paris (France).

PRESENTATION

| | | Contents | Ref. |
|----------------------|---------------------------|-----------|-------|
| Chromatic™ MRSA Agar | 90 mm ready-to-use plates | 20 plates | 10599 |

TABLE OF SYMBOLS

| | | | | | |
|-----------------------------|--|---|--|---|---|
| LOT Batch code | IVD <i>In vitro</i> Diagnostic Medical Device |  Manufacturer |  Use by |  Fragile, handle with care |  Keep away from sunlight |
| REF Catalogue number |  Temperature limitation |  Contains sufficient for <n> tests |  Caution, consult Instruction For Use |  Do not reuse | |

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