MacConkey Agar No.2

For detection of enterococci from clinical and nonclinical specimens with mixed microbial flora.

INTENDED PURPOSE

Selective medium used for the identification of enterococci in the presence of non-lactose fermenters and coliform bacteria in water, food and clinical specimens. This medium is intended as an aid in the diagnosis, requiring additional tests to complete the diagnostic results.

DESCRIPTION

MacConkey Agar No.2 is a modification of MacConkey Agar (also known as MacConkey Agar No.3).

This medium contains Bile Salts No.2 which are less inhibitory allowing the growth and recognition of enterococci. Such bacteria are frequently sought as an index of faecal pollution.

TYPICAL FORMULA*	(g/litre)
Peptone	20.0
Lactose	10.0
Sodium Chloride	5.0
Bile Salts No.2	1.5
Crystal Violet	0.001
Neutral Red	0.05
Agar	13.5
Final pH 7.2 ± 0.2 at 25°C	

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptone provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Lactose is the fermentable carbohydrate. Sodium Chloride maintains the osmotic balance of the medium. Bile salts and crystal violet are the selective agents, inhibiting Gram-positive organisms other than enterococci and allowing Gram-negative bacteria to grow. Neutral red is the pH indicator. Agar is the solidifying agent.

PREPARATION

Dehydrated medium

Suspend 50.0 g of the powder in 1 litre of distilled or deionized water. Heat to boil for 1 minute shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: Autoclave, sterile Petri plates, test tubes, inoculating loops, swabs, incubator, quality control organisms.

SPECIMENS

Clinical specimens, such as urine and faeces, should be sampled at the acute stage, before antimicrobial therapy (where possible) and examined as soon as possible after collection. Good laboratory practices for collection, transport and storage of the clinical specimens should be applied. Refer to specific guidelines for more information about specimen collection and preparation.

TEST PROCEDURE

Ensure there is no visible moisture on the plates before use.

Inoculate the plates by directly streaking the specimen onto the agar surface in order to obtain well-isolated colonies.

Incubate aerobically at 35 ± 2°C for 18-24 h. It may be necessary to extend incubation for further 24 h.

For more details, consult appropriate guidances.

INTERPRETING RESULTS

Colourless or pink to red colonies are produced on this medium depending upon the ability of the isolate to ferment lactose.

Enterococci appear as small, intensely coloured red colonies surrounded by a zone of pale red precipitate. Lactose-fermenting bacteria, such as *E. coli* and *Klebsiella* spp., grow as pink-red colonies with or without a zone of bile precipitate.

Non-lactose fermenting organisms, such as *Salmonella* spp, *Shigella* and *Proteus* spp, form colourless colonies, whilst staphylococci and non-faecal streptococci, are partially or completely inhibited.

Note: Further testing should be conducted to confirm the presumptive identification of organisms isolated on this medium.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Ready-to-use plates: 6 months.

QUALITY CONTROL

Appearance of Dehydrated Medium: Free-flowing, homogeneous, beige-pink.

Appearance of Prepared Medium: Slightly opalescent, red-violet.

Expected Cultural Response:

Control strain	Inoculum	Incubation	Criteria	Specification	
<i>Enterococcus faecalis</i>	WDCM 00087 (ATCC® 29212; NCTC 12697)	50-100 CFU	18-24 h / 35 ± 2°C	Good growth (P _R ≥ 0.5)	Red colonies
<i>Salmonella</i> Enteritidis	WDCM 00030 (ATCC® 13076; NCTC 12694)				Colourless colonies
<i>Escherichia coli</i>	WDCM 00013 (ATCC® 25922; NCTC 12241)				Pink-red colonies (bile precipitate)
<i>Staphylococcus aureus</i>	WDCM 00034 (ATCC® 25923; NCTC 12981)	10 ⁴ -10 ⁶ CFU		Inhibition	—

A productivity ratio (P_R) of 0.5 is equivalent to a recovery rate of 50%.

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of MacConkey Agar No.2 was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Invalid results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology. The operator should understand the principles of the procedures, including its performance limitations, in advance of operation to avoid potential mistakes.

Due to nutritional variation, some strains may result in poor growth or fail to grow on this medium. A nonselective medium, such as Columbia Agar with 5% sheep blood, should also be inoculated with the sample for primary isolation.

MacConkey Agar No.2 is intended as an aid in the diagnosis of infectious diseases, requiring further tests to complete the diagnostic results. All identification tests should ideally be performed from non-selective agar.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.

- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
MacConkey Agar No.2	Plate 90 mm	20 plates	10603
	Dehydrated media	500 g	610057

Significant changes from previous version:

Document	Release Date	Change Summary
10603_IFU-0	2022-10-26	Document creation (layout and content in compliance with IVDR 2017/746)

This IFU document and the SDS are available from the online Support Center:

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
liofilchem.com/ifu-sds

References / Riferimenti

1. Public Health England (2021) UK SMI ID 4 issue 4: Identification of *Streptococcus* species, *Enterococcus* species and morphologically similar organisms. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/1019677/ID_4i4.pdf
2. EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
3. Public Health England (2019) UK SMI B 41 (issue 8.7): investigation of urine. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/770688/B_41i8.7.pdf
4. Public Health England (2017) Inoculation of culture media for bacteriology. UK Standards for Microbiology Investigations. Q 5 Issue 2. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/583859/Q_5i2.pdf
5. Murray, Baron, Jorgensen, Landry and Pfaller ed. (2007) Manual of clinical microbiology, 9th ed. American Society for Microbiology, Washington, D.C.
6. Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). 2003. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C. USA.
7. Farmer III, J.J. (2003) Enterobacteriaceae: introduction and identification. In: Murray, P. R., E. J. Baron, J.H. Jorgensen,
8. Downes, F.P., and K. Ito. (2001) Compendium of methods for the microbiological examination of foods. 4th edition. American Public Health Association (APHA). Washington, D.C. USA.
9. Mazura-Reetz, G., T.R. Neblett, and J.M. Galperin. 1979. MacConkey agar: CO₂ vs. ambient incubation, abstr. C 179, p. 339. Abstr. 79th Annu. Meet. Am. Soc. Microbiol. 1979.
10. MacConkey, A. (1905) Lactose-fermenting bacteria in faeces. J. Hyg. 5:333–379.
11. MacConkey, A.T. (1900) Note on a new medium for the growth and differentiation of the *Bacillus coli* communis and the *Bacillus typhi* abdominalis. The Lancet, Part II:20.

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