## Dermatophyte Agar (DTM)

Selective medium with chloramphenicol for detection of pathogenic dermatophytic fungi from clinical and veterinary specimens.

### DESCRIPTION

Dermatophyte Agar (DTM) is a selective medium used for the isolation and differentiation of dermatophytes from clinical and veterinary specimens including those contaminated with other microorganisms.

### TYPICAL FORMULA

	(g/l)
Enzymatic Digest of Soybean Meal	10.0
Glucose	10.0
Phenol Red	0.2
Cycloheximide	0.5
Gentamicin	0.1
Chloramphenicol	0.1
Agar	20.0
Final pH 5.5 ± 0.2 at 25°C	

### METHOD PRINCIPLE

Enzymatic digest of soybean meal provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Glucose is an energy source. Phenol red is the pH indicator. Cycloheximide inhibits most saprophytic moulds. Gentamicin suppresses Gram-negative bacteria including *Pseudomonas* spp. Chloramphenicol is included to replace chlortetracycline used in the original formulation made by Taplin in 1969. It is a broad-spectrum antibiotic, inhibiting a wide range of Gram-positive and Gram-negative bacteria. Agar is the solidifying agent.

### PREPARATION

**Medium in bottles** Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into final containers.

### TEST PROCEDURE

Inoculate either plates or slant tubes by streaking directly the sample onto the agar surface. Streak the specimen as soon as possible after it is received in the laboratory. A control medium, Sabouraud Dextrose Agar (ref. 10035), may be inoculated in parallel. Incubate media aerobically at room temperature (15-30°C) for up to 14 days.

### INTERPRETING RESULTS

Examine media daily and observe for development of a red color change in the medium. Pathogenic dermatophytes (*Epidermophyton*, *Microsporum*, and *Trichophyton* spp) typically produce alkaline metabolites and most of them will produce a color change in 3-6 days.

Growth, without a color change to red, indicates that the organism is probably not a dermatophyte. Further biochemical and/or serological testing is recommended for complete identification.

If growth appears on the control medium (Sabouraud Dextrose Agar) and no growth appears on DTM, the organism is not a dermatophyte. Colonies with green or black hyphae are not typical of dermatophytes even though the media may turn red.

Transfer of growth from slants to plated media may be required in order to obtain pure cultures of fungi.

### APPEARANCE

Clear, yellowish-orange.

### STORAGE

Store at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

**SHELF LIFE**

Medium in bottles: 2 years.

Slant tubes: 1 year.

Ready-to-use plates: 6 months.

**QUALITY CONTROL**

The medium is inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Inoculum for selectivity: 10<sup>4</sup>-10<sup>6</sup> CFU.

Incubation conditions: aerobically at 25 ± 2°C for 2-7 days.

**QC Table.**

Microorganism	ATCC®	Specification
<i>Trichophyton mentagrophytes</i>	ATCC® 9533	Good growth, red medium
<i>Candida albicans</i>	ATCC® 10231	Good growth, no color change in the medium
<i>Aspergillus brasiliensis</i>	ATCC® 16404	Inhibition
<i>Staphylococcus aureus</i>	ATCC® 25923	Inhibition
<i>Escherichia coli</i>	ATCC® 25922	Inhibition

**WARNING AND PRECAUTIONS**

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

**DISPOSAL OF WASTE**

Disposal of waste must be carried out according to national and local regulations in force.

**BIBLIOGRAPHY**

- Larone (1995) Medically important fungi: a guide to identification, 3rd ed. American Society for Microbiology, Washington, D.C.
- MacFaddin, J.D. (1985) Media for isolation-cultivation-identification-identification-maintenance of medical bacteria, p. 695-699, vol. 1. Williams & Wilkins, Baltimore, MD.
- Taplin D., N. Zaias, N. Rebell and H. Blank (1969) Isolation and recognition of dermatophytes on a new medium (DTM). Arch Dermatol. 99:203.
- Sabouraud (1892) Ann. Dermatol. Syphil. 3:1061.

PRESENTATION		Contents	Ref.
Dermatophyte Agar (DTM)	90 mm ready-to-use plates	20 plates	11052
Dermatophyte Agar (DTM)	90 mm ready-to-use plates	100 plates	11052*
Dermatophyte Agar (DTM)	140 mm ready-to-use plates	10 plates	10244
Dermatophyte Agar (DTM)	Slant tubes	10 x 7 ml tubes	33086
Dermatophyte Agar (DTM)	Slant tubes	20 x 7 ml tubes	31086
Dermatophyte Agar (DTM)	Bottles	6 x 100 ml bottles	402200

**TABLE OF SYMBOLS**

<b>LOT</b> Batch code	<b>IVD</b> <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
<b>REF</b> Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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