Enterococco Agar

Selective medium for detection of group D streptococci from clinical and nonclinical specimens.

INTENDED PURPOSE

Enterococco Agar is a selective medium used for the isolation and enumeration of faecal streptococci (group D) from clinical and nonclinical specimens. This medium is intended as an aid in the diagnosis, requiring further identification tests to complete the diagnostic results.

DESCRIPTION

This medium is based on a standard formulation (Bile Aesculin Azide Agar) outlined in ISO 7899-2 for confirmation and enumeration of intestinal enterococci in water. With a higher concentration of sodium azide Enterococco Agar is suitable for isolation of pure cultures from clinical specimens containing mixed flora.

TYPICAL FORMULA*	(g/litre)
Tryptone	17.0
Peptone	3.0
Yeast Extract	5.0
Sodium Chloride	5.0
Ox Bile, dehydrated	10.0
Aesculin	1.0
Ferric Ammonium Citrate	0.5
Sodium Azide	0.25
Agar	15.0
Final pH 7.1 ± 0.2 at 25°C	

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Tryptone and peptone provide amino acids, carbon, nitrogen, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium. Bile inhibits Gram-positive bacteria other than enterococci. The glycoside aesculin is hydrolyzed from enterococci to give aesculetin and glucose. Aesculetin reacts with iron ions into the medium to form a dark brown or black complex. Sodium azide suppresses the accompanying Gram-negative microbial flora. Agar is the solidifying agent.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: Autoclave, water bath, sterile Petri plates, test tubes, inoculating loops, swabs, incubator, quality control organisms.

SPECIMENS

All types of clinical specimens. Specimens should be sampled at the acute stage, before antimicrobial therapy (where possible) and examined as soon as possible after collection. Good laboratory practices for collection, transport and storage of the clinical specimens should be applied. Refer to specific guidelines for more information about specimen collection and preparation.

TEST PROCEDURE

Inoculate the plates by directly streaking the specimen on the agar surface or spread the sample from an enrichment culture to obtain well-isolated colonies.

Incubate aerobically at 35 ± 2°C for 24-48 h.

For more details, consult appropriate guidances.

INTERPRETING RESULTS

Regard all typical colonies showing a brown to black colouration in the surrounding medium as giving a positive aesculin reaction and count as group D streptococci (*Enterococcus* spp. and *Streptococcus bovis*).

Note: Further testing should be conducted to confirm the presumptive identification of organisms isolated on this medium.

STORAGE

Store plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

6 months.

QUALITY CONTROL

Appearance of the Medium: Slightly opalescent, dark amber.

Expected Cultural Response:

Control strain	Inoculum	Incubation	Criteria	Specification	
<i>Enterococcus faecalis</i>	WDCM 00087 (ATCC® 29212; NCTC 12697)	50-100 CFU	24-48 h / 35 ± 2°C	Good growth (P _R ≥ 0.5)	Blackening around colonies
<i>Enterococcus faecalis</i>	WDCM 00009 (ATCC® 19433; NCTC 775)				
<i>Streptococcus pyogenes</i>	ATCC® 19615 (NCTC 12696)	10 ⁴ -10 ⁶ CFU	Inhibition (partial to complete)	Colourless colonies	
<i>Escherichia coli</i>	WDCM 00013 (ATCC® 25922; NCTC 12241)				

A productivity ratio (P_R) of 0.5 is equivalent to a recovery rate of 50%.

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of Enterococco Agar was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Invalid results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology. The operator should understand the principles of the procedures, including its performance limitations, in advance of operation to avoid potential mistakes.

Due to nutritional variation, some strains may result in poor growth or fail to grow on this medium. Additional agar media, such as Columbia agar with 5% sheep blood, should be used alongside the Enterococco Agar for primary isolation.

Organisms other than enterococci and *S. bovis* may be aesculin positive and may grow on this medium (e.g., *Pediococcus* and *Lactococcus* species).

Enterococco Agar is intended as an aid in the diagnosis of infectious diseases, requiring further tests to complete the diagnostic results. All identification tests should ideally be performed from non-selective agar.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.

- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Enterococco Agar	Plate 90 mm	20 plates	11057

Significant changes from previous version:

Document	Release Date	Change Summary
11057_IFU-0	2022-10-12	Document creation (layout and content in compliance with IVDR 2017/746)

This IFU document and the SDS are available from the online Support Center:

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liofilchem.com/ifu-sds

References / Riferimenti

1. Public Health England (2021) UK SMI ID 4 issue 4: Identification of *Streptococcus* species, *Enterococcus* species and morphologically similar organisms. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/1019677/ID_4i4.pdf
2. EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
3. Public Health England (2017) Inoculation of culture media for bacteriology. UK Standards for Microbiology Investigations. Q 5 Issue 2. <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
4. Clinical and Laboratory Standards Institute (2004) Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard, 3rd ed. M22-A3. CLSI, Wayne, PA.
5. ISO 7899-2:2000. Water Quality – Detection and enumeration of intestinal enterococci. Part 2: Membrane filtration method.
6. Facklam, R.R. (1972) Recognition of group D strptococcal species of human origin by biochemical and physiological test. Appl. Microbiol., 23; 1131- 1139.
7. Isenberg, H.D., Goldber D. (1970) Laboratory studies with a selective Enterococcus medium – Appl. Microbiol. vol. 20, p. 433-436.

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