Sabouraud CAF Agar + Actidione

Selective medium for isolation and cultivation of dermatophytes and other fungi.

DESCRIPTION

Sabouraud CAF Agar + Actidione is a selective medium used for the isolation and cultivation of dermatophytes and other pathogenic fungi from clinical specimens with mixed flora.

TYPICAL FORMULA (g/l)

Enzymatic Digest of Casein	5.0
Enzymatic Digest of Animal Tissue	5.0
Glucose	40.0
Chloramphenicol	0.05
Actidione (Cycloheximide)	0.5
Agar	15.0
Final pH 6.8 ± 0.2 at 25°C	

METHOD PRINCIPLE

Enzymatic digests of casein and enzymatic digest of animal tissue provide amino acids, carbon, nitrogen, vitamins and minerals for the growth of fungi. The high glucose concentration makes this medium particularly well suited for cultivating fungi. Chloramphenicol is a broad-spectrum antibiotic inhibitory to a wide range of Gram-negative and Gram-positive bacteria. Cycloheximide suppresses the growth of saprophytic fungi. Agar is the solidifying agent.

TEST PROCEDURE

Inoculate either plates or slant tubes by streaking directly the sample onto the agar surface. Streak the specimen as soon as possible after it is received in the laboratory. Incubate aerobically at 25-30°C for 2-7 days.

INTERPRETING RESULTS

Examine containers for fungal colonies exhibiting typical color and morphology. Biochemical tests and serological procedures should be performed to confirm findings.

Transfer of growth from slants to plated media may be required in order to obtain pure cultures of fungi.

APPEARANCE

Slightly opalescent, amber.

STORAGE

Store tubes and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Medium in tubes: 1 year.

Ready-to-use plates: 6 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Inoculum for selectivity: 10⁴-10⁶ CFU.

Incubation conditions: aerobically at 25-30°C for 2-7 days.

QC Table.

Microorganism		Growth
<i>Candida albicans</i>	ATCC® 10231	Good
<i>Trichophyton mentagrophytes</i>	ATCC® 9533	Good
<i>Aspergillus niger</i>	ATCC® 16404	Partially to completely inhibited
<i>Saccharomyces cerevisiae</i>	ATCC® 9763	Inhibited
<i>Escherichia coli</i>	ATCC® 8739	Inhibited

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. Sabouraud. 1892. Ann. Dermatol. Syphil. 3:1061.
2. Larone. 1995. Medically important fungi: a guide to identification, 3rd ed. American Society for Microbiology, Washington, D.C.
3. Wehr and Frank (ed.). 2004. Standard methods for the examination of dairy products, 17th ed. American Public Health Association, Washington, D.C.

PRESENTATION

		Contents	Ref.
Sabouraud CAF Agar + Actidione	90 mm ready-to-use plates	20 plates	11236
Sabouraud CAF Agar + Actidione	Slant tubes	10 x 9 ml tubes	30024
Sabouraud CAF Agar + Actidione	Slant tubes	20 x 9 ml tubes	31024

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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