Chromatic Detection

Non-selective medium for the isolation and differentiation of microorganisms.

INTENDED PURPOSE

Chromogenic non-selective medium for the enumeration and the presumptive identification of microorganisms directly from clinical and nonclinical specimens. This medium is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Chromatic Detection is a chromogenic medium used for the isolation and differentiation of microorganisms directly from clinical and nonclinical specimens.

Concomitant cultures are necessary to recover organisms for further microbiological testing or epidemiological typing.

The medium allows to carry out indole test for confirmation of *Escherichia coli*.

TYPICAL FORMULA*

(g/litre)

Peptone	14.0
Tryptone	6.0
Yeast Extract	3.0
Sodium Chloride	5.0
Chromogenic Mix	13.1
Agar	15.0

Final pH 7.2 ± 0.2 at 25°C

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptone and tryptone provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium. Chromogenic mix allows to identify microorganisms based on the color and morphology of the colonies. Agar is the solidifying agent.

PREPARATION

Dehydrated medium

Suspend 56.1 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil for 1 minute shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.

Medium in bottles

Melt the content of the bottle in a water bath at 100°C (losing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: Autoclave, water bath, sterile Petri plates, test tubes, inoculating loops, swabs, incubator, quality control organisms.

SPECIMENS

Chromatic Detection can be used for enumeration and the presumptive identification of microorganisms from all types of specimens. Good laboratory practices for collection, transport and storage of the clinical specimens should be applied.

Refer to specific guidelines for more information about specimen collection and preparation.

TEST PROCEDURE

Ensure there is no visible moisture on the plates before use.

Inoculate the plates by directly streaking the specimen on the agar surface or spread the sample from an enrichment culture to obtain well-isolated colonies.

Incubate at $35 \pm 2^{\circ}\text{C}$ for 18-24 h.

For more details, consult appropriate guidance.

INTERPRETING RESULTS

After incubation, observe the color and the morphology of the colonies for presumptive identification according to the following table.

Microorganism	Typical Colony
<i>Escherichia coli</i>	Pink-reddish-mauve
<i>Klebsiella</i> spp., <i>Enterobacter</i> spp., <i>Serratia</i> spp.	Green-blue
<i>Proteus</i> spp.	Brown
<i>Pseudomonas</i> spp.	Yellowish-green
<i>Staphylococcus aureus</i>	Cream
<i>Enterococcus faecalis</i>	Green-turquoise
<i>Staphylococcus saprophyticus</i>	Light pink

See Appendix I for example pictures.

Note: Further testing should be conducted to confirm the presumptive identification of organisms isolated on this medium.

STORAGE

The powder is very hygroscopic, store the powder at $10-30^{\circ}\text{C}$, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at $2-8^{\circ}\text{C}$ away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 2 years.

Medium in bottles: 1 years.

Ready-to-use plates: 4 months.

QUALITY CONTROL

Appearance of Dehydrated Medium: Free-flowing, homogeneous, beige.

Appearance of Prepared Medium: Slightly opalescent, amber.

Expected Cultural Response:

Control strain	Inoculum	Incubation	Criteria	Specification
<i>Escherichia coli</i>	50-100 CFU	18-24 h/ $35 \pm 2^{\circ}\text{C}$	Good growth ($P_R \geq 0.7$)	Pink-reddish-mauve colonies
<i>Klebsiella pneumoniae</i>				Green-blue colonies
<i>Proteus mirabilis</i>				Brown colonies
<i>Pseudomonas aeruginosa</i>				Yellowish-green colonies
<i>Staphylococcus aureus</i>				Cream colonies
<i>Enterococcus faecalis</i>				Green-turquoise colonies

A productivity ratio (P_R) of 0.7 is equivalent to a recovery rate of 70%.

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of Chromatic Detection was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Invalid results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology. The operator should understand the principles of the procedures, including its performance limitations, in advance of operation to avoid potential mistakes.

Chromatic Detection is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

The final identification must be confirmed by biochemical tests, immunological tests. They can be done directly from the suspicious colonies observed on the medium.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the various configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Chromatic Detection	Plate 90 mm	20 plates	11611
	Bottle	6 x 100 ml	481130
	Dehydrated media	500 g	610612
		100 g	620612
		5 kg	6106125

Revision History

Revision	Release Date	Change Summary
1	2023-11-09	Updated layout and content in compliance with IVDR 2017/746

This IFU document and the SDS are available from the online Support Center:

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liofilchem.com/ifu-sds

References / Riferimenti

1. EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. CLSI 2004 Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard -Third Edition. CLSI document M22-A3. Wayne, PA: Clinical and Laboratory Standards Institute.
3. Samra Z., et al. 1998 Evaluation of use of a new chromogenic Agar in detection of urinary tract pathogens. J. Clin. Microbiol. 36: 990-994.
4. Merlino J., S. Siarakas, G.J. Robertson, G.R. Funnel, T. Gottlieb, and R. Bradbury 1996 Evaluation of Colorex Orientation for differentiation and presumptive identification of Gram-negative bacilli and Enterococcus species. J. Clin. Microbiol. 34: 1788-1793.

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
liofilchem@liofilchem.com

Chromatic™ Detection

Chromogenic medium for the enumeration and identification of microorganisms directly from clinical and nonclinical specimens.

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