Chromatic Salmonella

Chromogenic medium for the isolation and differentiation of *Salmonella* spp.

INTENDED PURPOSE

Chromogenic selective medium for the detection and presumptive identification of *Salmonella* spp from clinical and nonclinical specimens. This medium is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Salmonella is ubiquitous in animal populations and is generally isolated from the intestinal tract of animals and humans. It is one of the most prevalent organisms associated with food-borne illnesses, which is often linked to animal origin. Illnesses caused by *Salmonella* have been associated with poultry, beef, chocolate, dairy and vegetable products.

Chromatic Salmonella is a chromogenic medium used with supplements for the selective isolation and differentiation of *Salmonella* species from different materials.

Clinical specimens such as stools and rectal swabs can be inoculated directly onto the agar plate.

This medium can also be used as second isolation agar in the four-step procedure described in EN ISO 6579-1 for detection of *Salmonella* in food, animal feed and environmental samples from the food production area, including samples from the primary production stage such as animal faeces, dust and swabs.

TYPICAL FORMULA* Per Litre of Purified Water

Base medium

Proteose Peptone	7.0 g
Meat Extract	1.0 g
Yeast Extract	3.0 g
Sodium Chloride	5.0 g
Chromogenic Mix	3.7 g
Agar	15.0 g
Final pH 7.5 ± 0.2 at 25°C	

*Adjusted and/or supplemented as required to meet performance specifications.

Prepared medium

Base

Chromatic Salmonella

+
+

Supplement

Tween 20 Supplement

Chromatic Salmonella Selective Supplement (optional)
--

METHOD PRINCIPLE

Proteose peptone and meat extract provide amino acids, nitrogen, carbon, minerals and vitamins essential for bacterial growth. Yeast extract is a source of vitamins, particularly of group B. Sodium chloride maintains the osmotic balance of the medium. The chromogenic mix allows the identification of bacteria on the basis of the colony colour. Agar is the solidifying agent.

The following supplements are added to the base medium:

- Tween 20 Supplement which enhances microbial growth.
- Chromatic Salmonella Selective Supplement containing selective agents to inhibit growth of competing microbiota such as species of *Proteus*, *Pseudomonas* and *Aeromonas*.

PREPARATION

Selective supplement

Aseptically reconstitute the contents of one vial of Chromatic Salmonella Selective Supplement with 5 ml sterile distilled or deionized water and swirl to complete dissolution (1 vial per 500 ml of medium).

Dehydrated medium

Suspend 34.7 g of the powder in 1 liter of distilled or deionized water. Add 3 ml of Tween 20 Supplement. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 100°C for 5 minutes.

If desired, aseptically add rehydrated contents of 2 vials (10 ml) of Chromatic Salmonella Selective Supplement to the melted agar medium, cooled to 45-50°C. Mix well and pour into sterile Petri dishes.

Medium in bottles

Melt the content of the bottle in a water bath at 100°C (losing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: Autoclave, water bath, sterile Petri plates, test tubes, inoculating loops, swabs, incubator, quality control organisms.

SPECIMENS

Chromatic Salmonella can be used for detection and the presumptive identification of *Salmonella* spp. from all types of specimens.

Faecal specimens should be sampled at the acute stage, before antimicrobial therapy (where possible) and examined as soon as possible after collection. Ensure that specimen is placed in appropriate package (clean, leak-proof container) and, in case the laboratory processing is delayed (>1 h), use a suitable transport medium, such as Cary-Blair for transportation/storage.

Good laboratory practices for collection, transport and storage of the clinical specimens should be applied. Refer to specific guidelines for more information about specimen collection and preparation.

TEST PROCEDURE

Ensure there is no visible moisture on the plates before use.

Inoculate the plates by directly streaking the specimen on the agar surface or spread the sample from an enrichment culture (e.g. Rappaport Vassiliadis Soy Broth or Selenite Broth) to obtain well-isolated colonies.

For testing of food samples, Chromatic Salmonella can be used according to EN ISO 6579-1 alongside XLD agar after pre-enrichment in Buffered Peptone Water (BPW), followed by enrichment in RVS broth or MSRV agar and Muller Kauffmann Tetrathionate Novobiocin (MKTTn) broth.

For detailed instructions refer to the technical sheet of the medium being used.

Incubate inoculated plates of Chromatic Salmonella aerobically at $35 \pm 2^\circ\text{C}$ for 24 ± 3 h. An additional 24 h incubation may be required.

For more details, consult appropriate guidance.

INTERPRETING RESULTS

Salmonella spp. (including *S. Typhi*, *S. Paratyphi*, lactose-positive, sucrose-positive) form typically light mauve to mauve colonies.

E. coli, *Enterobacter* spp and *Klebsiella* spp, grow as blue-green colonies.

Other organisms, if not inhibited, will appear as colourless.

Notes: Further testing should be conducted to confirm the presumptive identification of organisms isolated on this medium. For example, oxidase test allows to distinguish between *Salmonella* (oxidase negative) and rare strains of *Pseudomonas* (oxidase positive) which are able to growth on this medium with mauve colonies. Subculture to a non-selective Nutrient Agar before performing appropriate biochemical and serological tests.

STORAGE

Store the freeze-dried supplement at 2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature.

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed.

Store bottles and prepared plates at 2-8°C away from direct light in their original pack until just prior to use. Avoid quick temperature shifts of plated medium to prevent condensation.

Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Freeze-dried Supplement: 2 years

Dehydrated medium: 2 years.

Medium in bottles: 1 year

Ready-to-use plates: 6 months

QUALITY CONTROL

Appearance of Selective Supplement: Whitish lyophilized pellet, limpid solution once reconstituted.

Appearance of Dehydrated Medium: Free-flowing, homogeneous, beige.

Appearance of Prepared Medium: Clear, beige.

Expected Cultural Response:

Control strain		Inoculum	Incubation	Criteria	Colony color
<i>Salmonella</i> Typhimurium	WDCM 00031 (ATCC® 14028; NCTC 12023)	50-100 CFU	24 ± 3 h / 35 ± 2°C	Good growth ($P_R \geq 0.5$)	Mauve
<i>Salmonella</i> Enteritidis	WDCM 00030 (ATCC® 13076; NCTC 12694)				
<i>Escherichia coli</i>	WDCM 00013 (ATCC® 25922; NCTC 12241)	10 ³ -10 ⁴ CFU		Growth	Blue-green
<i>Enterobacter cloacae</i>	WDCM 00082 (ATCC® 23355; NCTC 13380)			Growth	Blue-green
<i>Shigella flexneri</i>	WDCM 00126 (ATCC® 12022; NCTC 12698)			Growth	Colorless
<i>Proteus mirabilis</i> *	WDCM 00023 (ATCC® 29906; NCTC 11938)	10 ⁴ -10 ⁶ CFU		Partial to complete inhibition	---
<i>Pseudomonas aeruginosa</i> *	WDCM 00025 (ATCC® 27853; NCTC 12903)				---

(*) The strain will be able to grow (colourless colonies) if the medium is used without the selective supplement.

A productivity ratio (P_R) of 0.5 is equivalent to a recovery rate of 50%.

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of Chromatic Salmonella agar was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Invalid results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology. The operator should understand the principles of the procedures, including its performance limitations, in advance of operation to avoid potential mistakes.

Due to nutritional variation, some *Salmonella* strains may produce a weak coloration, result in poor growth or fail to grow on this medium. Additional agar media, such as MacConkey agar, should be used alongside Chromatic Salmonella for primary isolation.

Certain strains of Gram-negative bacteria other than *Salmonella* may produce colonies resembling *Salmonella*. Chromatic Salmonella is intended as an aid in the diagnosis of infectious diseases, requiring further tests to complete the diagnostic results. All identification tests should ideally be performed from non-selective agar.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

See ordering info below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Chromatic Salmonella	Plate 90 mm	20 plates	11614
	Bottle	6 x 100 ml	481140
	Dehydrated media	100 g	620611
500 g		610611	
Chromatic Salmonella Selective Supplement	Lyophilized supplement	10 vials	81102
Tween 20 Supplement	Liquid supplement	2 x 50 ml	80032
Tween 20 Supplement	Liquid supplement	4 x 50 ml	80432

Revision History

Revision	Release Date	Change Summary
0	2024-01-16	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0


In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.


*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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

References / Bibliografia

1. US Food and Drug Administration. Bacteriological Analytical Manual (BAM) Chapter 5: Salmonella. Rev 2022-03-18.
2. Public Health England (2021) UK SMI ID 24 issue 4: identification of *Salmonella* species. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/969498/ID_24i4.pdf
3. EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
4. EN ISO 6579-1:2017+Amd1:2019. Microbiology of the food chain – Horizontal method for the detection, enumeration and serotyping of Salmonella – Part 1: Horizontal method for the detection of Salmonella spp.
5. Public Health England (2017) Inoculation of culture media for bacteriology. UK Standards for Microbiology Investigations. Q 5 Issue 2. <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
6. D'Aoust, Mauer and Bailey (2001) In Doyle, Beuchat, and Montville (ed.) Food microbiology: fundamentals and frontiers, 2nd ed. American Society for Microbiology, Washington, DC.
7. Bopp, Brenner, Wells and Strockbine (1999) In Murray, Baron, Pfaller, Tenover and Tenover (ed.). Manual of clinical microbiology, 7th ed American Society for Microbiology, Washington, DC.

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