Chromatic™ MH

Chromogenic Mueller Hinton for the microbial differentiation simultaneously with susceptibility testing directly from clinical specimens.

DESCRIPTION

Chromatic™ MH is a chromogenic medium (Mueller Hinton base) used for the microbial differentiation simultaneously with susceptibility testing of various microorganisms* directly from clinical specimens. AST can be performed on this medium by using either antibiotic discs (Kirby-Bauer method) or Liofilchem® MIC Test Strip(s) (Gradient test method).

* AST optimized for Gram-negative bacteria.

TYPICAL FORMULA

	(g/l)
Beef Extract	2.0
Acid Hydrolysate of Casein	17.5
Starch	1.5
Chromogenic Mix	0.5
Agar	17.0
Final pH 7.3 ± 0.2 at 25°C	

METHOD PRINCIPLE

Acid hydrolysate of casein and beef extract supply amino acids, nitrogen, carbon, minerals, vitamins and other nutrients which support the growth of microorganisms. Starch acts as a protective colloid against toxic molecules which can be present in the medium. Hydrolysis of starch during autoclaving supplies a little amount of glucose, which is a source of energy. Chromogenic mix allows the identification of microorganisms on the basis of the colony color and morphology. Agar is the solidifying agent.

The Kirby-Bauer method is based on the diffusion, through the agar, of the antimicrobial substance which soaks the paper disc. Each disc has a single concentration of the antimicrobial agent that inhibits the microorganism growth showing an halo around the disc. The diameter of the inhibiting halo is correlated with the Minimal Inhibitory Concentration (MIC).

The Gradient Strip method is a quantitative assay for determining the Minimum Inhibitory Concentration (MIC) of an antimicrobial agent against a microorganism and for detecting the resistance mechanisms. Each MIC Test Strip contains a preformed exponential gradient of the antimicrobial agent to apply onto an inoculated agar surface. The MIC value is read where the edge of the inhibition ellipse intersects the strip.

Chromatic™ MH Agar is manufactured to contain low levels of thymine and thymidine and controlled levels of calcium and magnesium ions.

PREPARATION

Dehydrated medium Suspend 38.5 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes. Cool to 45-50°C. Aseptically dispense in Petri dishes.

TEST PROCEDURE

1. Inoculate the plate by streaking directly 10µl specimen onto the agar surface. If required dilute the specimen with saline solution and mix with vortex before spreading.
2. Apply antimicrobial discs or MIC Test Strip(s) onto the surface of the inoculated agar plate.
3. Incubate aerobically at 37°C for 18-24 hours.

INTERPRETING RESULTS

Preliminary identification

After incubation observe the color of the colonies and interpret the results as indicated in the ID table.

ID Table.

Microorganism	Typical colony color
<i>E. coli</i>	Pink-reddish-mauve
<i>Klebsiella</i> spp, <i>Enterobacter</i> spp, <i>Serratia</i> spp	Green-blue
<i>Proteus</i> spp	Brown
<i>Pseudomonas</i> spp	Yellowish-green
<i>S. aureus</i>	Cream
<i>E. faecalis</i>	Green-turquoise

Antimicrobial Susceptibility Testing

For each microorganism measure the diameter of the zone of complete inhibition around the disc or read the MIC value. Interpret zone diameters or MIC values by referring to the current standards and report the organism as susceptible, intermediate or resistant to the agents that have been tested.

See pictures in Appendix I.

APPEARANCE

Dehydrated medium: fine, dry, homogeneous, free of extraneous material, beige.
Prepared medium: slightly opalescent, amber.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed.
Store prepared plates at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.
Ready-to-use plates: 4 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table. Antibiotic discs are then promptly placed onto the agar surface.

Inoculum: 0.5 Mc Farland.

Incubation conditions: aerobically at 35±2°C for 18-24 h.

QC Table.

Microorganism		Specification
<i>Escherichia coli</i>	ATCC® 25922	Pink-reddish-mauve colonies, inhibition halo diameters as recommended by CLSI M100 for Mueller Hinton agar
<i>Pseudomonas aeruginosa</i>	ATCC® 27853	Yellowish-green colonies, inhibition halo diameters as recommended by CLSI M100 for Mueller Hinton agar

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use.
The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. Bauer, Kirby, Sherris and Turck. Am. J. Clin. Pathol. 1966, 45:493.
2. Bouza E., Burillo A. and Muñoz P. Role of the Microbiology Laboratory in the Diagnosis of Ventilator-Associated Pneumonia, in Nosocomial Pneumonia: Strategies for Management (ed J. Rello), John Wiley & Sons, Ltd, Chichester, UK, 2007.
3. Iztok Strumbelj. New broth - chromogenic Mueller Hinton agar procedure for urine samples - next-day result of Enterobacteriaceae antimicrobial susceptibility testing. ECCMID 2013, poster no. 1547
4. F.H. M'Zali, L. Stecken, X. Berard, A. Quinart and F. Sztark. A novel, direct susceptibility testing method of sonicated vascular prosthetic graft samples by combination of Liofilchem® MIC Test Strip and chromogenic agar plates. ECCMID 2013, poster no. 1562.

PRESENTATION		Contents	Ref.
Chromatic™ MH	Dehydrated medium	500 g of powder	611618
Chromatic™ MH	Dehydrated medium	100 g of powder	621618
Chromatic™ MH	90 mm ready-to-use plates	20 plates	11618
Chromatic™ MH	140 mm ready-to-use plates	10 plates	10246

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care	 Keep away from sunlight
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse	

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