Chromatic CRE

Chromogenic medium for the detection of carbapenem-resistant Enterobacteriaceae, directly from clinical specimens.

INTENDED USE

Chromatic CRE is a selective and differential chromogenic medium used for the qualitative and presumptive detection of carbapenem-resistant Enterobacteriaceae (CRE) directly from clinical specimens.

This medium is intended as an aid in the detection of the following bacteria:

Escherichia coli and KES (*Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Enterobacter cloacae* complex, and *Serratia marcescens*) that are non-susceptible to carbapenems.

Chromatic CRE is not intended to diagnose infection or guide therapy. Results can be interpreted after incubation for 24 hours. Subculture to a non-selective medium is required for confirming identification, antimicrobial susceptibility testing and epidemiological typing.

SUMMARY AND EXPLANATION

CRE are multidrug-resistant organisms that can cause serious infections and require interventions in healthcare settings to prevent spread.

Carbapenemase enzyme production is one of the mechanisms, probably the most alarming, that can lead to carbapenem resistance in Enterobacteriaceae. Carbapenemases represent an emerging threat to healthcare as they significantly reduce treatment options. The prevalence of carbapenemase producing microorganisms remains low, but this phenomenon is increasing and spreading worldwide.

With no enrichment steps required, Chromatic CRE allows for early detection of colonization by carbapenemase-producing Enterobacteriaceae (CPE), which is essential for control of CPE infections and to prevent the spread of resistance.

TYPICAL FORMULA*

	(g/litre)
Peptones	30.0
Selective Mix	11.0
Chromogenic Mix	1.0
Agar	15.0

Final pH 7.0 ± 0.2 at 25°C

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptones provide amino acids, nitrogen, minerals, vitamins, and nutrient factors for growth of bacteria. The chromogenic and selective mix facilitate the identification of bacteria on the basis of the color and colony morphology, inhibiting most of the non-Enterobacteriaceae. Agar is the solidifying agent.

TEST PROCEDURE

Inoculate the specimen with a swab or loop onto the medium and streak for isolation, using a loop. Incubate the agar plates in an inverted position (medium side up) aerobically at 35 ± 2°C for 24-48 hours.

Notes:

Do not inoculate more than one specimen per plate. Chromatic CRE is primarily used with rectal and perianal swabs but can be used with specimens from other body sites suspected to contain CPE.

Avoid exposure to light during incubation as this might destroy the chromogens.

INTERPRETING RESULTS

At the end of incubation, observe the appearance of the colonies and interpret the results as indicated in the ID Table.

ID Table.

Microorganism	Typical colony color
Carbapenem non-susceptible <i>Escherichia coli</i>	Pink to red
Carbapenem non-susceptible KES	Blue-violet to blue-green
Other organisms (if not inhibited)	Colorless or natural color

See Appendix I for example pictures.

Notes:

A lack of growth or the absence of pink to red or blue-violet to blue-green colonies on Chromatic CRE does not preclude the presence of *Escherichia coli* and KES that are non-susceptible to carbapenems.

Appropriate susceptibility tests, molecular methods or phenotypic methods must be performed to confirm the presence of CPE isolates.

APPEARANCE

Light amber, slightly opalescent.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

3 months.

QUALITY CONTROL

The medium is inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Inoculum for selectivity: 10⁴-10⁶ CFU.

Incubation conditions: aerobically at 35 ± 2°C for 24-48 hours.

QC Table.

Microorganism		Growth	Colony color
<i>Klebsiella pneumoniae</i>	ATCC® BAA-1705	Good	Blue-violet
<i>Klebsiella pneumoniae</i>	ATCC® 700603	Inhibited	—
<i>Escherichia coli</i>	ATCC® 25922	Inhibited	—

LIMITATIONS

Chromatic CRE is considered for **screening**, and **not for final identification** of carbapenemase producers.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

Product Description	Format	Packaging	Ref.
Chromatic CRE	Plate 90 mm	20 plates	11619

This document is available from the online Support Center:

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liofilchem.com/ifu-sds

BIBLIOGRAPHY / BIBLIOGRAFIA

- Centers for Disease Control and Prevention. <https://www.cdc.gov/hai/organisms/cre/technical-info.html>
- Kochar S, Sheard T, Sharma R, Hui A, Tolentino E, Allen G, Landman D, Bratu S, Augenbraun M, Quale J. (2009) Success of an infection control program to reduce the spread of carbapenem-resistant *Klebsiella pneumoniae*. Infect. Control Hosp. Epidemiol. 30(5):447-52.
- Podschun R, Ullman U (1998). *Klebsiella spp* as Nosocomial Pathogens: Epidemiology, Taxonomy, Typing Methods, and Pathogenicity Factors. Clinical Microbiology Reviews 11 (4): 589–603.

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
Chromatic CRE





Chromogenic medium for detection of carbapenem-resistant
Enterobacteriaceae directly from clinical specimens.

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