Chromatic™ ESB+AmC

ENGLISH

Chromogenic media for detection of ESBs and AmC in Enterobacteriaceae directly from clinical specimens.

DESCRIPTION

Chromatic™ ESB+AmC is a chromogenic medium used for detection of ESB and AmC producers. ESBs are enzymes that hydrolyze most penicillins and cephalosporins, inhibited by β -lactamase inhibitors such as clavulanic acid, sulbactam and tazobactam. AmC-type enzymes hydrolyze penicillins, cephalosporins and monobactam and are poorly inhibited by the classical ESB inhibitors, especially clavulanic acid.

ESB-producing Enterobacteriaceae have become one of the most important causes of nosocomial community-acquired infections caused by *Escherichia coli* and *Klebsiella pneumoniae*, but also other Gram-negative bacteria. Patients infected by ESB producers often received inadequate empirical therapy until the pathogen's resistance is recognized. Early detection is, therefore, essential in order to limit the spread of these pathogens. Organisms producing enough AmC β -lactamase will typically give a positive ESB screening test but fail the confirmatory test involving increased sensitivity with clavulanic acid.

TYPICAL FORMULA

	(g/l)
Peptone Mix	43.2
Chromogenic Mix	1.0
Selective Mix	0.5
Agar	15.0

Final pH 7.2 \pm 0.2 at 25°C

METHOD PRINCIPLE

Peptones supply amino acids, nitrogen, carbon, minerals, vitamins and other nutrients which support the growth of microorganisms. Chromogenic mix allows the identification of microorganisms on the basis of the colony color and morphology. Selective mix inhibits the ESB-non-producing organisms, but allows the growth of AmC producers. Agar is the solidifying agent.

TEST PROCEDURE

Inoculate the plates by streaking directly the specimen onto the agar surface. Incubate aerobically at 37°C for 18-24 hours.

INTERPRETING RESULTS

After incubation observe the color of the colonies and interpret the results as indicated in the ID table.

ID Table.

Microorganism*	Typical colony color
<i>E. coli</i>	Pink-reddish-mauve
<i>Klebsiella</i> spp, <i>Enterobacter</i> spp, <i>Serratia</i> spp	Green-blue
<i>Proteus</i> spp	Brown

*Notice that only ESB-producing organisms grow on this medium, irrespective whether they are or not AmC producers too.

See pictures in Appendix I.

APPEARANCE

Slightly opalescent, amber.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

4 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 10-100 CFU/ml.

Inoculum for selectivity: 10⁴-10⁵ CFU/ml.

Incubation conditions: aerobically at 35±2°C for 18-24 h.

QC Table.

Microorganism		ESBL phenotype	AmpC phenotype	Growth	Specification
<i>Enterobacter cloacae</i>	ATCC® BAA-1143	Negative	Positive	Good	Green-blue colonies
<i>Klebsiella pneumoniae</i>	ATCC® 700603	Positive	Negative	Good	Green-blue colonies
<i>Klebsiella pneumoniae</i>	ATCC® BAA-1144	Negative	Positive	Good	Green-blue colonies
<i>Escherichia coli</i>	ATCC® 25922	Negative	Negative	Inhibited	---

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 1.0, 2013.
2. Podschun R, Ullman U. *Klebsiella* spp as Nosocomial Pathogens: Epidemiology, Taxonomy, Typing Methods, and Pathogenicity Factors. *Clinical Microbiology Reviews*. 1998; 11 (4): 589–603.
3. Geiss H.K. Comparison of two test kits for rapid identification of *Escherichia coli* by a beta-glucuronidase assay. *European Journal of Clinical Microbiology & Infections Diseases*. 1990; 9 (2):151-152.

PRESENTATION

		Contents	Ref.
Chromatic™ ESBL+AmpC	90 mm ready-to-use plates	20 plates	11629

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care	 Keep away from sunlight
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse	

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BIBLIOGRAFIA

1. EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 1.0, 2013.
2. Podschun R, Ullman U. *Klebsiella* spp as Nosocomial Pathogens: Epidemiology, Taxonomy, Typing Methods, and Pathogenicity Factors. *Clinical Microbiology Reviews*. 1998; 11 (4): 589–603.
3. Geiss H.K. Comparison of two test kits for rapid identification of *Escherichia coli* by a beta-glucuronidase assay. *European Journal of Clinical Microbiology & Infections*. 1990; 9 (2):151-152.

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
liofilchem@liofilchem.net

Chromatic™ ESB�+AmpC

Chromogenic medium for detection of ESB� and AmpC in Enterobacteriaceae directly from clinical specimens.

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