Chromatic Colistin

Chromogenic medium for detection of colistin-resistant organisms from clinical specimens and other materials.

DESCRIPTION

Chromatic Colistin is a selective chromogenic medium used for the isolation and differentiation of Colistin-resistant Gram-negative bacteria.

The medium is intended for screening of faecal samples and can also be used to confirm colistin resistance of isolates from non-selective and non-chromogenic culture media.

TYPICAL FORMULA	(g/l)
Meat Peptone	14.0
Casein Peptone	6.0
Yeast Extract	3.0
Sodium Chloride	5.0
Chromogenic and Selective Mix	14.0
Agar	15.0

Final pH 7.3 ± 0.2 at 25°C

METHOD PRINCIPLE

Peptones provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of group B. Sodium chloride maintains the osmotic balance of the medium. Chromogenic and selective mix allows to identify microorganisms on the basis of the color and morphology of the colonies while inhibiting most of Gram-positive bacteria, yeasts and moulds as well as Gram-negative bacteria that are sensitive to colistin. Agar is the solidifying agent.

TEST PROCEDURE

For screening use, inoculate the medium with rectal swabs or fresh stools (1 g solid sample or 50 µl liquid) by spread plating technique (100 µl inoculum) after selective enrichment in appropriate media. For example, place a Colistin sulfate 10 µg disc (ref. 9023) into a 9 ml tube of Brain Heart Infusion broth (ref. 24480), inoculate with the sample and incubate at 35 ± 2°C for 4-5 hours.

For confirmation of resistance detection, inoculate by streaking a pure strain obtained by isolation on a non-selective and non-chromogenic plate.

Incubate aerobically at 35 ± 2°C for 18-24 hours. An additional 24 h incubation may be required.

INTERPRETING RESULTS

After incubation observe the color and the morphology of the colonies and interpret the results as indicated in the ID table.

ID Table.

Microorganism	Typical colony appearance
<i>E. coli</i>	Pink-reddish-mauve
<i>Klebsiella</i> , <i>Enterobacter</i> , <i>Serratia</i> spp.	Green-blue
<i>Pseudomonas</i> spp.	Yellowish-green
Other colistin-resistant organisms	Colourless

See pictures in Appendix I.

NOTES:

- Only colistin-resistant Gram-negative microorganisms may develop on this medium, including bacteria naturally resistant to colistin such as *Serratia marcescens* and *Proteus mirabilis*. Gram-positive bacteria, yeasts, moulds and other Colistin-sensitive microorganisms are inhibited.
- Isolates obtained on this medium should be confirmed for colistin resistance by additional tests, such as through broth microdilution method with ComASP Colistin (ref. 75001).

APPEARANCE

Slightly opalescent, light amber.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

4 months.

QUALITY CONTROL

The plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Inoculum for selectivity: 10⁴-10⁶ CFU.

Incubation conditions: aerobically at 35 ± 2°C for 18-24 hours.

QC Table.

Microorganism		Growth	Colony color
<i>Escherichia coli</i> (MCR-1 pos)	NCTC 13846	Good	Pink-reddish-mauve
<i>Escherichia coli</i>	ATCC® 25922	Inhibited	---
<i>Enterococcus faecalis</i>	ATCC® 29212	Inhibited	---

PERFORMANCE

The sensitivity of the medium was evaluated using stool specimens spiked with colistin resistant target organisms, i.e. 16 *E. coli* and 2 *Salmonella* spp. The specificity was evaluated using 50 fresh clinical samples. The medium demonstrated an overall sensitivity of 85% and a specificity of 100%. The limit of detection (LoD) was 10³-10⁴ CFU.

A total of 67 characterised pure strains (from clinical origin) were tested on both Chromatic Colistin and ComASP Colistin, which was considered as MIC testing reference method. The tested strains included 34 strains resistant to colistin, including 12 *mcr-1* positive strains, and 33 strains susceptible to colistin. These were the results: 97%(65/67) category agreement with 5.8%(2/34) very major discrepancy rate. The two very major discrepancies (false susceptible on Chromatic Colistin) were observed to *S. marcescens* with MIC of 16 and 8 µg/l. No major discrepancies (false resistant on Chromatic Colistin) were found.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

- Zackery P, Bulman et al. (2017) Polymyxin Combinations Combat *Escherichia coli* Harboring *mcr-1* and *bla*_{NDM-5}: Preparation for a Postantibiotic Era. *mBio*. 8(4): e00540-17.
- Liu YY., Wang Y., Walsh TR. et al. (2016) Emergence of plasmid-mediated colistin resistance mechanism MCR-1 in animals and human beings in China: a microbiological and molecular biological study. *Lancet Infect Dis*, 16(2): 161-8.
- Merlino J. et al. (1996) Evaluation of Colorex Orientation for differentiation and presumptive identification of Gram-negative bacilli and *Enterococcus* species. *J. Clin. Microbiol.* 34: 1788-1793.

PRESENTATION	Format	Packaging	Ref
Chromatic Colistin	90 mm Plate	20 plates	11640

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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BIBLIOGRAFIA

- Zackery P. Bulman et al. (2017) Polymyxin Combinations Combat *Escherichia coli* Harboring mcr-1 and bla_{NDM-5}: Preparation for a Postantibiotic Era. mBio. 8(4): e00540-17.
- Liu YY., Wang Y., Walsh TR. et al. (2016) Emergence of plasmid-mediated colistin resistance mechanism MCR-1 in animals and human beings in China: a microbiological and molecular biological study. Lancet Infect Dis, 16(2): 161-8.
- Merlino J. et al. (1996) Evaluation of Colorex Orientation for differentiation and presumptive identification of Gram-negative bacilli and *Enterococcus* species. J. Clin. Microbiol. 34: 1788-1793.

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