Chromatic Super CAZ/AVI

Chromogenic medium for the detection of
Ceftazidime-avibactam resistance in Gram-negative bacteria.

Instructions For Use

ENGLISH

INTENDED USE

Chromatic Super CAZ/AVI is a selective and differential chromogenic medium used for screening ceftazidime-avibactam resistance in *Enterobacteriales* and *Pseudomonas aeruginosa*.

Chromatic Super CAZ/AVI is not intended to diagnose infection or guide therapy. Results can be interpreted after incubation for 18-24 hours. Subculture to a non-selective medium is required for confirming identification, antimicrobial susceptibility testing and epidemiological typing.

SUMMARY AND EXPLANATION

Ceftazidime-avibactam (CZA) is a novel cephalosporin/beta-lactamase inhibitor combination with activity against multidrug-resistant gram-negative bacteria, mainly among *Enterobacteriales* and *Pseudomonas aeruginosa* strains. This includes isolates producing a variety of Ambler class A and C beta-lactamases, i.e. AmpC, ESBLs, and KPC, as well as some class D OXA-48-like enzymes. In contrast, CZA is ineffective against class B MBL producers. In addition, acquisition of resistance has been emerging shortly after introduction of CZA in clinical use, especially among carbapenem-resistant strains. This is a cause for concerns, explaining the need for the surveillance of CZA-resistant pathogens.

TYPICAL FORMULA*	(g/litre)
Peptone	14.0
Tryptone	6.0
Yeast Extract	3.0
Sodium Chloride	5.0
Chromogenic Mix	13.1
Ceftazidime	0.006
Avibactam	0.004
Daptomycin	0.01
Amphotericin B	0.005
Zinc Sulfate	0.07
Agar	15.0
Final pH 7.2 ± 0.2 at 25°C	

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptone and tryptone provide amino acids, nitrogen, minerals, vitamins, and nutrient factors for growth of bacteria. Yeast extract is a source of vitamins, particularly of group B. Sodium chloride maintains the osmotic balance of the medium. The chromogenic mix facilitates the identification of bacteria on the basis of the color and colony morphology, while the antimicrobial agents inhibit most of yeasts and bacteria with the exception of CZA-resistant isolates. Zinc sulfate is included to enhance the activity of MBL producers. Agar is the solidifying agent.

TEST PROCEDURE

Inoculate the specimen with a swab or loop directly onto the medium, or after enrichment in a broth, and streak for isolation, using a loop. Incubate the agar plates in an inverted position (medium side up) aerobically at 35 ± 2°C for 18-24 hours. If no colonies are seen, reincubate for a total of 48 hours.

Notes:

Do not inoculate more than one specimen per plate. Avoid exposure to light during incubation as this might destroy the chromogens.

INTERPRETING RESULTS

At the end of incubation, observe the appearance of the colonies and interpret the results as indicated in the ID Table.

ID Table.

Microorganism	Typical colony color
CZA-resistant <i>Escherichia coli</i>	Pink-reddish-mauve
CZA-resistant <i>Klebsiella, Enterobacter, Serratia</i> spp.	Green-blue
CZA-resistant <i>Proteus</i> spp.	Brown
CZA-resistant <i>Pseudomonas</i> spp.	Yellowish-green
Other organisms (if not inhibited)	Colorless or natural color

See Appendix I for example pictures.

Notes:

- Growth depends on the requirements of each individual microorganism. A lack of growth on Chromatic Super CAZ/AVI does not preclude the presence of bacteria with acquired CZA resistance.
- Appropriate susceptibility tests, molecular methods or phenotypic methods must be performed to confirm the presence of CZA-resistant isolates.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

Avoid quick temperature shifts to prevent condensation.

SHELF LIFE

4 months.

QUALITY CONTROL

Appearance of Medium: Slightly opalescent, amber.

Expected Cultural Response:

Microorganism		Inoculum	Incubation	Specification
<i>Klebsiella pneumoniae</i> (MBL/NDM-1 pos.)	ATCC® BAA-2146	50-100 CFU	18-24 h/ 35 ± 2°C	Good growth, green colonies ($P_R \geq 0.5$)
<i>Klebsiella pneumoniae</i>	ATCC® 700603	$10^4\text{-}10^6$ CFU	24-48 h/ 35 ± 2°C	Inhibition
<i>Escherichia coli</i>	ATCC® 25922			Inhibition
<i>Staphylococcus aureus</i>	ATCC® 29213			Inhibition
<i>Candida albicans</i>	ATCC® 10231			Inhibition

A productivity ratio (P_R) of 0.5 is equivalent to a recovery rate of 50%.

Please refer to the actual batch related Certificate of Analysis (CoA).

LIMITATIONS

Chromatic Super CAZ/AVI is considered for **screening**, and **not for final identification** of CZA-resistant isolates.

WARNING AND PRECAUTIONS

For professional use only. Operators must be trained and have certain experience in the laboratory methods. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

Product	Format	Packaging	Ref.
Chromatic Super CAZ/AVI	Plate 90 mm	20 plates	11641

This IFU document and the SDS are available from the online Support Center:

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