Sabouraud CAF Agar (Irradiated)

Instructions For Use
ENGLISH

Ready-to-use plated medium for detection of yeasts and moulds.

DESCRIPTION

Sabouraud CAF Agar (Irradiated) is a selective medium used for isolation and enumeration of fungi from sanitized surfaces as well as for monitoring fungal contamination by air sampling.

This medium is based on Sabouraud Dextrose Agar (SDA) with the addition of Chloramphenicol to inhibit contaminating bacteria while permitting successful selective isolation of yeasts and moulds. SDA is formulated according to the requirements of the harmonized USP/EP/JP method for determination of total combined yeast/mould count (TYMC) and detection of *Candida* spp in non-sterile pharmaceutical products. It also complies with EN ISO 11133 for microbiological examination of food, animal feed and water.

An irradiation dose of 9-20 kGy is used to sterilize the plated medium in its final packaging, which includes the cardboard box. Gamma irradiation of the product is indicated by an orange to red color of the irradiation indicator dot on the inner bag.

TYPICAL FORMULA*

	(g/litre)
Enzymatic Digest of Casein	5.0
Enzymatic Digest of Animal Tissue	5.0
Glucose	40.0
Chloramphenicol	0.5
Agar	15.0
Final pH 5.6 ± 0.2 at 25°C	

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Enzymatic digests of casein and enzymatic digest of animal tissue provide nitrogen and vitamins for the growth of fungi. The high glucose concentration along with the acid pH make this medium particularly well suited for cultivating fungi. Chloramphenicol is a broad-spectrum antibiotic inhibitory to a wide range of Gram-negative and Gram-positive bacteria. Agar is the solidifying agent.

TEST PROCEDURE

For surface hygiene assessment, take a swab sample and streak directly onto a 90 mm plate. Swabs can be used to sample irregular non-flat materials as well as a defined area of regular surfaces, for example by means of the sampling template 10x10 (ref. 96762).

Contact plates are used for active air sampling as well as for surface testing. Selected surfaces are sampled by firmly pressing the agar medium against the test area for about 10 seconds. Residues of the medium should be subsequently removed from the area tested.

For personnel monitoring, both contact plates (gown sampling) and 90 mm plates (finger prints) may be used.

After sampling, incubate the plates aerobically at 20-25°C for 2-7 days. Other temperatures and times of incubation may be chosen to suit specific requirements of target organisms but such conditions should be validated under the specific protocol used for environmental monitoring testing.

INTERPRETING RESULTS

Examine plates for fungal colonies exhibiting typical color and morphology.

Biochemical tests and serological procedures should be performed to confirm findings.

STORAGE

Store at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

Avoid quick temperature shifts to prevent condensation.

SHELF LIFE

90 mm plates: 6 months.

Contact plates: 9 months.

QUALITY CONTROL

Appearance of Medium: Slightly opalescent, amber.

Expected Cultural Response:

Control strain		Inoculum	Incubation	Specification
<i>Saccharomyces cerevisiae</i>	WDCM 00058 (ATCC® 9763, NCTC 10716)	≤ 100 CFU	44-74 h/ 22.5 ± 2.5°C	Good growth, 50-200% recovery
<i>Aspergillus brasiliensis</i>	WDCM 00053 (ATCC® 16404, NCPF 2275)			
<i>Candida albicans</i>	WDCM 00054 (ATCC® 10231, NCPF 3179)			
<i>Staphylococcus aureus</i>	WDCM 00032 (ATCC® 6538, NCTC 10788)	10 ⁴ -10 ⁶ CFU	70-74 h/ 32.5 ± 2.5°C	Inhibition
<i>Escherichia coli</i>	WDCM 00012 (ATCC® 8739, NCTC 12923)			

A productivity ratio (P_R) of 0.5 is equivalent to a recovery rate of 50%.

Please refer to the actual batch related Certificate of Analysis (CoA).

WARNING AND PRECAUTIONS

For professional use only. Operators must be trained and have certain experience in the laboratory methods. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the various configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Sabouraud CAF Agar (Irradiated)	Contact Plate	20 plates	15380S
Sabouraud CAF Agar (Irradiated)	90 mm Plate	20 plates	11035S

This IFU document and the SDS are available from the online Support Center:

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liofilchem.com/ifu-sds

BIBLIOGRAPHY / BIBLIOGRAFIA

1. European Pharmacopoeia 10th Ed. (2020) 2.6.12. Microbiological examination of non-sterile products: Microbial enumeration tests.
2. European Pharmacopoeia 10th Ed. (2020) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
3. United States Pharmacopoeia 42 NF 37 (2019) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
4. EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
5. ISO 18593:2018. Microbiology of the food chain - Horizontal methods for surface sampling.
6. Japanese Pharmacopoeia 17th Ed. (2017) 4.05 Microbiological examination of non-sterile products: Test for specified microorganisms.
7. Wehr and Frank (ed.) (2004) Standard methods for the examination of dairy products, 17th ed. American Public Health Association, Washington, D.C.
8. Larone (1995) Medically important fungi: a guide to identification, 3rd ed. American Society for Microbiology, Washington, D.C.
9. Sabouraud (1892) Ann. Dermatol. Syphil. 3:1061.

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