m-Faecal Coliform Agar

Selective medium for detection of faecal coliform bacteria.

DESCRIPTION

m-Faecal Coliform Agar is a selective medium used for the detection and enumeration of faecal coliforms in water samples, food contact surfaces and other materials of sanitary importance.

This medium, which is recommended by the APHA, was formulated to be used with the membrane filter technique and to differentiate faecal from non-faecal coliforms.

| TYPICAL FORMULA | (g/l) |
|-------------------------------------|-------|
| Tryptose | 10.0 |
| Proteose Peptone | 5.0 |
| Yeast Extract | 3.0 |
| Sodium Chloride | 5.0 |
| Bile Salts | 1.5 |
| Lactose | 12.5 |
| Methyl Blue (formerly Aniline Blue) | 0.1 |
| Agar | 15.0 |

Final pH 7.4 ± 0.2 at 25°C

METHOD PRINCIPLE

Tryptose and proteose peptone provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium. Bile salts and methyl blue are selective agents effective against the accompanying Gram-positive flora. Methyl blue acts also as growth indicator. Lactose is the fermentable carbohydrate. Agar is the solidifying agent.

TEST PROCEDURE

1. Filter duplicate water samples through separate membrane filters.
2. Place the membranes on the surface of the agar avoiding formation of air bubbles between the filter and the agar surface.
3. Incubate a set of plates at 35 ± 2°C and the other one at 44.5 ± 0.5°C for 24 hours under aerobic atmosphere.

For hygiene monitoring in food handling and manufacturing environments, use a swab and the sampling template 10x10 (ref. 96762) to sample a well defined area of the test surface. If sample areas are irregular develop a standard sampling procedure which is documented and used consistently. After swabbing the test area, inoculate the medium by streaking the swab onto the agar surface. Alternatively, contact plates can be directly used for surface sampling by firmly pressing the agar medium against the test area.

INTERPRETING RESULTS

Colonies produced by fecal coliforms form blue colonies. Coliforms from other sources should be inhibited; however, if present they may appear, as other organisms, as gray to cream-colored colonies.

APPEARANCE

Clear, blue to violet.

STORAGE

Store at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

6 months.

QUALITY CONTROL

The plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Inoculum for selectivity: 10⁴-10⁶ CFU.

Incubation conditions: 24 ± 2 h at 44.5 ± 0.5°C.

QC Table.

| Microorganism | Specification | |
|-------------------------------|---------------|---|
| <i>Escherichia coli</i> | ATCC® 25922 | Good growth, blue to dark blue colonies |
| <i>Salmonella</i> Typhimurium | ATCC® 14028 | Good growth, pink to red colonies |
| <i>Enterobacter aerogenes</i> | ATCC® 13048 | Good growth, gray to gray-blue colonies |
| <i>Enterococcus faecalis</i> | ATCC® 19433 | Inhibition |

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. E.E. Geldreich, H.F. Clark, C.B. Huff, L.C. Best, J. Am. (1965) Water Works Assoc., 57, 208.
2. A.E. Greenberg, R.R. Trussel, L.S. Clesceri (Eds.), Standard Methods for the Examination of Water and Wastewater, 16th ed., APHA Washington D.C.

| PRESENTATION | | Contents | Ref. |
|------------------------|---|-----------|--------|
| m-Faecal Coliform Agar | 90 mm ready-to-use plates | 20 plates | 11210 |
| m-Faecal Coliform Agar | 60 mm ready-to-use plates | 20 plates | 163442 |
| m-Faecal Coliform Agar | 55 mm ready-to-use contact plates (in blister packs) | 20 plates | 15368 |

TABLE OF SYMBOLS

| | | | | |
|-----------------------------|---|---|---|---|
| LOT Batch code |  Keep away from sunlight |  Manufacturer |  Use by |  Fragile, handle with care |
| REF Catalogue number |  Temperature limitation |  Contains sufficient for <n> tests |  Caution, consult Instruction For Use |  Do not reuse |

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