MSA / Chromatic™ MRSA

Bi-plate for selective isolation of pathogen staphylococci and detection of methicillin-resistant *S. aureus* in clinical specimens.

DESCRIPTION

MSA / Chromatic™ MRSA is a ready-to-use plate containing two distinct selective media used for screening of *Staphylococcus aureus*, including methicillin-resistant strains (MRSA), in nasal swabs and soft tissue infections.

Staphylococcus aureus is a Gram positive, coagulase positive cocci with well documented pathogenicity for humans. *S. aureus* has often been implicated in nosocomial infections and food poisoning. Infections affecting MRSA continue to increase especially in hospital settings and represent a real threat to public health.

Mannitol Salt agar (MSA) is the medium recommended by the American Public Health Association (APHA) for the identification, differentiation and enumeration of staphylococci based on their ability to ferment mannitol.

Chromatic™ MRSA contains specific inhibitors and chromogenic substrates that allow the determination of MRSA colonies while the growth of methicillin-sensitive *S. aureus* (MSSA) strains is inhibited.

TYPICAL FORMULA

Mannitol Salt Agar	(g/l)	Chromatic™ MRSA	(g/l)
Pancreatic Digest of Casein	5.0	Peptone and Yeast Extract	30.0
Peptic Digest of Animal Tissue	5.0	Sodium Chloride	10.0
Meat Extract	1.0	Dibasic Sodium Phosphate	2.5
D-Mannitol	10.0	Selective and Matting agents	16.5
Sodium Chloride	75.0	Chromogenic and Antibiotic Mix	0.8
Phenol Red	0.025	Agar	15.0
Agar	15.0	Final pH 7.2 ± 0.2 at 25°C	

Final pH 7.4 ± 0.2 at 25°C

METHOD PRINCIPLE

Mannitol Salt Agar

Pancreatic digest of casein, peptic digest of animal tissue and meat extract provide amino acids, nitrogen, carbon, vitamins and minerals for the growth of microorganisms. Mannitol is the fermentable carbohydrate. The presence of 7.5% sodium chloride inhibits most bacteria except staphylococci. Phenol red is the pH indicator. Agar is the solidifying agent.

Chromatic™ MRSA

Peptone and yeast extract provide amino acids, nitrogen, carbon, minerals, vitamins and other nutrients that support the growth of microorganisms. Sodium chloride has a positive effect on the growth of *S. aureus*. Sodium phosphate is the buffer. Selective agents inhibit the growth of yeasts and most Gram negative and Gram positive bacteria with the exception of methicillin-resistant staphylococci. Matting agents improve the contrast of colonies on the medium. Agar is the solidifying agent. The chromogenic and antibiotic mixture allows the optimal recovery of MRSA and the identification based on the mauve or orange-mauve coloration of colonies.

TEST PROCEDURE

Inoculate the plate by direct streaking over the agar surface. Incubate aerobically at 35 ± 2°C for 18-24 hours.

NOTE: If there is no visible growth on MSA incubate for further 24 hours. ATTENTION, results on Chromatic MRSA must be interpreted after the first 24 h of incubation.

RESULTS INTERPRETATION

On MSA, *S. aureus* produces yellow or white colonies surrounded by a yellow halo. Coagulase-negative staphylococci form tiny colorless to red colonies with no changes in the color of the medium.

Observe the colonies on Chromatic™ MRSA and interpret the results as indicated in the ID table.

ID Table.

Microorganism	Typical colony color
MRSA	Mauve or orange-mauve
Other Gram-positive bacteria (if not inhibited)	White or blue

APPEARANCE

MSA: Slightly opalescent, pinkish-red.

Chromatic™ MRSA: Opaque, whitish.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

4 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Inoculum for selectivity: 10⁴-10⁶ CFU.

Incubation conditions: 35 ± 2°C for 18-24 hours.

QC Table.

Microorganism		Growth on Mannitol Salt Agar	Growth on Chromatic™ MRSA
<i>Staphylococcus aureus</i> (MRSA)	ATCC® 43300	Good, yellow colonies with yellow halo	Good, mauve colonies
<i>Staphylococcus aureus</i> (MSSA)	ATCC® 25923	Good, yellow colonies with yellow halo	Inhibited
<i>Staphylococcus epidermidis</i>	ATCC® 12228	Good, white colonies on red medium	Inhibited
<i>Escherichia coli</i>	ATCC® 25922	Inhibited	Inhibited

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used by properly trained operators only.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulation in force.


BIBLIOGRAPHY

1. Evaluation of CHROMagar Staph aureus, a new chromogenic medium, for isolation and presumptive identification of *Staphylococcus aureus* from human clinical specimens. Gaillot O. et al. (2001) Journal of Clinical Microbiology, 38 : 1587-1591.
2. Dépistage nasal de *Staphylococcus aureus*. Nécessité de standardiser les protocoles. Laudat P. et al. (2000) Poster 343/P2 presented at RICAI in Paris (France).
3. Kloos, W.E., and T.L. Bannerman (1995) *Staphylococcus* and *Micrococcus*. In Manual of clinical microbiology, 6th ed.
4. Chapman, G.H. (1945) The significance of sodium chloride in studies of staphylococci. J. Bacteriol. 50:201-203.

PRESENTATION

		Packaging	Ref.
MSA / Chromatic™ MRSA	90 mm ready-to-use plates	20 plates	18024

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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