Columbia CNA / Chocolate

Bi-plate for the isolation, cultivation and differentiation of Gram-positive and Gram-negative bacteria from clinical specimens.

DESCRIPTION

Columbia CNA / Chocolate is a ready-to-use plate containing two distinct media used for the detection of fastidious and non-fastidious microorganisms from clinical specimens and other materials of sanitary importance.

Columbia CNA Agar is a selective medium for the isolation and differentiation of Gram-positive bacteria, especially Streptococci and Staphylococci.

Chocolate Agar is a non-selective medium for the cultivation of fastidious organisms, primarily *Neisseria* and *Haemophilus* spp.

TYPICAL FORMULA

Columbia CNA Agar	(g/l)	Chocolate Agar	(g/l)
Enzymatic Digest of Casein	12.0	Enzymatic Digest of Casein	7.5
Enzymatic Digest of Animal Tissues	5.0	Enzymatic Digest of Animal Tissues	7.5
Beef Extract	3.0	Corn Starch	1.0
Yeast Extract	3.0	Dipotassium Phosphate	4.0
Corn Starch	1.0	Monopotassium Phosphate	1.0
Sodium Chloride	5.0	Sodium Chloride	5.0
Agar	15.0	Agar	17.0
Colistin	0.01	Horse Blood, Defibrinated	50.0 ml
Nalidixic Acid	0.015	Vitalex Growth Supplement	10.0 ml
Sheep Blood, Defibrinated	50.0 ml	Final pH 7.2 ± 0.2 at 25°C	
Final pH 7.3 ± 0.2 at 25°C			

METHOD PRINCIPLE

Enzymatic digest of casein, enzymatic digest of animal tissues and beef extract provide amino acids, nitrogen, carbon, vitamins and minerals required for organisms growth. Yeast extract is a source of vitamins, particularly of group B. Corn starch neutralizes any toxic metabolites. Sodium Chloride maintains the osmotic balance of the medium. The phosphates are buffering agents. Agar is the solidifying agent.

Colistin and nalidixic acid incorporated into Columbia CNA Agar suppress the growth of *Proteus*, *Klebsiella* and *Pseudomonas* while Staphylococci, haemolytic Streptococci and Enterococci still growing. Sheep blood allows detection of hemolytic reactions

Horse blood and Vitalex Growth Supplement in Chocolate Agar, enhance recovery of *Haemophilus* and *Neisseria* spp.

TEST PROCEDURE

Inoculate the plates by direct streaking or spread method. Incubate at 35 ± 2°C for 18-48 hours in an aerobic or CO₂-enriched atmosphere.

RESULTS INTERPRETATION

Examine for bacterial growth and interpret the results on Columbia CNA Agar as indicated in the ID table.

ID Table.

Microorganism	Appearance of colonies
Streptococci (non-group D)	Small, white to grayish, beta or alpha hemolysis
Enterococci (group D)	Small, but larger than group A Streptococci, blue-gray, beta or alpha hemolysis
Staphylococci	Large, white to gray or cream to yellow, with or without hemolysis
<i>Listeria monocytogenes</i>	Blue-gray with beta hemolysis

Alfa hemolysis: greenish coloration around the colony

Beta hemolysis: clear zone around or under the colony

APPEARANCE

Columbia CNA Agar: opaque, red.

Chocolate Agar: opaque, chocolate brown.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

SHELF LIFE

2 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU

Inoculum for selectivity: 10⁴-10⁶ CFU

Incubation conditions: 35 ± 2°C for 18-48 hours in aerobic atmosphere supplemented with 5-10% carbon dioxide.

QC Table.

Microorganism		Growth on Columbia CNA Agar	Growth on Chocolate Agar
<i>Staphylococcus aureus</i>	ATCC® 25923	Good	Good
<i>Streptococcus pneumoniae</i>	ATCC® 6305	Good, α-hemolysis	Good
<i>Streptococcus pyogenes</i>	ATCC® 19615	Good, β-hemolysis	Good
<i>Proteus mirabilis</i>	ATCC® 12453	Inhibited (partially)	Good
<i>Haemophilus influenzae</i>	ATCC® 10211	Inhibited	Good
<i>Neisseria gonorrhoeae</i>	ATCC® 43069	Inhibited	Good

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulation in force.

BIBLIOGRAPHY

- Clinical Microbiology Procedures Handbook, 3rd ed. (2010) ASM. Washington, D.C. 2,540 pages.
- Murray, P. R., E. J. Baron, J. H. Jorgensen, M. L. Landry, and M. A. Pfaller (ed.) Manual of Clinical Microbiology, 9th ed. (2007) ASM Press, Washington, D. C
- Isenberg, H.D. (ed.). 2004. Clinical microbiology procedures handbook, vol. 1, 2 and 3, 2nd ed. American Society for Microbiology, Washington, D.C.
- Chapin, C.K., G.V. Doern (1983). J. Clin. Microbiol. 17: 1163-1165.
- Thayer J.D., R.L. Jackson (1975) Public Health Rep. 82:361.
- Martin, J.E., Armstrong J.H., Smith P.B. (1974). Appl. Microbiol. 27: 802-805.
- Seth, A. (1970). Brit. J. Vener. Dis. 46: 201-202.

PRESENTATION		Packaging	Ref.
Columbia CNA / Chocolate	90 mm ready-to-use plates	20 plates	18507
Columbia CNA / Chocolate	90 mm ready-to-use plates	100 plates	18507*

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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