Muller Kauffmann Tetrathionate Novobiocin Broth (MKTTn)

Selective enrichment medium for detection of *Salmonella* spp,
according to ISO 6579-1.

DESCRIPTION

Muller Kauffmann Tetrathionate Novobiocin Broth (MKTTn) is a liquid medium used for the selective enrichment of *Salmonella* spp, including *S. Typhi* and *S. Paratyphi*.

This medium is recommended by ISO 6579-1 for the examination of food, animal feed and environmental samples from the food production area.

TYPICAL FORMULA	(g/l)
Meat Extract	4.3
Enzymatic Digest of Casein	8.6
Sodium Chloride	2.6
Calcium Carbonate	38.7
Sodium Thiosulfate	30.5*
Ox Bile	4.78
Brilliant Green	0.0096
Iodine	4.0
Potassium Iodide	5.0
Novobiocin	0.04

Final pH 8.0 ± 0.2 at 25°C

(*) Equivalent to 47.8 g of sodium thiosulfate pentahydrate.

METHOD PRINCIPLE

Meat extract and enzymatic digest of casein provide amino acids, carbon, nitrogen, vitamins and minerals for organisms growth. Sodium chloride maintains the osmotic balance of the medium. Calcium carbonate neutralizes toxic metabolites. Sodium thiosulfate react with iodine and potassium iodide to form tetrathionate, which inhibit the contaminant bacterial flora. Ox bile and brilliant green are additional selective agent used to suppress coliform bacteria and Gram-positive bacteria. Novobiocin is mainly effective against *Proteus* spp.

TEST PROCEDURE

Dilute food sample 1:10 in Buffered Peptone Water (ref. 414020), e.g. a 25 g test portion is mixed with 225 ml of BPW. Incubate at 34-38°C for 18 ± 2 h.

Transfer 1 ml of the pre-enrichment culture to a tube containing 10 ml of MKTTn broth. Incubate at 37 ± 1°C for 24 ± 3 h. An additional 24 h incubation may be required.

INTERPRETING RESULTS

Turbidity indicates microbial growth.

Use this enrichment culture to inoculate XLD Agar plates (ref. 10056) and a second different selective solid medium. Incubate aerobically at 37 ± 1°C for 24 ± 3 h. Examine for typical colonies of *Salmonella* spp which have a black centre and a lightly transparent zone of reddish color on XLD agar.

APPEARANCE

Very pale green, opaque, heterogeneous with precipitate.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Medium in bottles: 2 years.

Medium in tubes: 1 year.

QUALITY CONTROL

The medium is inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: ≤ 100 CFU.

Inoculum for selectivity: $> 10^3$ CFU.

Incubation conditions: $37 \pm 1^\circ\text{C}$ for 24 ± 3 hours.

QC Table.

Microorganism	Specifications
<i>Salmonella</i> Typhimurium	WDCM 00031 >10 characteristic colonies on XLD agar
+ <i>Escherichia coli</i>	WDCM 00013
+ <i>Pseudomonas aeruginosa</i>	WDCM 00025
<i>Escherichia coli</i>	WDCM 00013 Partial inhibition, ≤ 100 colonies on Tryptic Soy Agar
<i>Enterococcus faecalis</i>	WDCM 00009 <10 colonies on TSA

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

- ISO 6579-1:2017. Microbiology of the food chain – Horizontal method for the detection, enumeration and serotyping of *Salmonella* – Part 1: Detection of *Salmonella* spp.
- EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- Kauffmann F. (1935) Weitere Erfahrungen mit dem kombinierten Anreicherungsverfahren für Salmonellenbacillen. Z. Hyg. Infekt. 117:26-32.
- Muller L. (1923) Un nouveau milieu d'enrichissement pour la recherche du bacille typhique et des paratyphiques. Comp. rend. Soc. biol. 89:434-437.

PRESENTATION	Category	Packaging	Ref.
Muller Kauffmann Tetrathionate Novobiocin Broth (MKTTn)	Tubes - Bottles	20 x 10 ml tubes	20072
Muller Kauffmann Tetrathionate Novobiocin Broth (MKTTn)	Tubes - Bottles	100 x 10 ml tubes	26072
Muller Kauffmann Tetrathionate Novobiocin Broth (MKTTn)	Tubes - Bottles	6 x 100 ml bottles	403110

TABLE OF SYMBOLS

LOT Batch code	 Do not reuse	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	

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